


Weight Bearing Protection



Our new Weight Bearing Protective Cover has been designed to accommodate Digital DR Panels and CR Cassettes.

Our covers are strong yet lightweight, offering a high degree of protection during weight bearing procedures.

Our latest cover is constructed of a high density polyethylene frame and a carbon fiber top, offering protection up to 1300 lbs. Our easy lift handle makes removal and carrying a breeze.

The Weight Bearing Cover has been tested by an accredited independent laboratory.



Custom X-Ray Digital
Equipment Sales & Service

2120 W. Encanto Blvd., Phoenix, AZ 85009
[602] 439-3100 | [800] 230-XRAY | Fax: [602] 252-5950
www.customxray.com

Weight Bearing Protection



Smoked Polycarbonate



Carbon Fiber

- Available in Polycarbonate or Carbon Fiber
- Rated up to 1300 lbs.
- Easy lift handle
- Non-skid rubber feet
- Multiple sizes available
- Printed target area

This product is perfect for bariatric, equine, and podiatry procedures.



Easy Lift Handle

Protect your *investment* use a **COVER!**

Custom X-Ray Digital
Equipment Sales & Service

2120 W. Encanto Blvd., Phoenix, AZ 85009
[602] 439-3100 | [800] 230-XRAY | Fax: [602] 252-5950
www.customxray.com