



INNOVATIVE CARE SOLUTIONS


Custom X-Ray Digital
Equipment Sales & Service

2019-2020
Physical Therapy
Equipment Catalog

customxray.com
800.230.XRAY

ELECTRIC TREATMENT TABLES

6 TOP CONFIGURATIONS



PT1101 PT1102 PT1103 PT1107 PT1105 PT1106



Interested in branding
your equipment
with your school or
company logo?

Contact us today for your
FREE mockup.

800.743.7738



PT1103

FEATURES:

- 77" L x 27" W
- Height: 17" - 35"
- Lifting capacity: 400 lbs.
- 2" high density firm foam with rounded corners
- Single foot pedal
- Low profile locking casters
- White base powder coat
- Whisper quiet lift system

AVAILABLE OPTIONS:


- 30"W cushion
- Dual foot bar
- 650 lb. lifting options

THERA-P TABLE

Ease of use, multiple options and a sleek but durable design make the **THERA-P** line of electric treatment tables the choice of therapists everywhere. The **THERA-P** electric treatment table is constructed of precision-cut flat and tube steel for durability and style. It also has a whisper quiet lift system (with elevating foot pedal) and easy-to-operate caster system for better ergonomics.

BARIATRIC VERSION

available with 40" wide cushion



PT1113

LAST TABLE™

LEG & SHOULDER THERAPY TABLE™

The **LAST Table™** is the first ergonomically designed treatment/exercise table that permits bilateral training and specific anatomical positioning necessary to achieve desired results. The distinct Hi-Lo elevating base allows individuals with limited mobility access to the table in any setting.

Patent Pending.

FEATURES:

- 84" L x 30" W
- Height: 19" - 37"
- Lifting capacity: 400 lbs.
- Adjustable top with standard elevating base
- Elevating back 0 - 70°
- Elevating legs 0 - 90°
- Leg height of 17"
- Drop legs 0 - 90°
- 2" high density firm foam with rounded corners
- Single foot pedal
- Low profile locking casters
- White powder coat base



CAB-100-LST
with options

LEG & SHOULDER TREATMENT CABINET

Visit our website for more details:
PHSmedicalsolutions.com



PT9060
with options

ME2000

ELEVATING TREATMENT TABLE

ME2000
with options

This dynamic table was designed and built with the highest standards of quality, integrity and innovation. With its quiet and smooth elevation, the **ME2000 Elevating Treatment Table** eliminates the constant squatting and bending while treating patients, leaving you with more energy at the end of the work day. This table offers maximum storage with access on each side of the table.

FEATURES:

- 77" L x 28" W
- Height: 24" - 36"
- Lifting capacity: 400 lbs.
- Complete cabinet with 2 sliding doors on each side of the table and 1 foot pedal
- Flat 2" high density cushion with rounded corners and nose hole
- Lift back option available

ME2001

BARIATRIC ELEVATING TREATMENT TABLE

FEATURES:

- 84" L x 40" W
- Height: 24" - 36"
- Lifting capacity: 400 lbs.

With a wider top, the ME2001 Bariatric Elevating Treatment Table will give you the flexibility to work with patients of all sizes without compromising comfort.



ME2001

TREATMENT TABLES

19407



FEATURES:

- 30"W x 73"L x 31"H
- Load Capacity: 450 lbs.
- Solid hardwood legs
- H-brace construction
- Shelf • Rounded corners and nose hole
- 2" High density foam cushion
- 70° Manual multi-position elevating lift back

BASIC WOOD TREATMENT TABLE

Meet the next generation of cost-effective exam tables. Our **Basic Wood Treatment Table** is built to meet your needs and budget. Each table comes standard with H-brace construction, shelf, comfortable 2" high density foam cushion, 70° manual lift back with multi-position elevation and softly rounded corners. This table is the perfect fit for any treatment space. Some assembly required.

ESSENTIAL WOOD TREATMENT TABLE

Each table comes standard with a poplar hardwood frame, H-brace construction, comfortable 2" high density foam cushion, softly rounded corners and dual-bolted folding leg assembly for easy shipping. Make it your own by adding length, width, shelves or a cabinet that can be mounted on either end.



TEWT
with options

FEATURES:

- 27" or 30"W x 72" or 78"L x 31"H
- Load Capacity: 500 lbs.
- Solid poplar hardwood frame
- H-brace construction
- 2" High density foam cushion
- Rounded corners
- Dual-bolted folding leg assembly for convenient shipping
- Lift back option available



CLASSIC WOOD TREATMENT TABLE



TESS
with options

FEATURES:

- 27" or 30"W x 72" or 78"L x 31"H
- Load Capacity: 550 lbs.
- Solid hardwood H-brace construction
- 2" High density foam cushion
- Rounded corners
- Dual-bolted leg assembly that folds down for convenient shipping
- Lift back option available

Custom build your wood treatment table to meet the specific needs of your practice. Make it your own by adding length or width, shelves or a cabinet that can be mounted on either end. Patient-friendly customizations include a nose hole and 70° manual lift back with 6 position elevation.



FESS-6L

Designed for greater stability, our **6 Leg Classic Wood Treatment Table** offers a rich variety of standard options as well as customizations. Tailor your table to your needs with an additional cabinet that can be mounted on either end or include a nose hole, 70° manual lift back with 6 position elevation and 45° split leg elevation.

CLASSIC 6 LEG WOOD TREATMENT TABLE

FEATURES:

- 30"W x 78"L x 31"H
- Load Capacity: 800 lbs.
- Solid hardwood 6 Leg construction
- 2" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly
- Bottom shelf for added storage
- Lift back and split leg options available



PT9080

HIP & KNEE FLEXION TREATMENT TABLE

The flexible split leg design of this solid hardwood treatment table makes it extremely functional for multiple treatment positions. Each padded leg cushion can be manually adjusted to the desired hip and knee flexion positions making treatments easier and more comfortable for the patient and therapist.

FEATURES:

- 30"W x 81.5"L x 32"H
- Load Capacity: 800 lbs.
- Solid hardwood construction
- Dual brace support
- 2" High density foam cushion
- Rounded corners
- Manual lift back
- Elevating Back 0-70°
- Elevating Legs 0-90°
- Hip/Knee flexion
- Bottom shelf for added storage

**A4327LB
with options**



ALUMA ELITE BASIC TREATMENT TABLE

The dependable **Aluma Elite Basic Treatment Table** is the perfect centerpiece for every therapy or rehab facility. The extraordinary load capacity of 650lbs. and the welded aluminum frame make the Aluma Elite Basic Treatment Table the perfect table for your office, clinic, gym, or facility. Add length or width, laminate shelves or a cabinet that can be mounted on either end to meet your specific needs. Other customizations include a nose hole and 70° manual lift back.

FEATURES:

- 27" or 30"W x 72" or 78"L x 31"H
- Load Capacity: 650 lbs.
- Welded aluminum frame
- H-brace construction
- 2" High density foam cushion
- Rounded corners
- Tough powder coat finish
- Lift back option available

**BRANDING
AVAILABLE**



COMBINATION LAB TABLE/DESK

FEATURES:

- 27" W x 72" L x 30" H
- Load capacity: 450 lbs.
- 4 leg, solid hardwood construction
- Flat, laminate top
- 24" clearance for stools



FESS-LAB

The **Combination Lab Table/Desk** is a desk or treatment table ideal for classroom training. The table is constructed of solid hardwood with a flat, laminate top. The table can be built using oak or maple wood. Choose from nine standard colors for the laminate top. Some on site assembly required.

CLASSROOM/LAB TREATMENT TABLE

FEATURES:

- 27" W x 78" L x 32" H
- Load capacity: 600 lbs.
- 6 leg, solid hardwood construction
- Flat, laminate top
- 24" clearance for stools
- 2" Vinyl/foam pad available



FESS-LAB-6L with options

The **Classroom/Lab Treatment Table** is a multi-purpose desk and treatment table perfect for a clinic or classroom setting. The table is constructed of solid hardwood with a flat, laminate top. The table can be built using oak or maple wood. Add storage to your table by adding 2 open cubbies. An optional removable 2" vinyl/foam pad attaches to the laminate table top using hook and loop fastener. Some on site assembly required.

WALL MOUNT FOLDING TREATMENT TABLE



FEATURES:

- 24" W x 72" L x 31" H
- Load capacity: 450 lbs.
- Solid hardwood construction
- Flat, seamless 2" high density foam cushion
- Durable wall mounting plate

The **Wall Mount Folding Treatment Table** is a fold-out treatment table ideal for small spaces. The table is constructed of solid hardwood with a flat, seamless 2" high density foam cushion. Three legs attach to the wall, while the two front legs and cushion fold out away from the wall when needed. A durable wall mounting plate is included (wall fasteners not included).

WMFTT-2472



FESS-LAB-6L with options

CUSHION BRANDING
available



FESS-LAB with options

CUSTOM LAMINATE COLORS &
CUSTOM PRINTED LAMINATE

available on orders of 10
or more tables.

MAT TABLES



PT2000

PT2000 BARIATRIC HILO MAT TABLE

FEATURES:

- Available in 4'x7', 5'x7' and 6'x8' configurations
- 18"-33" H
- Load Capacity: 1000 lbs.
- Single foot pedal
- Drop down casters
- Durable steel base with powder coat finish
- 1.5" High density foam cushion
- Rounded corners

The **PT2000 Bariatric Hi-Lo Mat Table** is ideal for rehabilitation settings and is available in a range of sizes to provide your facility with the option of delivering optimum care for larger patients. The table features 18-33" height range, single foot pedal, 1.5" high density foam cushion and convenient drop down casters.



**ALSO AVAILABLE:
PT2001 BARIATRIC
MAT TABLE**
(without casters)



WALL MOUNT MAT TABLE

available in
4'x7' or 5'x7'
configurations

WMMT



CLASSIC WOOD MAT TABLE

MAT-4872-4L

FEATURES:

- Sizes available: 4'x6', 4'x7', 5'x7' and 6'x8'
- Table height: 19"
- Load capacity: 450 lbs.
- Solid hardwood frame and 4 legs
- 3/4" Plywood top
- 1.5" High density foam cushion
- Rounded corners
- Triple-bolted leg assembly

The **Classic Wood Mat Table** is built for strength, stability exercise and other treatments. Our 4 Leg Mat table comes standard with solid hardwood on the base frame, four legs and supports with a thick plywood top construction. The mat table also features a layer of comfortable 1.5" high density foam with softly rounded corners and triple-bolted leg assembly for extra durability.

6 leg option available.



A4520

ALUMA ELITE MAT TABLE

FEATURES:

- Available in 4'x7', 5'x7' and 6'x8' configurations
- Table height: 18" H
- Load Capacity: 4'x7' - 800 lbs.
5'x7' & 6'x8' - 1200 lbs.
- Welded frames & leveling feet
- Tough powder coat finish
- 1.5" foam, semi-firm

The **Aluma Elite Mat Table** is one of the most dependable Aluma Elite products with a load capacity of 800 lbs. to 1,200 lbs. Built with durable aluminum, the Mat Table can handle the harshest environments and withstand the rigors physical therapy entails.



SMC-002
with options

SMC-001

SMDC-003
with options

STOR-EDGE™ MODALITY CARTS

Small with ample storage, our **Stor-Edge Modality Carts** come in standard and deluxe versions, with multiple drawer and shelf options that keep tools and treatment supplies close so therapists can work more efficiently.

FEATURES:

- Available in wood or laminate finishes
- Black wire pull handles
- Carts ship fully assembled, except for casters
- Full extension glides on all drawers
- 2" Swivel nylon casters



SC-009D

SC-003

STOR-EDGE™ STATIONARY CABINETS

The stylish **Stor-Edge Stationary Cabinet** system provides multi-purpose storage in a variety of options, including tall cabinets, shelves or drawers and more - ideal for storing gear and other supplies inside a treatment space or office.

FEATURES:

- Full extension glides
- Wood and laminate finishes available
- Cabinets available with or without exterior double doors
- Weight racks and slatwall hooks provided where applicable
- Black wire pull handles



STOR-EDGE™ MEDICAL CABINETS

Modular, functional, durable and aesthetically pleasing, our lines of **Stor-Edge Medical Wall & Base Cabinets** are cost-effective and suited for multiple applications.

FEATURES:

- Base cabinets: 18" or 24"D x 36"H
- Wall cabinets: 12" or 18"D x 24" or 30"H
- Must specify left (L) or right (R) hinge on single door cabinets
- Countertop included
- Cost effective, yet durable for multiple applications
- Available in wood or laminate finishes
- Full extension drawer glides on all drawers
- All cabinets come standard with (1) adjustable shelf
- Black wire pull handles



CAB-090N

CAB-010
with options

CAB-120

TREATMENT CABINETS

Your practice will shine with the addition of a new **treatment cabinet**. Our line of cabinets meet the functional needs of your practice, with ample storage, 55° lift back with manual gas shock and 45° split leg options.

STOR-EDGE™

MULTI-PURPOSE CARTS

Crafted of wood or laminate, our **Stor-Edge Multi-Purpose Carts** are custom built. Choose from many model choices for the right solution to organizing your weights, resistance bands and other tools and equipment.

FEATURES:

- Available in wood or laminate finishes
- Full extension glides on all drawers
- Black wire pull handles
- Durable slatwall construction
- Weight racks and J-hooks provided where applicable
- 4" heavy duty casters

**Carts shown with accessories (not included).
Weight kits are available.**



SR-009
with options



SR-005
with options



SR-011



SR-001



PTST-005-WHR

CONVERTIBLE STAIRCASES

FEATURES:

- 350 lbs. Load rating
- Stairs can be configured as a straight or corner unit
- Platform: 36"W x 36"D x 24"H
- (4) 36"W x 9"D x 6"H steps and (6) 36"W x 9"D x 4"H steps
- 9' ceiling height required

PTST-002-WHR



STRAIGHT STAIRCASES

FEATURES:

- 350 lbs. Load rating
- Platform: 30"W x 24"D x 24"H
- (4) 30"W x 9"D x 6"H steps
- 9' ceiling height required

PTST-001-WHR-ADA



PASS THROUGH STAIRCASES

FEATURES:

- 350 lbs. Load rating
- Platform: 30"W x 24"D x 12"H
- (2) 30"W x 9"D x 6"H steps and (3) 30"W x 9"D x 4"H steps
- 8' ceiling height required

ALL STAIRCASES FEATURE:

- Handrails made from solid hardwood (shipped ready to assemble)
- Anti-slip safety tread on platform and all steps for increased safety
- ADA compliant staircases are available with 11"D steps and handrails with 12" of overhang
- Uprights made from 2" steel with a black powder coat finish
- Stairs available in wood or laminate finishes



Hands Free Electric Table

CREATE YOUR PIVOTAL PRACTICE

Every PHS Medical product is made to help your facility succeed. From our treatment tables to our staircases, Hi-Lo tables and modality carts, **our product innovation doesn't end there.**



Other brands of Pivotal Health Solutions:



DOC Decompression Table



The Sport Portable Aluminum Treatment Table



Be the Difference.
 pivotalhealthsolutions.com

Innovative Care Solutions - PHSmedicalsolutions.com



AP2-DT-5000-C

APOLLO DESKTOP LASER UNIT

2 Channel Desktop with choice of
3000, 4000 or 5000mW probe

9"W x 11.25"D x 4.5"H

FEATURES:

- Heavy-duty construction for office and clinic environments
- Probe holders on both sides of the unit that can accommodate two probes and both right and left handed practitioners
- Elevated face for easy viewing and control

ROLLING STOOLS




Our rolling stools combine comfort with control thanks to a 4" thick deluxe seat cushion, a five-star base with heavy-duty casters and pneumatic controls that adjust seat height from 16"-22".

FEATURES:

- 16"W x 16"D x 16-22"H
- Load Capacity: 300 lbs.
- Padded vinyl seat
- 4" deluxe foam cushion
- Heavy duty casters
- Swivel stool

**VISIT OUR WEBSITE
FOR BOLSTERS &
OTHER ACCESSORIES**



FMST TOOLS



Fluid Motion Soft Tissue Tools are the soft tissue mobilization tools that can take myofascial manipulation and mobilization to the next level in your practice. The design of the two tools available enables you, the practitioner, to work over all surfaces of the body quickly and easily. The tools help reduce scar tissue problems from old injuries/surgeries and achieve greater results with pain control, joint and muscle flexibility.

FEATURES:

- Custom molded, non-slip thumb grips create maximum angular pressure without over use by the therapist.
- Ambidextrous design: no more switching through multiple tools to find the right edge or angle. Both tools are designed to be used on both sides and on every surface for increased compatibility of treatment on multiple body areas.

Interested in branding
your equipment
with your school or
company logo?

Contact us today for your
FREE mockup.

800.743.7738



We are continuously developing, evolving
and improving products so we can deliver
unmatched value to our customers
in both product and service. It's our
priority to be the difference in providing
innovative care solutions.



phsmedicalsolutions.com