


CONVERT YOUR GE AMX4 AND AMX4+ TO A **WIRELESS DR UNIT**

with Medlink's AMX Wireless DR Upgrade Kit

INCREASE PATIENT THROUGHPUT AT A REASONABLE PRICE

- Fraction of the cost compared to the expense of a new system
- Low-Cost Service & Maintenance
- Minutes-to-Install
- Superior Image Quality
- Completely Self Contained
- Utilizing Award-Winning Viewing Software

“Medlink’s AMX Wireless Upgrade Kit is the fastest, most reliable way to convert your aging mobile unit into a modern Wireless DR System.”



AMX UPGRADE SYSTEM

Software and Hardware Features

Save Time | Save Money | Manage Assets
WIRELESS DR UPGRADE

Medlink's AMX Upgrade Kit is designed to quickly, easily, and inexpensively turn your current portable X-ray unit into a wireless solution allowing to maximize patient thorough-put.

- RECEIVE IMMEDIATE IMAGE VISUALIZATION
- EASY INSTALLATION
- ERGONOMICALLY DESIGNED
- ADVANCED AED TECHNOLOGY (AUTOMATIC EXPOSURE DETECTION)
- DICOM 3.0 NETWORK INTERFACE
- AWARD-WINNING VXVUE VIEWING SOFTWARE



Mast Mount

MAST MOUNT

- Workstation Bracket Assembly with 270 degrees of rotation
- Dual Band Access Point



All-In-One Touchscreen PC

ALL-IN-ONE TOUCHSCREEN PC

- 17-inch (All-in-One) Medical Grade PC, (1280 x 1024 resolution)
- 15-Wire Resistive Touch Screen
- Intel Dual Band Wireless-AC 7260
- Touchscreen Stylus Pen
- 8GB DD R4-2400MHz SO-DIMM
- **Three (3) Year Warranty**

VIVIX WIRELESS DR PANEL with ACQUISITION SOFTWARE



VIVIX Wireless Panel and VXVUE Acquisition Software

- Scintillator: CsI
- Pixel Matrix: 2560 x 3072 pixels
- Pixel Pitch: 140µm
- Image Size: 14" x 17"
- Grayscale: 14bit
- Image Acquisition Time: 2 seconds
- Data Interface: Gigabit Ethernet | IEEE 802.11n
- Battery: Lithium Ion 4000 mAh - 4 h(capturing)
- Communication: Inside AP™ Wireless
- Generator Interface: Auto Trigger AED (Automatic Exposure Detection) Mode
- **Five (5) Year Warranty**