



# Apron Racks

2020  
PRODUCT  
GUIDE

## Mobile Apron Z-Rack

Z-Base mobile rack nests with others when not in use. Heavy duty casters on 12G steel base. Chromed-steel posts and hang bar. 60" long x 60" tall bar allows hanging without creasing. Ships @ 35 lbs.

**Rack ONLY**  
#ARM-UZNH  
\$280.00

OR

**Rack with Ten Hangers**  
#ARM-UZ10  
\$566.00



## NEW Mobile Apron Z-Rack Middle Bar Add-On

Middle Bar addition to hang vests and skirts! 63" bar will hold up to 150 lbs. (Mobile apron z-rack sold separately)


**Bar Add-on**  
#ARM-UBAR  
\$58.00



## NEW Mobile Apron Z-Rack Storage Base Add-On

Bottom shelf add-on for storing misc. protective accessories! 62" shelf will hold up to 150 lbs. (Mobile apron z-rack sold separately)

**Shelf Add-on**  
#ARM-UBASE  
\$195.00



## Apron Hangers

Designed for standard aprons, panoramic, vests, skirts, etc. Top hook swivels for easy hanging. 21" long chromed-steel with anti-slip tips. Use on ARM-UZNH rack, in lockers, or behind a door. Single ships @ 2 lbs.

**Single Hanger**  
#ARH-01  
\$29.00



**Set of 10 Hangers**  
#ARH-10  
\$286.00



## Mobile Apron Rack

Ten (10) swing-arm mobile garment stand. Durable all-steel construction. 24" x 24" H-Base for stability and easy cleaning. E-Z-Roll swivel locking casters. Convenient handle for easy mobility. White enamel finish. Overall height is 55". Some assembly required. Ships @ 110 lb.

### Standard Mobile Rack

#ARM-100  
\$1165.00



E-Z Roll, 3" medical grade, swivel, locking casters!

### Standard WITH Glove Rack

#ARM-G60  
\$1240.00



Arms swing inward to conserve space!

Easy-Clean H-Base!



(Aprons sold separately)

## 5 Arm Apron Rack

Five (5) swing-arm mobile garment stand. Durable all-steel construction with white enamel finish. 20" x 20" H-Base for stability and easy cleaning. E-Z roll swivel locking casters. Convenient handle for easy mobility. Arms extend to 29.5" outward. We offer the ability to choose whether you need right or left swing arms. Ships @ 80 lbs. motor freight.

#ARM-105  
\$880.00



You choose the right or left side!



Arms swing inward to conserve space!

E-Z Roll, 3" medical grade, swivel, locking casters!

Easy-Clean H-Base!



(Aprons sold separately)

### Multi-Tubular Apron Rack

The perfect solution safe apron storage. Durable all-steel construction. A 24" deep by 32" wide" H-base offers stability and easy cleaning. Features E-Z Roll, 3" medical grade swivel locking casters. Convenient handle for easy mobility from room to room. White baked enamel finish. Overall height is 49" and width is 32". Rods are 27" long and spaced 4" apart for easy clearance and storage. Stores up to 5 aprons in tight spaces without crinkling or damaging aprons.

#ARM-500 | \$950



E-Z Roll, 3" medical grade, swivel, locking casters!



Convenient handle for easy mobility and storage!

Easy-Clean H-Base!



Holds up to 5 Aprons at once!

(Aprons sold separately)

#### Available Steel Colors

White



Default

Black



Textured

Platinum



Textured

Almond



Textured

## NEW Mobile Combination Protection Rack

The ultimate option for serious radiation protection storage. This versatile rack holds 6 aprons, 3 pairs of gloves, and 6 pairs of radiation protection glasses. Comes with 6 heavy duty apron hangers for quick apron access. Durable all-steel construction. A 24" deep by 28 wide" H-base offers stability and easy cleaning. Features E-Z Roll, 3" medical grade swivel locking casters. White baked enamel finish. Overall height is 79", width is 30" and depth is 24". The perfect height to hang aprons avoiding damage from crinkling. Soft rubber rings and 6 high grade mounts create the ultimate protection for radiation glasses at all times.

#ARM-600 | \$950



Available Steel Colors			
White	Black	Platinum	Almond
			
Default	Textured		Textured



Holds up to 6 aprons, 3 pairs of gloves and 6 pairs of protective glasses at once!

### NEW Mobile Lead Apron Locker

These industrial strength storage cabinets will feature heavy-duty casters for mobility and have 22 gauge welded steel construction. Includes a full width adjustable shelf and a full width hanger rod. The 230 lb. capacity shelf meets ANSI/BIFMA shelf capacity standards. Double doors lock with 3 point latching and recessed handles that prevent banging and catching with a flush surface. Powder coat finish. The 5" polypropylene swivel casters (2 with brakes) are attached with a 6"W heavy-gauge steel plate that bolts through to 1-1/4" angle iron inside the cabinet, both running the length of the cabinet. Steel caster components are zinc treated to resist corrosion. Total caster capacity is 1100 lbs. Rod capacity 500 lbs. Locks with key. Available in charcoal and gray. Item will shipped assembled.



Features	
Width	46
Depth	24
Height	78
Adjustable Shelves	1
Construction	20/22 Gauge Steel
Door Quantity	2
Handle Type	Recessed
Shelf Capacity	230 lbs.
Wheel Type	Polypropylene 4 Swivel, 2 Brake

#### Charcoal Locker Only

MLALC | \$1330

#### Gray Locker Only

MLALG | \$1330

#### Locker & Hanger Bundles

MLALG-10 | \$1595

MLALG-10 | \$1595

Locker in Gray



\*Aprons sold separately



### Wall Mounted Apron Racks

Ten or five swing-arm wall-mount apron racks. Mounting holes on 16" centers. All-steel durability with white enamel finish.



10-Arm #ARW-10 \$436.00



5-Arm, Swing Right #ARW-R5 \$260.00



5-Arm, Swing Left #ARW-L5 \$260.00

### Glove Adapter

Retrofits to Mobile Rack. Holds six (6) 18" gloves.

#ARG-06  
\$103.00



### Apron/Glove Racks

Designed for one apron, but will hold more. All steel, heavy duty construction and a white finish. 16" & 24" holes match most wall studs. Glove holders mounted onto back panel. 28.5" long x 1.75" diameter bar prevents creasing.

#### WITHOUT GLOVE HOLDERS

#AGR-105  
\$187.00



#### WITH TWO GLOVE HOLDERS

#AGR-110  
\$251.00



#### WITH FOUR GLOVE HOLDERS

#AGR-115  
\$272.00



## Apron Peg Racks

Help extend the life of your aprons. Our apron racks are just the thing for you! 6" angled pegs keep garments from slipping. All-steel white enamel finish.

Custom sized racks may be available!

**2-Peg**  
(for 1 apron)  
**#ARP-02**  
\$72.00



**4-peg**  
(for 2 aprons)  
**#ARP-04**  
\$124.00



6-peg and 8-peg racks: mounting holes on 16" and 24" centers.

**6-peg**  
(for 3 aprons)  
**#ARP-06**  
\$145.00



**8-peg**  
(for 4 aprons)  
**#ARP-08**  
\$156.00



**16-peg**  
(for 8 aprons)  
**#ARP-16**  
\$398.00



**20-peg**  
(for 10 aprons)  
**#ARP-20**  
\$461.00








**FREE** Protective Apparel Tracking and Inventory System!

## Exclusive from Techno-Aide!

All new Techno-Aide garments include a scanable QR code for automatic tracking. Upload garment data into your unique ApronHub portal for permanent inventory management. Receive alerts for inspection dates through a customizable format. Retrofit QR code tags onto your existing Techno-Aide or non-Techno-Aide garments. Track inspection history, location, user, etc., through a customizable reports menu. 2D Scanners are available through Techno-Aide or you can use your own 2D-scanner. ApronHub has everything you need to quickly and easily meet the highest QA standards. For more information visit [www.apronhub.com](http://www.apronhub.com) today!





## Accessories

### Tools for Scanning and Labeling Protective Apparel

Stay one step ahead of the inspectors and enhance your QA! All new garments manufactured by Techno-Aide include a scanable QR Code Serial Number label for automated or manual tracking. However, we realize that you may have older, unlabeled garments or garments from another manufacturer. To update all your protective garments, we have designed this retrofit labeling kit so that your garments can be tracked utilizing the most up-to-date tracking technology available.

CAT #	Description	Price
AHSCB	2D Scanner with Cable	call
AHSXW	2D Scanner Wireless	call
AHETL	Eyelet Tool	call
AHES50	Silver Eyelets (50)	\$25.00
AHQRC	QR Code Retrofit Labels (100)	\$68.00
AHDT19	Date Tags: 2019 (50)	\$ 14.00
AHDT20	Date Tags: 2020 (50)	\$ 14.00

**#AHETL**  
call for pricing



**#AHQRC**  
\$68.00



**#AHSCB/AHSXW**  
call for pricing



**#AHES50**  
\$25.00



**#AHDT19/AHDT20**  
\$14.00








ADMINISTRATION/MANUFACTURING

7117 Centennial Blvd. Nashville, TN 37209

Orders can be placed online or sent to [orders@techno-aide.com](mailto:orders@techno-aide.com)

For official Quotes: [quotes@techno-aide.com](mailto:quotes@techno-aide.com)

Any questions can be directed to [customerservice@techno-aide.com](mailto:customerservice@techno-aide.com)

P: 800.251.2629 or 615.350.7030

F: 800.554.6275 or 615.350.7879

[techno-aide.com](http://techno-aide.com)