




Thyroid Collars

2020
PRODUCT
GUIDE

Thyro-Guards



Standard

Binding all around means affordable comfort and flexibility.



No Binding

No binding means extra comfort, flexibility and lighter weight.



Mid-Binding

No binding at the top means extra comfort under the chin.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TYR	Regular	\$44.00
TYL	Light Wt.	\$54.00
TYG	Lead-Free	\$64.00

CAT. #	Lead Type	Price
TNR	Regular	\$74.00
TNL	Light Wt.	\$84.00
TNG	Lead-Free	\$94.00

CAT. #	Lead Type	Price
TVR	Regular	\$59.00
TVL	Light Wt.	\$69.00
TVG	Lead-Free	\$79.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TYSR	Regular	\$34.00
TYSL	Light Wt.	\$39.00
TYSG	Lead-Free	\$48.00

CAT. #	Lead Type	Price
TNSR	Regular	\$53.00
TNSL	Light Wt.	\$57.00
TNSG	Lead-Free	\$66.00

CAT. #	Lead Type	Price
TVSR	Regular	\$44.00
TVSL	Light Wt.	\$49.00
TVSG	Lead-Free	\$58.00

Attached Standard

Attached No Binding

Attached Mid-Binding

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TAR	Regular	\$54.00
TAL	Light Wt.	\$64.00
TAG	Lead-Free	\$74.00

CAT. #	Lead Type	Price
TCR	Regular	\$84.00
TCL	Light Wt.	\$94.00
TCG	Lead-Free	\$104.00

CAT. #	Lead Type	Price
TWR	Regular	\$69.00
TWL	Light Wt.	\$79.00
TWG	Lead-Free	\$89.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TASR	Regular	\$46.00
TASL	Light Wt.	\$51.00
TASG	Lead-Free	\$59.00

CAT. #	Lead Type	Price
TCSR	Regular	\$63.00
TCSL	Light Wt.	\$67.00
TCSG	Lead-Free	\$76.00

CAT. #	Lead Type	Price
TWSR	Regular	\$54.00
TWSL	Light Wt.	\$59.00
TWSG	Lead-Free	\$68.00

Attached Thyroid Collars Info

All Attached Thyroid Collars include:

- 1) A snap ring for easy collar attachment.
- 2) Attach anytime with tethered strap.
- 3) Unique 2-in-1 pocket design.

[store collar in one pocket and use other pocket for accessories]

Thyroid collars are 6" tall x 6 1/2" wide, with a 24" strap.



Attached Thyroid option and our 2-in-1 pocket design!

Buckle Closure Add-On

If you prefer buckle closure to regular hook & loop, add "-B" to any of the collar numbers above and add \$5.00 to the price.

(Ex.: TYR-B would be \$49.00)



Buckle closure add-on available!

Specialty Thyro-Guards



Pediatric

Binding all around means affordable comfort and flexibility.



Mammography

No binding means extra comfort, flexibility and lighter weight.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TPR	Regular	\$55.00
TPL	Light Wt.	\$60.00
TPG	Lead-Free	\$73.00

CAT. #	Lead Type	Price
TMR	Regular	\$69.00
TML	Light Wt.	\$79.00
TMG	Lead-Free	\$89.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TPSR	Regular	\$53.00
TPSL	Light Wt.	\$57.00
TPSG	Lead-Free	\$66.00

CAT. #	Lead Type	Price
TMSR	Regular	\$53.00
TMSL	Light Wt.	\$57.00
TMSG	Lead-Free	\$66.00

Attached Thyroid Collars Info

All Attached Thyroid Collars include:

- 1) A snap ring for easy collar attachment.
- 2) Attach anytime with tethered strap.
- 3) Unique 2-in-1 pocket design.

[store collar in one pocket and use other pocket for accessories]

Thyroid collars are 6" tall x 6 1/2" wide, with a 24" strap.



Attached Thyroid option and our 2-in-1 pocket design!

Buckle Closure Add-On

If you prefer buckle closure to regular hook & loop, add "-B" to any of the collar numbers above and add \$5.00 to the price.

(Ex.: TYR-B would be \$49.00)




Buckle closure add-on available!


Protective Apparel Fabric Options

NYLON FABRICS

FEATHER RIPSTOP FABRIC Light Weight Nylon, Filament Reinforced Fabric



DIAMOND RIPSTOP FABRIC Medium Weight Nylon, Filament Reinforced Fabric




NYLON SOLIDS Liquid Repellent, Light & Supple



SPECIALTY FABRICS

COLORFUL PRINTS Liquid Repellent, Medium Weight



NEW 69-Skulls




70-Blue Plaid



71-Leopard



72-Lots•0•Dots



73-Sleepy Sheep



74-Playful Pups



83-Cool Waves



77-Life Vines



78-Waterfall



80-Red Paws



81-Black Paws



76-Blue Paws



82-Pink Camo




85-Digital Camo



86-Camo Night

Prints are shown
SMALLER than
actual size!

TEXTURED FABRIC Textured High Tech Fabrics



15-Diamondback



16-Crocogator



18-Bronze Tech

Premium Fabric Add-On


NEW Coolwear Fabric

Designed for optimum cooling comfort!

Our new Coolwear fabric add-on provides premium cooling comfort all day long! **This premium fabric can now be added to ANY garment!** The polyester mesh lining features air vented space and stretch-in-width for extreme comfort and fit. This innovative add-on is designed for ultimate performance and comfort.

CP25 - \$25 add-on for all thyroid collars

Features premium comfort with all day cooling!




Protective Apparel Add-ONS

FREE Embroidery Sewn-In


Select products come with the first line (up to 24 characters) of embroidery **FREE!** Add flexibility with an embroidered patch. Choose from a style below. Add flexibility with an embroidery patch! **Embroidered items are non-returnable & non-refundable.**

Available Styles



#EML2 | \$20 for **second embroidery line**, up to 24 letters.

Thyroid Embroidery Example



Embroidery Colors

- White* (default)
- Black** (default)
- Lt. Gray*
- Charcoal*
- Royal Blue*
- Red*
- Yellow*
- Orange*
- Tan*
- Brown*

Binding Colors

- Black** (default)
- Light Gray** (default)
- Navy**
- Royal Blue**
- Red**
- Yellow**
- Orange**
- Tan**
- Pink**
- Purple**
- Green**
- Burgundy**

Protective Media Guide


MEDIA	BENEFITS	PROTECTION	ATTENUATION	WEIGHT
Regular Lead	Most economical and recyclable	100% Lead	0.5mm Pb Equiv. Protection avg. 97% attenuation over 80, 100, 120 KVP	9 - lbs. (based on medium front protection apron)
Lightweight Lead	Lighter than regular lead and recyclable	Lead-Tin	0.5mm Pb Equiv. Protection avg. 96.8% attenuation over 80, 100, 120 KVP	8 - lbs. (based on medium front protection apron)
Lead Free	Lightest Weight, Lead-Free, Environmentally Friendly	Tungsten-Antimony	0.5mm Pb Equiv. Protection avg. 96.6% attenuation over 80, 100, 120 KVP	7 - lbs. (based on medium front protection apron)



FREE Protective Apparel Tracking and Inventory System!

Exclusive from Techno-Aide!

All new Techno-Aide garments include a scanable QR code for automatic tracking. Upload garment data into your unique ApronHub portal for permanent inventory management. Receive alerts for inspection dates through a customizable format. Retrofit QR code tags onto your existing Techno-Aide or non-Techno-Aide garments. Track inspection history, location, user, etc., through a customizable reports menu. 2D Scanners are available through Techno-Aide or you can use your own 2D-scanner. ApronHub has everything you need to quickly and easily meet the highest QA standards. For more information visit www.apronhub.com today!







Accessories

Tools for Scanning and Labeling Protective Apparel

Stay one step ahead of the inspectors and enhance your QA! All new garments manufactured by Techno-Aide include a scannable QR Code Serial Number label for automated or manual tracking. However, we realize that you may have older, unlabeled garments or garments from another manufacturer. To update all your protective garments, we have designed this retrofit labeling kit so that your garments can be tracked utilizing the most up-to-date tracking technology available.

CAT #	Description	Price
AHSCB	2D Scanner with Cable	call
AHSXW	2D Scanner Wireless	call
AHETL	Eyelet Tool	call
AHES50	Silver Eyelets (50)	\$25.00
AHQRC	QR Code Retrofit Labels (100)	\$68.00
AHDT19	Date Tags: 2019 (50)	\$ 14.00
AHDT20	Date Tags: 2020 (50)	\$ 14.00





ADMINISTRATION/MANUFACTURING

7117 Centennial Blvd. Nashville, TN 37209

Orders can be placed online or sent to orders@techno-aide.com

For official Quotes: quotes@techno-aide.com

Any questions can be directed to customerservice@techno-aide.com

P: 800.251.2629 or 615.350.7030

F: 800.554.6275 or 615.350.7879

techno-aide.com