



Protective Apparel

2020
PRODUCT
GUIDE

CHECK OUT WHAT'S NEW FOR 2020

Techno-Aide is ready to kick 2020 off with some amazing new products! Here are a few key highlights you can't afford to miss!

New Wide Belt Vest Guard Plus & Wide Belt Skirt & Kilt | PG.11

These new garments feature a 6" wide adjustable elastic belt for premium back support and comfort!

New Quick Release Wide Belt Apron | PG.8

This new revolutionary apron features a quick-release shoulder strap to easily slide off under sterile gowns!

New EZ Ship Apron Program | PG.4

Protection shouldn't be complicated, it should be EZ! Keep it simple and get it quick! Techno-Aide is proud to introduce the new EZ Ship Apron program. No fuss, no mess, no wait! Order one of our EZ Ship aprons and get it in 12 days!

Updated Quick Ship Program | PG.6

Need it fast? Techno-Aide's revolutionary Quick Ship program is back and better than ever. With our new default Anti-Microbial Sapphire fabric, these ready to order products ship on the next business day!

New Coolwear Fabric Option | PG.30

Our new Coolwear fabric add-on provides optimum cooling comfort! This premium fabric can now be added on to ANY garment! The polyester mesh lining features air vented space and stretch-in-width for extreme comfort and fit.

Custom Product Pricing Update | PG.2

Check out page two for an important update to our custom apparel process. **Avoid the headache and order a standard size!** Techno-Aide sizing charts are based on universal body sizes. Use our Best Fit Sizing Guide to find the perfect fit!

Techno-Aide Best Fit Sizing Guide

Techno-Aide Sizing Charts are based on universal body sizes. But let's face it... nobody is average! A great apron won't fit like a pair of jeans or a suit. **Aprons are meant to fit multiple body types and sizes!** So, in order to get the "Best Fit for You," please consider the following when measuring yourself for a standard garment:

PROTECTION IS PRIMARY: Why are you wearing an apron? Protection from radiation! At Techno-Aide, we understand that comfort is important, but never as important as protection. Our aprons are designed to keep the user safe above all else. But who says you can't add a little style and comfort along the way?

CIRCUMFERENCE CONSIDERATIONS: When measuring for a Front Protection Apron, measure from one side of the widest part of your body, around to the other side - the actual distance you want protected. Then, tell us if that measurement is your chest, tummy, or hips.

HEIGHT CONSIDERATIONS: In relation to your circumference, do you have a long or short torso? We want to make sure that you are adequately protected!

Front Protection Apron			
Size	Front Width	Length	Circumference
Petite	20"	32"	< 36"
Small	21"	34"	< 40"
Medium	22"	36"	< 44"
Large	24"	38"	< 48"
X-Large	26"	40"	< 52"
2X	28"	40"	< 56"
3X*	30"	40"	< 60"

Vest (Male) & Full Wrap-Around			
Size	Chest/Tummy	Vest Length	Full Wrap Length
Petite	35"-36"	19.5"	32"
Small	37"-38"	20"	34"
Medium	39"-41"	21"	36"
Large	42"-44"	22"	38"
X-Large	45"-49"	23"	40"
2X	50"-54"	24"	40"
3X*	55"-59"	24"	40"

Vest (Female)				
Size	Bust	Waist	Hip	Length
Petite	33"-34"	< 31"	35"-36"	19"
Small	35"-36"	< 33"	37"-38"	19.5"
Medium	37"-38"	< 35"	39"-40"	20"
Large	39"-41"	< 39"	41"-43"	21"
X-Large	42"-44"	< 41"	44"-46"	22"
2X	45"-49"	< 46"	47"-51"	23"
3X*	50"-54"	< 51"	52"-56"	23"

Skirt and Kilt			
Size	Waist	Hip	Length
Petite	25"-27"	32"-34"	19"
Small	28"-29"	35"-36"	19"
Medium	30"-33"	37"-40"	20"
Large	34"-37"	41"-44"	21"
X-Large	38"-41"	45"-48"	22"
2X	42"-45"	49"-52"	22"
3X	46"-50"	53"-57"	22"



Notes About Measuring!

Sizes shown for apparel (other than front protection) are **ACTUAL** body measurements.

All lengths are from the **sternal notch or waist** downward, depending on the garment.

One Size to Fit Them All!

Don't get hung up with complicated customization!

Standard sizes are designed to fit nearly every body type. **As shown to the right**, our standard Large Flex-Guard can be worn properly by three people with completely different body types!

Use our expert guides to find the perfect size today!



*Each person is wearing a standard (non-custom) large Flex-Guard apron appropriately. Photos provided to show standard sizing fit. Apron shown in default 61-Sapphire material.

EZ-Guard (Hook & Loop)

EZ-On / EZ-Off. Our Most Popular Front Protection apron with an easy criss-cross back and **Hook and Loop** closure on the front. Perfect for every day imaging protection. Double-fold binding for strength and tear resistance. Concealed shoulder pads for extra comfort. Free dual pocket to store an attached thyroid collar and a separate pocket for other items.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
EVR	Regular	\$229.00	\$250.00	\$289.00
EVL	Light Wt.	\$289.00	\$314.00	\$368.00
EVG	Lead-Free	\$374.00	\$435.00	\$474.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
EVSR	Regular	\$203.00	\$219.00	\$251.00
EVSL	Light Wt.	\$225.00	\$279.00	\$325.00
EVSG	Lead-Free	\$310.00	\$384.00	\$453.00

Our Most Popular Apron!



See our "Best-Fit" sizing guide on PG 2



EZ-Guard (Buckle)

EZ-On / EZ-Off. Our Most Popular Front Protection apron with an easy criss-cross back and **Buckle** closure on the front. Perfect for those who prefer a quick adjustable closure. Double-fold binding for strength and tear resistance. Concealed shoulder pads for extra comfort. Free dual pocket to store an attached thyroid collar and a separate pocket for other items.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
EBR	Regular	\$229.00	\$250.00	\$289.00
EBL	Light Wt.	\$289.00	\$314.00	\$368.00
EBG	Lead-Free	\$374.00	\$399.00	\$474.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
EBSR	Regular	\$203.00	\$219.00	\$251.00
EBSL	Light Wt.	\$225.00	\$279.00	\$325.00
EBSG	Lead-Free	\$310.00	\$384.00	\$453.00

TechnoAide
PRESENTS

NEW 12 DAY SHIPPING!

Protection shouldn't be complicated, it should be **EZ!** Simply choose from three **sizes** and three **media types**.

Keep it simple, get it quick!
We ship in **12 days!***

- Our most popular Apron
- Our most popular Sizes
- Our most popular Fabrics
- All media types available!
- 0.5mm Pb Equivalency

EZSHIP» APRONS



REGULAR

LIGHT WEIGHT

LEAD FREE

EZ Guard Medium (65)	EVRM	\$229	EVLM	\$289	EVGM	\$374
EZ Guard Large (61)	EVRL	\$250	EVL	\$314	EVGL	\$435
EZ Guard X-Large (64)	EVRX	\$250	EVLX	\$314	EVGX	\$435

**DON'T FORGET TO ADD A MATCHING THYRO-GUARD!
SHIPS AT THE SAME TIME!**

Standard Thyro-Guard	TYR	\$44	TYL	\$54	TYG	\$64
Attached Thyro-Guard	TAR	\$54	TAL	\$64	TAG	\$74

*Shipping dates are provided from date of order confirmation. "EZ Ship" is a limited time promotion. Techno-Aide reserves the right to cancel, change or discontinue any promotion at any time for any reason. Selections are limited to the provided/advertised CAT/SKUs, combinations, and products only. Quantities greater than 10 will require a manual quote and may not ship in advertised time. Please contact customer service at CustomerService@Techno-Aide.com for more information. Visit Techno-Aide.com for more product information.

Flex-Guard (Hook & Loop)

Front Protection apron with flexible stretch criss-cross back and **Hook and Loop** closure on the front. Excellent fit for full day imaging protection by providing less weight on the shoulders and snug fit at the back. Double-fold binding for strength and tear resistance. Free dual pocket to store an attached thyroid collar and a separate pocket for other items. Concealed shoulder pads for extra comfort.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
FVR	Regular	\$238.00	\$280.00	\$319.00
FVL	Light Wt.	\$300.00	\$344.00	\$383.00
FVG	Lead-Free	\$395.00	\$438.00	\$503.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
FVSR	Regular	\$211.00	\$249.00	\$281.00
FVSL	Light Wt.	\$260.00	\$309.00	\$355.00
FVSG	Lead-Free	\$333.00	\$376.00	\$416.00



See our "Best-Fit" sizing guide on PG 2



Flex-Guard (Buckle)

Front Protection apron with flexible stretch criss-cross back and **Buckle** closure on the front. Excellent fit for full day imaging protection by providing less weight on the shoulders and snug fit at the back. Double-fold binding for strength and tear resistance. Free dual pocket to store an attached thyroid collar and a separate pocket for other items. Concealed shoulder pads for extra comfort.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
FBR	Regular	\$238.00	\$280.00	\$319.00
FBL	Light Wt.	\$300.00	\$344.00	\$383.00
FBG	Lead-Free	\$395.00	\$438.00	\$503.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
FBSR	Regular	\$211.00	\$249.00	\$281.00
FBSL	Light Wt.	\$268.00	\$309.00	\$355.00
FBSG	Lead-Free	\$333.00	\$376.00	\$416.00




NEED IT FAST?

Keep it simple, get it quick!

Techno-Aide is proud to offer the fastest solution in protection! Our Quick Ship program keeps you protected when you need it most. Fast, easy and quality protection you can count on.

All Quick Ship Items ship by the next business day!*

Flex Guard - Small, Lightweight lead, in Sapphire (61)	QS-FVLS-61	\$300
Flex Guard - Medium, Lightweight lead, in Sapphire (61)	QS-FVLM-61	\$300
Flex Guard - Large, Lightweight lead, in Sapphire (61)	QS-FVLL-61	\$344
Flex Guard - X-Large, Lightweight lead, in Sapphire (61)	QS-FVLX-61	\$344
Gonadal Guard Set (S,M,L) - Regular lead, in Sapphire (61)	QS-GSR-61	\$120
Lap Guard Half Apron - Regular lead, in Sapphire (61)	QS-LGRL-61	\$110
Thyroid Collar - Regular lead, in Sapphire (61)	QS-TYR-61	\$44



A \$15 service fee applies to each quick ship item. Part number QKS-15 must be added for each item. Limit of 2 items per day, per order.

*Shipping dates are provided from date of order confirmation. Techno-Aide reserves the right to cancel, change or discontinue any promotion at any time for any reason. Selections are limited to the provided/advertised CAT/SKUs, combinations, and products only. Contact customer service for quantities greater than two (2). A \$15 fee applies to all quick ship items on any order regardless of quantity or type. Please contact customer service at CustomerService@Techno-Aide.com for more information. Visit Techno-Aide.com for more product information.

Surgery-Guard

This surgical drop off apron easily slides off under a sterile gown with no arm holes. No worry about compromising sterile fields in the OR. Crossover shoulder to **Hook and Loop** closure for quick release at the waist. Wide shoulders with concealed shoulder pads to optimize weight distribution and provide extra comfort. Double-fold binding for strength and tear resistance. Free dual pocket to store an attached thyroid collar and a separate pocket for other items.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
SRR	Regular	\$229.00	\$260.00	\$300.00
SRL	Light Wt.	\$290.00	\$325.00	\$379.00
SRG	Lead-Free	\$384.00	\$408.00	\$485.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
SRSR	Regular	\$192.00	\$229.00	\$262.00
SRS�	Light Wt.	\$236.00	\$289.00	\$305.00
SRSĠ	Lead-Free	\$327.00	\$395.00	\$415.00



See our "Best-Fit" sizing guide on PG 2



OR-Guard

Flexible-back apron that slides off easily under sterile gowns with break away shoulders and **Hook and Loop** closure at the hips. Double-fold binding for strength and tear resistance. Free dual pocket to store an attached thyroid collar and a separate pocket for other items. Concealed shoulder pads for extra comfort.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
ORR	Regular	\$270.00	\$291.00	\$324.00
ORL	Light Wt.	\$329.00	\$349.00	\$403.00
ORG	Lead-Free	\$412.00	\$472.00	\$520.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
ORSR	Regular	\$229.00	\$254.00	\$286.00
ORSL	Light Wt.	\$325.00	\$357.00	\$360.00
ORSG	Lead-Free	\$393.00	\$419.00	\$469.00

Ortho-Guard

Built-in Back-Guard for Maximum Support. Tethered belt can be adjusted vertically and laterally for effective back support. Back strap assures best fit and comfort. 0.5mm Pb equiv. front protection. Free dual pocket and shoulder pads.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
OGR	Regular	\$270.00	\$289.00	\$324.00
OGL	Light Wt.	\$329.00	\$349.00	\$403.00
OGG	Lead-Free	\$409.00	\$467.00	\$524.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
OGSR	Regular	\$229.00	\$254.00	\$286.00
OGSL	Light Wt.	\$325.00	\$351.00	\$360.00
OGSG	Lead-Free	\$398.00	\$419.00	\$488.00



Shoulder strap included!

See our "Best-Fit" sizing guide on PG 2



Quick release shoulder strap makes it easy to slide off under sterile gowns!

NEW Quick Release Wide Belt

Techno-Aide proudly presents our **NEW** Quick-Release Wide Belt Apron. Featuring a quick-release shoulder strap to easily slide off under sterile gowns. Upgrade your OR Guard to the Quick Release Wide Belt Guard with our added attached wide belt. Our 6" elastic supported belt guard helps reduce back stress in the Operating Room or during longer imaging procedures in surgery centers. 1" buckle closure for snug fit.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
OWR	Regular	\$320.00	\$341.00	\$374.00
OWL	Light Wt.	\$379.00	\$399.00	\$453.00
OWG	Lead-Free	\$462.00	\$522.00	\$570.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
OWSR	Regular	\$279.00	\$304.00	\$336.00
OWSL	Light Wt.	\$375.00	\$407.00	\$410.00
OWSG	Lead-Free	\$443.00	\$469.00	\$519.00

Tie-Guard Apron

Excellent protection and value prices. Comes with a 1" buckle closure that is long enough to tie. 0.5mm Pb equivalency protection. Double fold binding for strength and tear resistance. **Free dual pocket** to store an attached thyroid collar and a separate pocket for other items. Concealed shoulder pads or extra comfort.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
EPR	Regular	\$229.00	\$250.00	\$289.00
EPL	Light Wt.	\$289.00	\$314.00	\$368.00
EPG	Lead-Free	\$384.00	\$405.00	\$474.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
EPSR	Regular	\$193.00	\$219.00	\$251.00
EPSL	Light Wt.	\$225.00	\$279.00	\$325.00
EPSP	Lead-Free	\$303.00	\$384.00	\$412.00



See our "Best-Fit" sizing guide on PG 2



Econo-Guard

Excellent Protection and Value Priced. Comes with 1" buckle closure that is long enough to tie. Double-fold binding for strength and tear resistance. **There is no pocket for this garment. Not compatible with attached thyroid collar.** If pocket and attached collar are required, please refer to our EZ guard aprons.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
EGR	Regular	\$188.00	\$220.00	\$270.00
EGL	Light Wt.	\$256.00	\$288.00	\$348.00
EGG	Lead-Free	\$334.00	\$394.00	\$469.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
EGSR	Regular	\$175.00	\$194.00	\$223.00
EGSL	Light Wt.	\$220.00	\$248.00	\$290.00
EGSG	Lead-Free	\$304.00	\$321.00	\$373.00

Max-Wrap

Wrap around coat style garment offers 0.5mm Pb equiv. protection in front and 0.25mm in back. Solid panels provide no overlap meaning less bulk, better fit, better freedom of movement and secure closure. Double-fold binding for strength and tear resistance. Breakaway left shoulder and flexible back is fast and easy. Secure closure at the waist and shoulder for optimal fit. Flexible back support panels. Free dual pocket to store an attached thyroid collar and a separate pocket for other items. Concealed shoulder pads for extra comfort. **Hook and Loop** closure with flexible back is fast and easy to put on. Excellent protection for full coverage.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
MWL	Light Wt.	\$437.00	\$486.00	\$541.00
MWG	Lead-Free	\$542.00	\$610.00	\$649.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
MWSL	Light Wt.	\$358.00	\$385.00	\$427.00
MWSG	Lead-Free	\$464.00	\$494.00	\$551.00



See our "Best-Fit" sizing guide on PG 2



Max-Guard Plus

Overlap wrap around coat style garment offers 0.5mm Pb equiv. protection in front and 0.25mm in back. Distributes weight between shoulders and hips. **Hook and Loop** closure with inside elastic Back Guard for excellent back support. Double-fold binding for strength and tear resistance. Free dual pocket to store an attached thyroid collar and a separate pocket for other items. Concealed shoulder pads for extra comfort. Excellent protection for full coverage.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
MXL	Light Wt.	\$420.00	\$431.00	\$517.00
MXG	Lead-Free	\$525.00	\$560.00	\$620.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
MXSL	Light Wt.	\$367.00	\$431.00	\$517.00
MXSG	Lead-Free	\$525.00	\$560.00	\$620.00

NEW Male Wide Belt Vest Guard Plus

Full front overlap design distributes weight between shoulders and hips and provides full coverage. Provides .5mm Pb equiv protection in the front and .25mm Pb equiv. protection in the back. New inside elastic back guard and front belt assures a snug fit with maximum back support. Double-fold binding offers strength and tear resistance. Hook and Loop closure is fast and easy. Free dual pocket to store an attached thyroid collar and a separate pocket for other items. Concealed shoulder pads for extra comfort. 1" safety buckle at shoulder. **Male vest closes from right to left.**

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VPR	Regular	\$260.00	\$300.00	\$349.00
VPL	Light Wt.	\$320.00	\$350.00	\$414.00
VPG	Lead-Free	\$399.00	\$416.00	\$490.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VPSR	Regular	\$242.00	\$266.00	\$313.00
VPSL	Light Wt.	\$276.00	\$307.00	\$370.00
VPSG	Lead-Free	\$318.00	\$345.00	\$413.00



See our "Best-Fit" sizing guide on PG 2



NEW Wide Belt Kilt

Techno-Aide has added a 6" wide adjustable elastic belt to our standard vest and skirt/kilt for additional back support and comfort. Weight is distributed evenly around hips. Free pocket and waist pad. Double-fold binding for strength and tear resistance. **Wrap around protection: 0.5mm Pb Eq front, 0.25mm back.** **Male kilt closes from right to left.**

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
KWR	Regular	\$244.00	\$264.00	\$284.00
KWL	Light Wt.	\$304.00	\$324.00	\$344.00
KWG	Lead-Free	\$354.00	\$374.00	\$394.00

NEW Female Wide Belt Vest Guard Plus

Full front overlap design distributes weight between shoulders and hips and provides full coverage. Provides .5mm Pb equiv protection in the front and .25mm Pb equiv. protection in the back. New inside elastic back guard and front belt assures a snug fit with maximum back support. Double-fold binding offers strength and tear resistance. Hook and Loop closure is fast and easy. Free dual pocket to store an attached thyroid collar and a separate pocket for other items. Concealed shoulder pads for extra comfort. 1" safety buckle at shoulder. **Female vest closes from left to right.**

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VWR	Regular	\$260.00	\$300.00	\$349.00
VWL	Light Wt.	\$320.00	\$350.00	\$414.00
VWG	Lead-Free	\$399.00	\$416.00	\$490.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VWSR	Regular	\$242.00	\$266.00	\$313.00
VWSL	Light Wt.	\$276.00	\$307.00	\$370.00
VWSG	Lead-Free	\$318.00	\$345.00	\$413.00



See our "Best-Fit" sizing guide on PG 2



New 6" wide, adjustable elastic belt for additional back support!

NEW Wide Belt Skirt

Techno-Aide has added a 6" wide adjustable elastic belt to our standard vest and skirt/kilt for additional back support and comfort. Weight is distributed evenly around hips. Free pocket and waist pad. Double-fold binding for strength and tear resistance. **Wrap around protection: 0.5mm Pb Eq front, 0.25mm back.** **Female vest closes from left to right.**

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
SWR	Regular	\$244.00	\$264.00	\$284.00
SWL	Light Wt.	\$304.00	\$324.00	\$344.00
SWG	Lead-Free	\$354.00	\$374.00	\$394.00

Vest-Guard Male

Full front overlap design saves weight. Distributes weight between shoulders and hips. 0.5mm Pb equivalency protection front, 0.25mm back. **Hook and Loop** closure is fast and easy. Free dual pocket and shoulder pads. Double-fold binding for strength and tear resistance. Concealed shoulder pads for extra comfort. Free dual pocket to store an attached thyroid collar and a separate pocket for other items.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VMR	Regular	\$210.00	\$250.00	\$299.00
VML	Light Wt.	\$270.00	\$300.00	\$364.00
VMG	Lead-Free	\$349.00	\$366.00	\$440.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VMSR	Regular	\$192.00	\$216.00	\$263.00
VMSL	Light Wt.	\$226.00	\$257.00	\$320.00
VMSG	Lead-Free	\$268.00	\$295.00	\$363.00

(Vests do not include the two 1" safety buckles. If you would like these as a FREE Add-On, please specify at the time of ordering.)



See our "Best-Fit" sizing guide on PG 2

Kilt-Guard Male

All Day Comfort for the Lower Body

Designed for best comfort.. Weight is distributed evenly around hips. **Wrap around protection: 0.5mm Pb Eq front, 0.25mm back.** **Hook and Loop** closure with 2" safety buckle. Free pocket and waist pad. Double-fold binding for strength and tear resistance.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
KMR	Regular	\$210.00	\$236.00	\$255.00
KML	Light Wt.	\$270.00	\$293.00	\$319.00
KMG	Lead-Free	\$340.00	\$368.00	\$395.00



Vest-Guard Female

Full front overlap design saves weight. Distributes weight between shoulders and hips. 0.5mm Pb equivalency protection front, 0.25mm back. **Hook and Loop** closure is fast and easy. Free dual pocket and shoulder pads. Double-fold binding for strength and tear resistance. Concealed shoulder pads for extra comfort. Free dual pocket to store an attached thyroid collar and a separate pocket for other items.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VFR	Regular	\$210.00	\$250.00	\$299.00
VFL	Light Wt.	\$270.00	\$300.00	\$364.00
VFG	Lead-Free	\$349.00	\$366.00	\$440.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VFSR	Regular	\$192.00	\$216.00	\$263.00
VFSL	Light Wt.	\$226.00	\$257.00	\$320.00
VFSG	Lead-Free	\$268.00	\$295.00	\$363.00



(Vests do not include the two 1" safety buckles. If you would like these as a FREE Add-On, please specify at the time of ordering.)

See our "Best-Fit" sizing guide on PG 2

Skirt-Guard

All Day Comfort for the Lower Body

Designed for best comfort.. Weight is distributed evenly around hips. **Wrap around protection: 0.5mm Pb Eq front, 0.25mm back.** **Hook and Loop** closure with 2" safety buckle. Free pocket and waist pad. Double-fold binding for strength and tear resistance.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
SFR	Regular	\$210.00	\$236.00	\$255.00
SFL	Light Wt.	\$270.00	\$293.00	\$319.00
SFG	Lead-Free	\$344.00	\$354.00	\$364.00



Features an easy to use 2" safety buckle!

Vest-Wrap | Less Bulk, Better Fit

Separate Men and Women’s styles for **the best comfort and the ultimate freedom of movement!** Distributes weight between shoulders and hips. 0.25mm Pb equivalency protection on the back. **Hook and Loop** closure with flexible back. Generous armholes for EZ-fit. Concealed shoulder pads for extra comfort. Double-fold binding for strength and tear resistance. Dual pocket system stores an attached thyroid collar and other items such as markers in a separate slot.

Male Vest-Wrap

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VBR	Regular	\$215.00	\$263.00	\$307.00
VBL	Light Wt.	\$269.00	\$326.00	\$359.00
VBG	Lead-Free	\$348.00	\$400.00	\$446.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VBSR	Regular	\$192.00	\$216.00	\$263.00
VBSL	Light Wt.	\$226.00	\$257.00	\$320.00
VBSG	Lead-Free	\$268.00	\$295.00	\$363.00

Female Vest-Wrap

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VGR	Regular	\$215.00	\$263.00	\$307.00
VGL	Light Wt.	\$269.00	\$326.00	\$359.00
VGG	Lead-Free	\$348.00	\$400.00	\$446.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Petite/ Sm/ Med	LG / XL	2X / 3X
VGSR	Regular	\$192.00	\$216.00	\$263.00
VGSL	Light Wt.	\$226.00	\$257.00	\$320.00
VGGG	Lead-Free	\$268.00	\$295.00	\$363.00



Care-Guard Pregnancy Apron | Special Protection for Fetal Care

Front Details: Front protection with flexible back straps. 1.0mm Pb equivalency over fetal area, 0.5mm over remaining area. 22" x 24" pregnancy panel can be removed at completion of pregnancy to create a standard Flex-Guard Apron with 0.5mm Pb equivalency protection. **Hook and Loop** closure is fast and easy. Concealed shoulder pads for extra comfort. Double-fold binding for strength and tear resistance.

Back Details: Works with the front apron for full protection. 0.25mm Pb equivalency. 1" buckle closure (long enough to tie) in front.

Care-Guard Front Only

CAT. #	Lead Type	Regular	Tall/Full
CGL	Light Wt.	\$474.00	\$524.00
CGG	Lead-Free	\$524.00	\$572.00

Size	Dimensions
Reg	24" x 36"
Tall	24" x 40"
Full	28" x 38"


Care-Guard Back Only

CAT. #	Lead Type	Regular	Tall/Full
CBL	Light Wt.	\$288.00	\$348.00
CBG	Lead-Free	\$394.00	\$469.00

Attached Thyroid Collar Available. See Thyroid collars on pg. 16.



Thyro-Guards



Standard

Binding all around means affordable comfort and flexibility.



No Binding

No binding means extra comfort, flexibility and lighter weight.



Mid-Binding

No binding at the top means extra comfort under the chin.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TYR	Regular	\$44.00
TYL	Light Wt.	\$54.00
TYG	Lead-Free	\$64.00

CAT. #	Lead Type	Price
TNR	Regular	\$74.00
TNL	Light Wt.	\$84.00
TNG	Lead-Free	\$94.00

CAT. #	Lead Type	Price
TVR	Regular	\$59.00
TVL	Light Wt.	\$69.00
TVG	Lead-Free	\$79.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TYSR	Regular	\$34.00
TYSL	Light Wt.	\$39.00
TYSG	Lead-Free	\$48.00

CAT. #	Lead Type	Price
TNSR	Regular	\$53.00
TNSL	Light Wt.	\$57.00
TNSG	Lead-Free	\$66.00

CAT. #	Lead Type	Price
TVSR	Regular	\$44.00
TVSL	Light Wt.	\$49.00
TVSG	Lead-Free	\$58.00

Attached Standard

Attached No Binding

Attached Mid-Binding

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TAR	Regular	\$54.00
TAL	Light Wt.	\$64.00
TAG	Lead-Free	\$74.00

CAT. #	Lead Type	Price
TCR	Regular	\$84.00
TCL	Light Wt.	\$94.00
TCG	Lead-Free	\$104.00

CAT. #	Lead Type	Price
TWR	Regular	\$69.00
TWL	Light Wt.	\$79.00
TWG	Lead-Free	\$89.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TASR	Regular	\$46.00
TASL	Light Wt.	\$51.00
TASG	Lead-Free	\$59.00

CAT. #	Lead Type	Price
TCSR	Regular	\$63.00
TCSL	Light Wt.	\$67.00
TCSG	Lead-Free	\$76.00

CAT. #	Lead Type	Price
TWSR	Regular	\$54.00
TWSL	Light Wt.	\$59.00
TWSG	Lead-Free	\$68.00

Attached Thyroid Collars Info

All Attached Thyroid Collars include:

- 1) A snap ring for easy collar attachment.
- 2) Attach anytime with tethered strap.
- 3) Unique 2-in-1 pocket design.

[store collar in one pocket and use other pocket for accessories]

Thyroid collars are 6" tall x 6.5" wide, with a 24" strap.



Attached Thyroid option and our 2-in-1 pocket design!

Buckle Closure Add-On

If you prefer buckle closure to regular hook & loop, add "-B" to any of the collar numbers above and add \$5.00 to the price.

(Ex.: TYR-B would be \$49.00)



Buckle closure add-on available!

Specialty Thyro-Guards



Pediatric

Binding all around means affordable comfort and flexibility.



Default color:
Berry with Pink ribbon only

Mammography

No binding means extra comfort, flexibility and lighter weight.

0.5mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TPR	Regular	\$55.00
TPL	Light Wt.	\$60.00
TPG	Lead-Free	\$73.00

CAT. #	Lead Type	Price
TMR	Regular	\$69.00
TML	Light Wt.	\$79.00
TMG	Lead-Free	\$89.00

0.375mm Pb Equivalent Protection

CAT. #	Lead Type	Price
TPSR	Regular	\$53.00
TPSL	Light Wt.	\$57.00
TPSG	Lead-Free	\$66.00

CAT. #	Lead Type	Price
TMSR	Regular	\$53.00
TMSL	Light Wt.	\$57.00
TMSG	Lead-Free	\$66.00

Attached Thyroid Collars Info

All Attached Thyroid Collars include:

- 1) A snap ring for easy collar attachment.
- 2) Attach anytime with tethered strap.
- 3) Unique 2-in-1 pocket design.

[store collar in one pocket and use other pocket for accessories]

Thyroid collars are 6" tall x 6 1/2" wide, with a 24" strap.



Attached Thyroid option and our 2-in-1 pocket design!

Buckle Closure Add-On

If you prefer buckle closure to regular hook & loop, add "-B" to any of the collar numbers above and add \$5.00 to the price.

(Ex.: TYR-B would be \$49.00)



Buckle closure add-on available!

Breast-Guard

Female Breast Protection

Allows full viewing of anatomy EXCEPT sensitive breast tissue. Hook and Loop strap drapes around the neck. Adjust shields anywhere along the strap. 0.5mm Pb equiv. protection. Sets are two shields plus one strap

Two convenient sizes:

Adult: 6" wide x 7" high; Child: 4" wide x 4 1/2" high.

CAT #	Lead Type	Price Per Set	Size
BSR	Regular	\$44.00	A (Adult)
BSL	Light Wt.	\$54.00	C (Child)
BSG	Lead-Free	\$64.00	



Scolio-Guard

Unprotected neck for cervical view. Wrap around belt allows custom positioning of the open area. No-slip fabric underside. 0.5mm Pb equiv. protection.

Adult: 7.5" w x 26" L **Child:** 6.0" w x 20" L

CAT #	Lead Type	Price for Adult	Price for Child	Size
SCR	Regular	\$99.00	\$99.00	A, C
SCL	Light Wt.	\$119.00	\$119.00	A, C
SCG	Lead-Free	\$170.00	\$170.00	A, C



Scolio-Guard2

These breast guards offer the same as Scolio-Guards, but with the added **Hook and Loop** strip on the front to fold up out of the way when needed. Also available in Adult and Child sizes.

Adult: 7.5" w x 26" L **Child:** 6.0" w x 20" L

CAT #	Lead Type	Price for Adult	Price for Child	Size
SCR2	Regular	\$119.00	\$119.00	A, C
SCL2	Light Wt.	\$139.00	\$139.00	A, C
SCG2	Lead-Free	\$178.00	\$178.00	A, C



Semi-Guard

Extra Coverage Half Apronette

Convenient hook & loop closure. 0.5mm Pb equivalency protection. Excellent protection for patient and staff.



Default color is 61-Sapphire Blue

Size	W x L
S	12" x 12"
M	18" x 16"
L	24" x 20"
XL	24" x 24"
2XL	36" x 24"
Half	24" x 12"
Kids*	8" x 10"

*See Kids Lap-Guard PG.22

Lap-Guard

Buckle Closure Half Apronette

Unbeatable protection for patient and staff. Easy-to-use 1" buckle belt closure. 0.5mm Pb equivalency protection.



Semi-Guard and Lap-Guard Options

SEMI CAT #	LAP CAT #	Lead Type	Small	Med/Half	Large	X-Large	2XL
SGR	LGR	Regular	\$90.00	\$100.00	\$110.00	\$130.00	\$160.00
SGL	LGL	Light Wt.	\$110.00	\$120.00	\$130.00	\$150.00	\$190.00
SGG	LGG	Lead-Free	\$160.00	\$180.00	\$200.00	\$220.00	\$270.00

Semi-Guard Sets and Lap-Guard Sets

SEMI CAT #	LAP CAT #	Lead Type	Price Per Set
STR	LTR	Regular	\$448.00
STL	LTL	Light Wt.	\$518.00
STG	LTG	Lead-Free	\$748.00

Semi-Guard & Lap-Guard SETS include a 2-Peg Rack and the following apronettes:

Small (yellow) Large (royal blue)
Medium (berry) X-Large (green)



Pregnancy Semi & Lap-Guards

Half Apronette

Similar to Semi-Guard and Lap-Guard shown at the top of this page, but with 1.0mm Pb equivalency protection and a 2" belt on the Lap-Guard version.

Excellent protection for pregnant patients requiring x-ray.

1.0mm PREGNANCY SEMI-GUARD & LAP-GUARD

SEMI CAT #	LAP CAT #	Lead Type	Medium	Large	X-Large	2XL
SPL	LPL	Light Wt.	\$206.00	\$226.00	\$266.00	\$286.00
SPG	LPG	Lead-Free	\$324.00	\$344.00	\$364.00	\$384.00



Demi-Guard

Half Apron

Child to Adult gonadal protection. 0.5mm Pb equivalency front protection. Available with 1" buckle closure or convenient Waist Ring.

WAIST RINGS

Demi-Guard with Waist Rings 0.5mm

CAT. #	Lead Type	Small	Medium	Large	X-Large
DRR	Regular	\$126.00	\$136.00	\$146.00	\$156.00
DRL	Light Wt.	\$146.00	\$156.00	\$166.00	\$176.00
DRG	Lead-Free	\$166.00	\$186.00	\$196.00	\$206.00

Demi-Guard 0.5mm (Ring set of 4 sizes & rack)

CAT. #	Lead Type	Price
DSR	Regular	\$753.00
DSL	Light Wt.	\$823.00
DSG	Lead-Free	\$943.00



Size	Small	Medium	Large	X-Large
W x L	8" x 10"	12" x 12"	16" x 14"	20" x 16"

BUCKLE OR HOOK & LOOP CLOSURE

Demi-Guard with 0.5mm Pb Equivalent Protection

Buckle CAT#	Hook & Loop CAT#	Lead Type	Small	Medium	Large	X-Large
DBR	DVR	Regular	\$66.00	\$76.00	\$86.00	\$96.00
DBL	DVL	Light Wt.	\$76.00	\$87.00	\$113.00	\$144.00
DBG	DVG	Lead-Free	\$86.00	\$96.00	\$116.00	\$144.00

Demi-Guard 0.5mm (Buckle or Hook & Loop set of 4 sizes & rack)

Buckle CAT#	Hook & Loop CAT#	Lead Type	Price
DGR	DYR	Regular	\$523.00
DGL	DYL	Light Wt.	\$593.00
DGG	DYG	Lead-Free	\$713.00

Size	Small	Medium	Large	X-Large
W x L	8" x 10"	12" x 12"	16" x 14"	20" x 16"



ALSO AVAILABLE FOR ORDER!

Demi-Rack Only

#DDD | \$282.00



Gonadal-Guard

Adult Gonadal Shields

Attached belt with 1" buckle. 0.5mm Pb equivalency protection. Available in three sizes.



Gonadal-Guard (Individual)

CAT#	Lead Type	Small	Medium	Large
GNR	Regular	\$36.00	\$46.00	\$56.00
GNL	Light Wt.	\$46.00	\$56.00	\$66.00
GNG	Lead-Free	\$51.00	\$61.00	\$71.00

Gonadal-Guard (Set of all 3)

CAT#	Lead Type	Price	Size	(W x L)
GSR	Regular	\$120.00	S	4" x 4.5"
GSL	Light Wt.	\$140.00	M	6" x 7"
GSG	Lead-Free	\$170.00	L	7.5" x 9"

Diaper-Guard

Pediatric Protective Diaper

Specifically designed radiation protective diapers accommodate infants up to small children. Color-coded regular lead diapers: small (red), medium (green), large (yellow). 0.5mm Pb equivalency protection. **Regular lead.**

CAT#	Back x Length x Front	Size	Price
GBR-S	9.5" x 11" x 5"	Small	\$84.00
GBR-M	12.5" x 12" x 7.5"	Medium	\$94.00
GBR-L	15.5" x 13" x 9"	Large	\$104.00
GBS-3	Set of All 3	Set of 3	\$254.00



Gonadal-Guard

Pediatric Gonadal Shields

Specifically designed sets will accommodate neonatal up to small children. Gender appropriate gonadal shielding. 0.5mm Pb equivalency protection.

Male | Set of 4 Blue Gurads

CAT#	Lead Type	Price Per Set
GMR4	Regular	\$85.00
GML4	Light Wt.	\$105.00
GMG4	Lead-Free	\$135.00



Female | Set of 5 Pink Gurads

CAT#	Lead Type	Price Per Set
GFR5	Regular	\$95.00
GFL5	Light Wt.	\$115.00
GFG5	Lead-Free	\$155.00



Diaper-Guard/Brachytherapy

Adult Protective Diaper

Specifically designed radiation protective diapers for adults. 0.5mm Pb equivalency protection. Convenient tie strings at each hip. **Regular lead.**

#GDR | \$125.00



Kids Lap/Semi-Guard | Half Apronette

With a width of 8" and a length of 10", this semi-guard is excellent protection for small patients. LAP uses 1" buckle closure and SEMI uses hook & loop. 0.5mm Pb equivalency protection.

LAP CAT#	SEMI CAT#	Lead Type	Child
LKR	SKR	Regular	\$80.00
LKL	SKL	Light Wt.	\$90.00
LKG	SKG	Lead-Free	\$120.00



Different color and patten options on PGS 29-30!

Child-Guard | Front Protection For Kids

EZ Fit, Front Protection

Ideal for our younger patients, creating safer and more secure protection than lead drapes or shields. Front protection apron with buckle closure. 0.5mm Pb equiv. front protection.

Shoulder Pads concealed for comfort and durability!
Width x Length (18" x 24")

LAP CAT#	Lead Type	Child
KGR	Regular	\$189.00
KGL	Light Wt.	\$249.00
KGG	Lead-Free	\$299.00



See our "Best-Fit" sizing guide on PG 2!

Kiddie Covers | Kid-Friendly Protection

What makes Kiddie Covers Unique:

Decrease anxiety and fear with a “hug” from a friendly buddy. All purpose lead blanket protects children from radiation. Protective panel is 18” x 18” and safe for all modalities. [Regular lead.](#)




0.5mm Pb Equivalency

CAT#	Design Style	Price
KK-LBUG	Ladybug	\$180.00
KK-BBEE	Bumble Bee	\$180.00
KK-TURT	Turtle	\$180.00
KK-BFLY	Butterfly	\$180.00
KK-PUPP	Puppy	\$180.00
KK-FROG	Frog	\$180.00
KK-CATT	Cat	\$180.00
KK-MNKY	Monkey	\$180.00



0.375mm Pb Equivalency

CAT#	Design Style	Price
KS-LBUG	Ladybug	\$124.00
KS-BBEE	Bumble Bee	\$124.00
KS-TURT	Turtle	\$124.00
KS-BFLY	Butterfly	\$124.00
KS-PUPP	Puppy	\$124.00
KS-FROG	Frog	\$124.00
KS-CATT	Cat	\$124.00
KS-MNKY	Monkey	\$124.00



Omni-Guard

Protective Flexible Shield

Used over a patient or equipment, protects from radiation exposure. Used underneath a patient, prevents scatter radiation. All purpose covering for patients or equipment. 0.5mm Pb equivalency protection. Choose any color fabric.

Size	Width x Length
S	18" x 18"
M	24" x 24"
L	24" x 36"
X	24" x 48"



0.5mm Pb Equivalency

CAT. #	Lead Type	Small	Medium	Large	X-Large
OMR	Regular	\$124.00	\$144.00	\$184.00	\$214.00
OML	Light Wt.	\$154.00	\$203.00	\$284.00	\$351.00
OMG	Lead-Free	\$214.00	\$264.00	\$324.00	\$468.00

0.375mm Pb Equivalency

CAT. #	Lead Type	Small	Medium	Large	X-Large
OMSR	Regular	\$94.00	\$124.00	\$154.00	\$184.00
OMSL	Light Wt.	\$144.00	\$184.00	\$224.00	\$274.00
OMSG	Lead-Free	\$204.00	\$254.00	\$334.00	\$474.00



NEW Add-On Conversions

For All Omni-Guards

Make your Omni-Guard more versatile! Want to hang your omni-guard like a curtain or drape? Convert any omni-guard by adding grommets or loops!*

CAT. #	Type	QTY	Price
OM-GR	Grommets	Up to 5	\$25.00
OM-LP	Loops (Tabs)	Up to 5	\$25.00

Shop color and embroidery options on pages 28-30!



Custom Omni-Guards

Don't see the size you need? We are happy to provide a quote based on your needs! Simply fill out the information below and submit for a quote today!

Customer Information	
Name	
Phone Number	
Contact Email	
Address Line 1	
Address Line 2	


Omni-Guard Quote Information	
Length (in.)	
Width (in.)	
Grommets (Y/N)	
Tabs (Y/N)	
Drape Color (PG. 29-30)	

*Desired location and distance between grommets (maximum width of 48") must be noted on any Omni-Guard add-on order.

Belt-Guard | Back Support Belt

Helps reduce back stress. Can be used with any apron. Features a 6" elastic back support belt with 2" buckle closure. The same belt as used on the Ortho-Guard apron.

CAT#	Style/Size	Price
BGB	Regular	\$50.00
BGB-2	X-Large	\$60.00



Guard-Tote | Apparel Tote Bag

Convenient tote for all your protective wear. Durable rip-stop nylon fabric with zippered closure. 28" long to minimize folding your aprons. EZ-Carry handle and shoulder strap included. Default color for the Guard Tote is Anti-Microbial Black Onyx Diamond Ripstop but available in other colors at no additional charge. (shown in 64-Onyx)

#ACB-28 | \$110.00



Hat-Guard

Protection for the top of head. **Hook and Loop** fitted closure. 0.50mm Pb equiv. sides, 0.25mm Pb equiv. top.

CAT#	Lead Type	Price
HCL	Light Wt.	\$154.00
HCG	Lead-Free	\$174.00

Size	XX
S	21"
M	22"
L	23.5"



Palm Guard Radiation Protective Mittens

Full palm protection with thumb extension for grip. 0.5mm Pb equivalency dip molded core. Durable, soft vinyl covering.

Palm-Guard Slitted Mittens

CAT#	Quantity	Price
PSP	Pair (L & R)	\$531.00



Palm-Guard Mittens (no slit)

CAT#	Quantity	Price
PGP	Pair (L & R)	\$445.00



Hand Guard Radiation Protective Gloves

Five fingered protection for full dexterity. 0.5mm Pb equivalency dip molded core. Durable, soft navy vinyl covering. 15". Sold as a pair.

Hand-Guard/ 5 Finger

CAT#	Quantity	Price
HGP	Pair (L & R)	\$376.00



Finger-Guards | Now sold as a single pair!

Radiation reducing Finger-Guard offers extra sterile hand protection from secondary radiation. Multiple options to choose from! **NEW sold as a single pair!**

CAT#	Description	Glove Sizes	NEW QTY	Price
PLF-1	0.18mm Glove	[6.5] [7] [7.5] [8] [8.5] [9]	1 Pair	\$110.00
PFC-1	0.23mm Glove	[6.5] [7] [7.5] [8] [8.5] [9]	1 Pair	\$99.00
PRL-1	0.35mm Glove	[6.5] [7] [7.5] [8] [8.5] [9]	1 Pair	\$99.00

NOW SOLD AS A SINGLE PAIR!

Attenuation for Finger-Guard Broad Beam					
CAT#	Description	60 KvP	80 KvP	100 KvP	120 KvP
PLF	0.18mm Glove	56.5%	47.1%	41.6%	37.7%
PFC	0.23mm Glove	48.0%	40.0%	32.0%	28.0%
PRL	0.35mm Glove	60.0%	53.0%	46.0%	41.0%



PLF-1

Superior tactile sensitivity with excellent attenuation.

Lead-free/Powder-free/Latex-free

PFC-1

Excellent tactile sensitivity with good attenuation

Lead-free/Latex-free

PRL-1

Reasonable tactile sensitivity with the best attenuation

Lead-free/Powder-free/
Contains Natural Latex/ Beaded cuff/
Textured finger tips

("mm" is actual thickness of glove in millimeters, NOT lead equivalency. Test results in accordance with ASTM Test Method F2547-06.)

Custom Apparel Quote Form

Please fill out the following form and submit for a pricing quote on all custom apparel. All custom sized apparel will include an additional fee for each item ordered. Custom sized apparel is non-refundable, non-returnable.

Custom Sized Apparel

Please use the provided CAT numbers and pricing for each custom apparel item.

ACP - **\$125** fee for custom Aprons, Vests & Skirts

ACP2 - **\$45** fee for other custom apparel items (ex: thyroid guards, gonadal guards)

*3X is the maximum size for standard (listed) pricing. All sizes exceeding 3X will be considered custom products. Additional fees apply.

APRON CAT.# _____ DESCRIPTION _____ LEAD TYPE _____

FABRIC (pg 29-30) _____ BINDING COLOR (pg 28) _____


EMBROIDERY (OPTIONAL) LINE 1, UP TO 24-LETTERS (FREE) _____

EMBROIDERY LINE 2, UP TO 24-LETTERS (\$20.00) _____

EMBROIDERY THREAD COLOR (pg 28) _____ LETTER STYLE (pg 28) _____

ADD ATTACHED THYROID COLLAR? YES NO CAT.# _____ (COLOR WILL MATCH APRON)

Getting Accurate Measurements (all in inches)



Please submit this information with your order.

GENDER: _____

HEIGHT: _____

BUST/CHEST: _____

Circumference around the widest part of your chest.

SIDE-to-SIDE: _____

The widest part of your body that you want covered by a Front Protection Apron. Please specify the this area (circle): chest / tummy / hips

WAIST: _____

Standard Apron: measure your waistline.

Skirt/Kilt: measure where you want the skirt/kilt top (may be higher than your regular waistline).

HIP: _____

Circumference around the widest part of your hips.

LENGTH: _____

FRONT APRON or WRAP-AROUND: _____

Sternal notch to just above the knee.

VEST/TORSO: _____

Sternal Notch to desired length below your waist.

SKIRT/KILT: _____

Desired skirt/kilt top (may be above your regular waistline) to just above your knee.

Protective Apparel Add-Ons

FREE Embroidery Sewn-In

Select products come with the first line (up to 24 characters) of embroidery **FREE!** Add flexibility with an embroidered patch. Choose from a style below. Add flexibility with an embroidery patch! *Embroidered items are non-returnable & non-refundable.*

Available Styles

JOHN SMITH, MD

John Smith, MD

John Smith, MD

All Caps (default)

Block

Script

#EML2 | \$20 for second embroidery line, up to 24 letters.

Embroidery Add-On Patch

This flexible option has your information embroidered onto a patch. Hook & Loop on the back of the patch and on the garment secure the patch in place.

#EMP (Patch) | \$25.00



(Patch above pocket)



(Hook and loop closure)

Embroidery Colors

- White* (default)
- Black* (default)
- Lt. Gray*
- Charcoal*
- Royal Blue*
- Red*
- Yellow*
- Orange*
- Tan*
- Brown*


Binding Colors

- Black** (default)
- Light Gray** (default)
- Navy**
- Royal Blue**
- Red**
- Yellow**
- Orange**
- Tan**
- Pink**
- Purple**
- Green**
- Burgundy**

FREE Shoulder Strap

Simply ask for it with your order. Hook & loop Strap is attached to one shoulder and reaches across to the other.

#EMS | FREE



Helps keep shoulders high and comfortable.

Protective Media Guide

MEDIA	BENEFITS	PROTECTION	ATTENUATION	WEIGHT
Regular Lead	Most economical and recyclable	100% Lead	0.5mm Pb Equiv. Protection avg. 97% attenuation over 80, 100, 120 KVP	9 - lbs. <small>(based on medium front protection apron)</small>
Lightweight Lead	Lighter than regular lead and recyclable	Lead-Tin	0.5mm Pb Equiv. Protection avg. 96.8% attenuation over 80, 100, 120 KVP	8 - lbs. <small>(based on medium front protection apron)</small>
Lead Free	Lightest Weight, Lead-Free, Environmentally Friendly	Tungsten-Antimony	0.5mm Pb Equiv. Protection avg. 96.6% attenuation over 80, 100, 120 KVP	7 - lbs. <small>(based on medium front protection apron)</small>


Protective Apparel Fabric Options

NYLON FABRICS

FEATHER RIPSTOP FABRIC Light Weight Nylon, Filament Reinforced Fabric



DIAMOND RIPSTOP FABRIC Medium Weight Nylon, Filament Reinforced Fabric




NYLON SOLIDS Liquid Repellent, Light & Supple



SPECIALTY FABRICS


COLORFUL PRINTS Liquid Repellent, Medium Weight



NEW 69-Skulls




70-Blue Plaid



71-Leopard



72-Lots•0•Dots




73-Sleepy Sheep



74-Playful Pups



83-Cool Waves



77-Life Vines



78-Waterfall




80-Red Paws



81-Black Paws



76-Blue Paws



82-Pink Camo




85-Digital Camo



86-Camo Night

Prints are shown
SMALLER than
actual size!


TEXTURED FABRIC Textured High Tech Fabrics



15-Diamondback



16-Crocogator



18-Bronze Tech

New Premium Fabric Add-On


NEW Coolwear Fabric

Designed for optimum cooling comfort!

Our new Coolwear fabric add-on provides premium cooling comfort all day long! **This premium fabric can now be added to ANY garment!** The polyester mesh lining features air vented space and stretch-in-width for extreme comfort and fit. This innovative add-on is designed for ultimate performance and comfort.

- CP50** - \$50 add-on for aprons, vests, skirts, kilts, dental drapes, omni-blankets, & kiddie covers
- CP40** - \$40 add-on for half aprons and scolio-guards
- CP25** - \$25 add-on for thyroid collars, hats, & gonadal shields

Features premium comfort
with all day cooling!



Ultra-Guard

Ultra-light and super stylish.

Curved lens adds to lateral protection. Regular hinges with wrap around frame. Universal bridge support. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.
 Weight: 67 grams.
 Frame: TR-90 Nylon.

Model#	Color	Prescription Options
LER-9SPL	Silver	Single Vision \$320.00
LER-9KPL	Black	Bifocal \$340.00
LER-9TPL	Taupe	
LER-9PPL	Pink	
Plano Lens	\$195.00	



Avant-Guard

Two-tone colorful Ultra-Guards.

Curved lens adds to lateral protection. Regular hinges with wrap around frame. Universal bridge support. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.
 Weight: 67 grams.
 Frame: TR-90 Nylon.

Model#	Color	Prescription Options
LER-ACPL	Carbon/Gray	Single Vision \$320.00
LER-APPL	Purple/Gray	Bifocal \$340.00
LER-AYPL	Yellow/Blue	
Plano Lens	\$195.00	



NEW Elite Guard Wrap Around Glasses

Super Sleek curved lens provides all the great features of a High Quality Wrap Around Safety Frame. Techno-Aide Elite Guard-lightweight radiation protective glasses offer a stylish wrap-around style frame. Features rubberized soft grip temple bars for a secure fit. Rubberized soft grip nose pads. The Elite Guard fits best for Large size faces. Includes hard case, retainer cord and cleaning cloth

Lead Equiv.: 0.75mm lens.
 Weight: 58 grams.
 Frame: High-impact TR-90 nylon.

Model#	Color	Prescription Options
LER-EB-PL	Black	Single Vision \$320.00
LER-ES-PL	Silver	Bifocal \$340.00
Plano Lens	\$195.00	



Turbo-Guard

For that sharp and stream-lined look.

Large lens area for excellent coverage. Rubberized temple bars for comfortable grip. Rubber, ultra-soft, fixed nose pad. Available in almost all prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.
 Weight: 65 grams.
 Frame: Impact resistant nylon.

Model#	Color	Prescription Options
LER-TRPL	Red	Single Vision \$320.00
LER-TKPL	Black	Bifocal \$340.00
LER-TSPL	Silver	
Plano Lens	\$195.00	



Sure-Guard

Stylish design for medium faces.

Large lens area for excellent coverage. Regular hinges with wrap around frame. Rubber, ultra-soft, fixed nose piece. Available in almost all prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.
 Weight: 63 grams.
 Frame: TR-90 Nylon.

Model#	Color	Prescription Options
LER-6SPPL	Silver	Single Vision \$335.00
LER-6KPL	Black	Bifocal \$355.00
LER-6RPL	Red	
LER-6BPL	Blue	
Plano Lens	\$215.00	



Astro-Guard

Side Shield protection.

Attractive and sturdy women's frames. Flexible temples fit almost any face. Excellent visibility. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm front and sides.
 Weight: 66 grams.
 Frame: Impact resistant plastic.

Model#	Color	Prescription Options
LES-AFYK	Black/Yellow	Single Vision \$335.00
LES-AFPW	Purple/White	Bifocal \$355.00
LES-AFGP	Gray/Pink	
LES-AFOP	Brown/Orange	
Plano Lens	\$215.00	



Wrap-Around

The “original” with great side protection.

Generous lens area increases protection. Spring hinges and adjustable nose. Available in most prescriptions. Includes premium case. Prescription lens available.

Lead Equiv.: 0.75mm lens.	Model#	Color
Frame: Impact resistant plastic.	LES-53PLT	Tortoise
	LES-53PLB	Blue
	Plano Lens	\$215.00



Viva-Guard

Lively designs for small faces.

Smaller frame for excellent fit of smaller faces. Regular hinges with wrap around frame. Ultra soft bridge support. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.	Model#	Color	Prescription Options
Weight: 60 grams.	LER-VFPL	Fancy Flight	Single Vision \$365.00
Frame: Impact resistant polycarbonate.	LER-VVPL	Cherry Vines	
	LER-VSPL	Rad Skull	
	Plano Lens	\$235.00	



Power-Guard

Colortone accents increase comfort.

Large lens area for excellent coverage. Soft colored endpiece for comfort and fit. Padded bridge, temples, and endpiece. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.	Model#	Color	Prescription Options
Weight: 64 grams.	LER-PBPL	Blue/Orange	Single Vision \$365.00
Frame: Impact resistant polycarbonate.	LER-PSPL	Silver/Pink	Bifocal \$385.00
	LER-PKPL	Black/Gray	
	Plano Lens	\$235.00	



Secure-Guard

Vented goggle for extreme protection.

Adjustable headband keeps goggles securely on head and prevents slippage. Neoprene foam cushioning for facial comfort. Large lens pocket allows for maximum viewing area and vented to prevent fogging. Prescription lens available.

Lead Equiv.: 0.75mm lens.	Model#	Color	Prescription Options
Weight: 66 grams.	LES-RKPL	Black	Single Vision \$320.00
Frame: Impact resistant plastic.	Plano Lens	\$175.00	Bifocal \$340.00



Cover-Guard

Comfortable protection over regular glasses.

Excellent front and side visibility. Provides limited splash protection. Can be worn over prescription eyewear. Regular hinges with molded nose piece. Includes 0.75mm Pb glass side shields. Includes premium case.

Lead Equiv.: 0.75mm front and sides.	Model#	Color
Weight: 85 grams.	LES-CKPL	Black
Frame: Black rubberized nylon.	LES-CWPL	White
	LES-CTPL	Tortoise
	Plano Lens	\$150.00



Mobile Apron Z-Rack

Z-Base mobile rack nests with others when not in use. Heavy duty casters on 12G steel base. Chromed-steel posts and hang bar. 60" long x 60" tall bar allows hanging without creasing. Ships @ 35 lbs.

Rack ONLY
#ARM-UZNH
\$280.00

OR

Rack with
Ten Hangers
#ARM-UZ10
\$566.00



NEW Mobile Apron Z-Rack Middle Bar Add-On

Middle Bar addition to hang vests and skirts! 63" bar will hold up to 150 lbs. (Mobile apron z-rack sold separately)


Bar Add-on
#ARM-UBAR
\$58.00



NEW Mobile Apron Z-Rack Storage Base Add-On

Bottom shelf add-on for storing misc. protective accessories! 62" shelf will hold up to 150 lbs. (Mobile apron z-rack sold separately)

Shelf Add-on
#ARM-UBASE
\$195.00



Apron Hangers

Designed for standard aprons, panoramic, vests, skirts, etc. Top hook swivels for easy hanging. 21" long chromed-steel with anti-slip tips. Use on ARM-UZNH rack, in lockers, or behind a door. Single ships @ 2 lbs.

Single Hanger
#ARH-01
\$29.00



Set of 10 Hangers
#ARH-10
\$286.00



Mobile Apron Rack

Ten (10) swing-arm mobile garment stand. Durable all-steel construction. 24" x 24" H-Base for stability and easy cleaning. E-Z-Roll swivel locking casters. Convenient handle for easy mobility. White enamel finish. Overall height is 55". Some assembly required. Ships @ 110 lb.

Standard Mobile Rack

#ARM-100
\$1165.00



E-Z Roll, 3" medical grade, swivel, locking casters!

Standard WITH Glove Rack

#ARM-G60
\$1240.00



Arms swing inward to conserve space!

Easy-Clean H-Base!



(Aprons sold separately)

5 Arm Apron Rack

Five (5) swing-arm mobile garment stand. Durable all-steel construction with white enamel finish. 20" x 20" H-Base for stability and easy cleaning. E-Z roll swivel locking casters. Convenient handle for easy mobility. Arms extend to 29.5" outward. We offer the ability to choose whether you need right or left swing arms. Ships @ 80 lbs. motor freight.

#ARM-105
\$880.00



You choose the right or left side!



Arms swing inward to conserve space!

Easy-Clean H-Base!



(Aprons sold separately)

E-Z Roll, 3" medical grade, swivel, locking casters!

Multi-Tubular Apron Rack

The perfect solution safe apron storage. Durable all-steel construction. A 24" deep by 32" wide" H-base offers stability and easy cleaning. Features E-Z Roll, 3" medical grade swivel locking casters. Convenient handle for easy mobility from room to room. White baked enamel finish. Overall height is 49" and width is 32". Rods are 27" long and spaced 4" apart for easy clearance and storage. Stores up to 5 aprons in tight spaces without crinkling or damaging aprons.

#ARM-500 | \$950



E-Z Roll, 3" medical grade, swivel, locking casters!



Convenient handle for easy mobility and storage!

Easy-Clean H-Base!



Holds up to 5 Aprons at once!

(Aprons sold separately)

Available Steel Colors*			
White	Black	Platinum	Almond
			
Default	Textured		Textured




*Non-default color options are considered custom items. Custom items are non-returnable and non-refundable. 1-year limited warranty on materials and workmanship for most products offered. Please see full terms and conditions for more information and included products.

NEW Mobile Combination Protection Rack

The ultimate option for serious radiation protection storage. This versatile rack holds 6 aprons, 3 pairs of gloves, and 6 pairs of radiation protection glasses. Comes with 6 heavy duty apron hangers for quick apron access. Durable all-steel construction. A 24" deep by 28 wide" H-base offers stability and easy cleaning. Features E-Z Roll, 3" medical grade swivel locking casters. White baked enamel finish. Overall height is 79", width is 30" and depth is 24". The perfect height to hang aprons avoiding damage from crinkling. Soft rubber rings and 6 high grade mounts create the ultimate protection for radiation glasses at all times.

#ARM-600 | \$950



Available Steel Colors*			
White	Black	Platinum	Almond
			
Default	Textured		Textured

Holds up to 6 aprons, 3 pairs of gloves and 6 pairs of protective glasses at once!

*Non-default color options are considered custom items. Custom items are non-returnable and non-refundable. 1-year limited warranty on materials and workmanship for most products offered. Please see full terms and conditions for more information and included products.

NEW Mobile Lead Apron Locker

These industrial strength storage cabinets will feature heavy-duty casters for mobility and have 22 gauge welded steel construction. Includes a full width adjustable shelf and a full width hanger rod. The 230 lb. capacity shelf meets ANSI/BIFMA shelf capacity standards. Double doors lock with 3 point latching and recessed handles that prevent banging and catching with a flush surface. Powder coat finish. The 5" polypropylene swivel casters (2 with brakes) are attached with a 6"W heavy-gauge steel plate that bolts through to 1-1/4" angle iron inside the cabinet, both running the length of the cabinet. Steel caster components are zinc treated to resist corrosion. Total caster capacity is 1100 lbs. Rod capacity 500 lbs. Locks with key. Available in charcoal and gray. Item will shipped assembled.



Features	
Width	46
Depth	24
Height	78
Adjustable Shelves	1
Construction	20/22 Gauge Steel
Door Quantity	2
Handle Type	Recessed
Shelf Capacity	230 lbs.
Wheel Type	Polypropylene 4 Swivel, 2 Brake

Charcoal Locker Only

MLALC | \$1330

Gray Locker Only

MLALG | \$1330

Locker & Hanger Bundles

MLALG-10 | \$1595

MLALG-10 | \$1595

Locker in Gray



*Aprons sold separately

Wall Mounted Apron Racks

Ten or five swing-arm wall-mount apron racks. Mounting holes on 16" centers. All-steel durability with white enamel finish.



10-Arm #ARW-10 \$436.00



5-Arm, Swing Right #ARW-R5 \$260.00



5-Arm, Swing Left #ARW-L5 \$260.00

Glove Adapter

Retrofits to Mobile Rack. Holds six (6) 18" gloves. White enamel finish.

#ARG-06
\$103.00



Apron/Glove Racks

Designed for one apron, but will hold more. All steel, heavy duty construction and a white finish. 16" & 24" holes match most wall studs. Glove holders mounted onto back panel. 28.5" long x 1.75" diameter bar prevents creasing.

WITHOUT GLOVE HOLDERS

#AGR-105
\$187.00



WITH TWO GLOVE HOLDERS

#AGR-110
\$251.00



WITH FOUR GLOVE HOLDERS

#AGR-115
\$272.00



Apron Peg Racks

Help extend the life of your aprons. Our apron racks are just the thing for you! 6" angled pegs keep garments from slipping. All-steel, white enamel finish.

Custom sized racks may be available!

2-Peg
(for 1 apron)
#ARP-02
\$72.00



4-peg
(for 2 aprons)
#ARP-04
\$124.00



6-peg and 8-peg racks: mounting holes on 16" and 24" centers.

6-peg
(for 3 aprons)
#ARP-06
\$145.00



8-peg
(for 4 aprons)
#ARP-08
\$156.00



16-peg
(for 8 aprons)
#ARP-16
\$398.00



20-peg
(for 10 aprons)
#ARP-20
\$461.00






FREE Protective Apparel Tracking and Inventory System!

Exclusive from Techno-Aide!

All new Techno-Aide garments include a scanable QR code for automatic tracking. Upload garment data into your unique ApronHub portal for permanent inventory management. Receive alerts for inspection dates through a customizable format. Retrofit QR code tags onto your existing Techno-Aide or non-Techno-Aide garments. Track inspection history, location, user, etc., through a customizable reports menu. 2D Scanners are available through Techno-Aide or you can use your own 2D-scanner. ApronHub has everything you need to quickly and easily meet the highest QA standards. For more information visit www.apronhub.com today!







Accessories

Tools for Scanning and Labeling Protective Apparel

Stay one step ahead of the inspectors and enhance your QA! All new garments manufactured by Techno-Aide include a scanable QR Code Serial Number label for automated or manual tracking. However, we realize that you may have older, unlabeled garments or garments from another manufacturer. To update all your protective garments, we have designed this retrofit labeling kit so that your garments can be tracked utilizing the most up-to-date tracking technology available.

CAT #	Description	Price
AHSCB	2D Scanner with Cable	call
AHSXW	2D Scanner Wireless	call
AHETL	Eyelet Tool	call
AHES50	Silver Eyelets (50)	\$25.00
AHQRC	QR Code Retrofit Labels (100)	\$68.00
AHDT19	Date Tags: 2019 (50)	\$ 14.00
AHDT20	Date Tags: 2020 (50)	\$ 14.00





ADMINISTRATION/MANUFACTURING

7117 Centennial Blvd. Nashville, TN 37209

Orders can be placed online or sent to orders@techno-aide.com

For official Quotes: quotes@techno-aide.com

Any questions can be directed to customerservice@techno-aide.com

P: 800.251.2629 or 615.350.7030

F: 800.554.6275 or 615.350.7879

techno-aide.com