

radiology



IMAGING IN MOTION

 **ustom X-Ray Service, Inc.**
Digital Equipment Sales & Service

CONTENTS

PROTECTIVE APPAREL

2-25

SHIELDING / BARRIERS

26-31

GENERAL ACCESSORIES

32-45

ULTRASOUND / MRI

46-51

POSITIONING

52-73

TABLE PADS

74-83

PATIENT HANDLING

84-93

MARKERS

94-107

DARKROOM / VIEWING

108-109

FILING SYSTEMS

110-115

FUSION ERGO

116-119

INDEX ON PAGE 120-121

Copyright 2018, Techno-Aide, Inc., all rights reserved.

FREE STUFF on applicable garments:

Double-Fold Binding
Shoulder Pads
Double Pocket


Embroidery
Recycling

MEDIA	BENEFITS	PROTECTION	ATTENUATION	WEIGHT
Regular Lead	Most Economical and Recyclable	100% Lead	0.5mm Pb Equiv. Protection averages 97% Attenuation over 80,100,120 KVP	9-lbs (based on medium front protection apron)
Lightweight Lead	Lighter than Regular Lead and Recyclable	Lead-Tin	0.5mm Pb Equiv. Protection averages 96.8% Attenuation over 80,100,120 KVP	8-lbs (based on medium front protection apron)
"LMG"*	Lightest Weight, Lead-Free, Environmentally Friendly	Tungsten-Antimony	0.5mm Pb Equiv. Protection averages 96.6% Attenuation over 80,100,120 KVP	7-lbs (based on medium front protection apron)

*Media is IEC/EN 61331-1 / ASTM Approved.

FDA REGISTERED / Every batch of protective media is tested by an independent third-party health physics laboratory and actual attenuation figures reported.

AH APRON-HUB Accessories



Tools for Scanning and Labeling Protective Apparel

Stay one step ahead of the inspectors and enhance your QA!


All new garments manufactured by Techno-Aide include a scanable QR Code Serial Number label for automated or manual tracking. However, we realize that you may have older, unlabeled garments or garments from another manufacturer. To update all your protective garments, we have designed this retrofit labeling kit so that your garments can be tracked utilizing the most up-to-date tracking technology available.



cat. #	description	price
AHSCB	2D Scanner with Cable	call
AHSXW	2D Scanner Wireless	call
AHETL	Eyelet Tool	call
AHES50	Silver Eyelets (50)	\$24.00
AHQRC	QR Code Retrofit Labels (100)	\$65.00
AHDT17	Date Tags: 2017 (50)	\$12.50
AHDT18	Date Tags: 2018 (50)	\$12.50



Protective Apparel Tracking and Inventory System!



www.apronhub.com

The exclusive online apparel tracking and inventory system from Techno-Aide.

All new Techno-Aide garments include a scannable QR code for automatic tracking. Upload garment data into your unique ApronHub portal for permanent inventory management. Receive alerts for inspection dates through a customizable format. Retrofit QR code tags onto your existing Techno-Aide or non-Techno-Aide garments. Track inspection history, location, user, etc., through a customizable reports menu. 2D Scanners are available through Techno-Aide or you can use your own 2D-scanner.

ApronHub has everything you need to quickly and easily meet the highest QA standards. For more information go to: www.apronhub.com today!



Most Popular Apron!



Buckle

EZ-Guard

All Purpose Front Protection Apron

EZ-On / EZ-Off. Available with hook & loop or 2" buckle closure. 0.5mm Pb equivalency protection (0.375mm available). Free dual pocket and shoulder pads.

EZ-GUARD with Hook & Loop Closure

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
EVR	Regular	\$220.00	\$240.00	\$260.00
EVL	Light Wt.	\$280.00	\$300.00	\$320.00
EVG	LMG Lead-Free	\$330.00	\$350.00	\$370.00

EZ-GUARD with Buckle Closure

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
EBR	Regular	\$220.00	\$240.00	\$260.00
EBL	Light Wt.	\$280.00	\$300.00	\$320.00
EBG	LMG Lead-Free	\$330.00	\$350.00	\$370.00

Flex-Guard

Flexible-Back Apron

See Best-Fit Sizing Chart on page 22.

Stretch back provides greater comfort by placing more weight on the hips and less weight on the shoulders. Available with hook & loop or 2" buckle closure. 0.5mm Pb equivalency protection (0.375mm available). Free dual pocket and shoulder pads.



Buckle

FLEX-GUARD with Hook & Loop Closure

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
FVR	Regular	\$230.00	\$250.00	\$270.00
FVL	Light Wt.	\$290.00	\$310.00	\$330.00
FVG	LMG Lead-Free	\$340.00	\$360.00	\$380.00

FLEX-GUARD with Buckle Closure

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
FBR	Regular	\$230.00	\$250.00	\$270.00
FBL	Light Wt.	\$290.00	\$310.00	\$330.00
FBG	LMG Lead-Free	\$340.00	\$360.00	\$380.00

COLORS & FABRICS: page 25.
Default Fabric: **Royal Blue nylon.**
2-Year Limited Warranty.



Surgery-Guard

Surgical Drop-Off Apron

Easily Slides Off under Sterile Gown

No compromise of sterile field. Crossover shoulder-to-hip hook & loop closure. 0.5mm Pb equiv. front protection. Free dual pocket and shoulder pads.

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
SRR	Regular	\$220.00	\$240.00	\$260.00
SRL	Light Wt.	\$280.00	\$300.00	\$320.00
SRG	LMG Lead-Free	\$330.00	\$350.00	\$370.00



OR-Guard

Shoulder Drop-Off Apron

Slides Off under Sterile Gown

No compromise of sterile field. Flexible back, like the Flex-Guard. Break-away shoulders and hips. 0.5mm Pb equiv. front protection. Free dual pocket and shoulder pads.

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
ORR	Regular	\$260.00	\$280.00	\$300.00
ORL	Light Wt.	\$320.00	\$340.00	\$360.00
ORG	LMG Lead-Free	\$370.00	\$390.00	\$410.00



Ortho-Guard

Apron with attached Belt-Guard

Built-in Back-Guard for Maximum Support

Tethered belt can be adjusted vertically and laterally for effective back support. Back strap assures best fit and comfort. 0.5mm Pb equiv. front protection. Free dual pocket and shoulder pads.

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
OGR	Regular	\$260.00	\$280.00	\$300.00
OGL	Light Wt.	\$320.00	\$340.00	\$360.00
OGG	LMG Lead-Free	\$370.00	\$390.00	\$410.00

Shoulder Strap Included!



Econo-Guard

Economy Apron

Excellent Protection, Value Priced

1" buckle closure (long enough to tie). 0.5mm Pb equivalency front protection. Available in any fabric. See Best-Fit size chart on page 22.

ECONO-GUARD

Ineligible for attached thyroid collar or shoulder strap.

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
EGR	Regular	\$150.00	\$170.00	\$190.00
EGL	Light Wt.	\$210.00	\$230.00	\$250.00
EGG	LMG Lead-Free	\$260.00	\$280.00	\$300.00

COLORS & FABRICS: page 25.

Default Fabric: **Royal Blue nylon.**

2-Year Limited Warranty.



Vest-Wrap

Less Bulk, Better Fit for Freedom of Movement

Separate Men and Women's styles for best comfort. Distributes weight between shoulders and hips. 0.5mm Pb equivalency protection front, 0.25mm back. Hook & loop closure with flexible back. Generous armholes for EZ-fit. Free dual pocket and shoulder pads.

VEST-WRAP Female

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
VGR	Regular	\$200.00	\$220.00	\$240.00
VGL	Light Wt.	\$260.00	\$280.00	\$300.00
VGG	LMG Lead-Free	\$310.00	\$330.00	\$350.00

VEST-WRAP Male

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
VBR	Regular	\$200.00	\$220.00	\$240.00
VBL	Light Wt.	\$260.00	\$280.00	\$300.00
VBG	LMG Lead-Free	\$310.00	\$330.00	\$350.00

COLORS & FABRICS: page 25.
Default Fabric: **Royal Blue nylon.**
2-Year Limited Warranty.



**Newly Tailored
to assure the most
Comfortable Fit!**

**See Best-Fit
Sizing Chart
on page 22.**



Vest-Guard

New Design: Saves Weight, Better Fit

Separate Men and Women's styles for best comfort. Full front overlap design saves weight. Distributes weight between shoulders and hips. 0.5mm Pb equivalency protection front, 0.25mm back. Hook & loop closure is fast and easy. Free dual pocket and shoulder pads.

VEST-GUARD Female

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
VFR	Regular	\$200.00	\$220.00	\$240.00
VFL	Light Wt.	\$260.00	\$280.00	\$300.00
VFG	LMG Lead-Free	\$310.00	\$330.00	\$350.00

VEST-GUARD Male

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
VMR	Regular	\$200.00	\$220.00	\$240.00
VML	Light Wt.	\$260.00	\$280.00	\$300.00
VMG	LMG Lead-Free	\$310.00	\$330.00	\$350.00

(Vests do not include the two 1" safety buckles as shown. If you would like these as a FREE Add-On, please specify at the time of ordering.)

Skirt/Kilt-Guard

All Day Comfort for the Lower Body

Separate Men and Women's styles for best comfort. Weight is distributed evenly around hips. Wrap around protection: .5mm Pb Eq front, .25mm back. Hook & loop closure with 2" safety buckle. Free pocket and waist pad.

SKIRT-GUARD Female

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
SFR	Regular	\$200.00	\$220.00	\$240.00
SFL	Light Wt.	\$260.00	\$280.00	\$300.00
SFG	LMG Lead-Free	\$310.00	\$330.00	\$350.00

KILT-GUARD Male

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
KMR	Regular	\$200.00	\$220.00	\$240.00
KML	Light Wt.	\$260.00	\$280.00	\$300.00
KMG	LMG Lead-Free	\$310.00	\$330.00	\$350.00



Max-Wrap

Wrap around, full-body protection coat-style garment
Distributes weight between shoulders and hips. 0.5mm
Pb equivalency protection front, 0.25mm back. Hook
& loop closure with flexible back is fast and easy.
Generous armholes for EZ-fit. Free dual pocket and
shoulder pads.

Innovative Max Design!

Solid front and back panels provide no overlap
means less bulk, better fit, better freedom of move-
ment, secure closure at waist and shoulder, flexible
back support panels and lighter weight.

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
MWL	Light Wt.	\$400.00	\$420.00	\$440.00
MWG	LMG Lead-Free	\$450.00	\$470.00	\$490.00



Max-Guard Plus

Overlap wrap around offers full-body protection
Distributes weight between shoulders and hips. 0.5mm
Pb equivalency protection front, 0.25mm back. Hook
& loop closure with inside elastic Back-Guard. Front
belt assures snug fit and maximum support. Generous
armholes for EZ-fit. Free dual pocket and shoulder pads.

cat. #	lead type	petite/ sm/ med	lg / xl	2X / 3X
MXL	Light Wt.	\$400.00	\$420.00	\$440.00
MXG	LMG Lead-Free	\$450.00	\$470.00	\$490.00

Inside:
built-in 6"
elastic back
support.

See Best-Fit
Sizing Chart
on page 22.

Thyro-Guards



Standard

Binding all around means affordable comfort and flexibility.

cat. #	lead type	price
TYR	Regular	\$52.00
TYL	Light Wt.	\$61.00
TYG	LMG Lead-Free	\$72.00



No Binding

No binding means extra comfort, flexibility and lighter weight.

cat. #	lead type	price
TNR	Regular	\$67.00
TNL	Light Wt.	\$76.00
TNG	LMG Lead-Free	\$87.00



Mid-Binding

No binding at the top means extra comfort under the chin.

cat. #	lead type	price
TVR	Regular	\$60.00
TVL	Light Wt.	\$69.00
TVG	LMG Lead-Free	\$80.00

Attached Standard

cat. #	lead type	price
TAR	Regular	\$64.00
TAL	Light Wt.	\$73.00
TAG	LMG Lead-Free	\$84.00

Attached No Binding

cat. #	lead type	price
TCR	Regular	\$80.00
TCL	Light Wt.	\$89.00
TCG	LMG Lead-Free	\$100.00

Attached Mid-Binding

cat. #	lead type	price
TWR	Regular	\$73.00
TWL	Light Wt.	\$82.00
TWG	LMG Lead-Free	\$93.00



(Berry with Pink ribbon only)

Mammography

Shorty patient collar offers protection without compromising breast imaging.

cat. #	lead type	price
TMR	Regular	\$64.00
TML	Light Wt.	\$73.00
TMG	LMG Lead-Free	\$84.00



Pediatric

Small patient collar offers protection during pediatric imaging.

cat. #	lead type	price
TPR	Regular	\$47.00
TPL	Light Wt.	\$56.00
TPG	LMG Lead-Free	\$67.00

Attached Pediatric

cat. #	lead type	price
TTR	Regular	\$59.00
TTL	Light Wt.	\$68.00
TTG	LMG Lead-Free	\$79.00



Buckle Closure

If you prefer buckle closure to regular hook & loop, add "-B" to any of the collar numbers above and add \$5.00 to the price. (Ex.: TYR-B would be \$57.00)

Unique 2-pockets in one!



All Attached Thyroid Collars include: Snap ring for easy collar attachment. Attach anytime with tethered strap. Unique 2-in-1 pocket design: store collar in one pocket and use other pocket for accessories.

Standard, No-Binding, and Mid-Binding collars are 6" tall x 6½" wide x 24" strap.

Select COLOR and FABRIC on page 25. Default color is Royal Blue except for Mammo Collar (Berry only).



CARE-GUARD Front Only

cat. #	lead type	regular	tall/full
CGL	Light Wt.	\$470.00	\$500.00
CGG	LMG Lead-Free	\$520.00	\$550.00

Care-Guard

Pregnancy Apron

Special Protection for Fetal Care

Front Details: Front protection with flexible back straps. 1.0mm Pb equivalency over fetal area, 0.5mm over remaining area. 22" x 24" pregnancy panel can be removed at completion of pregnancy to create a standard Flex-Guard Apron with 0.5mm Pb equivalency protection. Hook & loop closure is fast and easy. Free dual pocket and shoulder pads.

Back Details: Complements the Front Apron for full protection. 0.25mm Pb equivalency. 1" buckle closure (long enough to tie) in front.

CARE-GUARD Back Only

cat. #	lead type	regular	tall/full	size	(width x length)
CBL	Light Wt.	\$160.00	\$180.00	Tall	(24" x 40")
CBG	LMG Lead-Free	\$210.00	\$230.00	Full	(28" x 38")



Breast-Guard

Female Breast Protection

Allows full viewing of anatomy EXCEPT sensitive breast tissue. Hook and Loop strap drapes around the neck. Adjust shields anywhere along the strap. 0.5mm Pb equiv. protection.

Two convenient sizes:

Adult: 6" wide x 7" high; **Child:** 4" wide x 4½" high.

Set = two shields + one strap

cat. #	lead type	price/set	size
BSR	Regular	\$60.00	A (Adult)
BSL	Light Wt.	\$70.00	C (Child)
BSG	LMG Lead-Free	\$80.00	



Pregnancy Semi & Lap-Guards

Half Apronette

Similar to Semi-Guard and Lap-Guard shown at the top of page 11, but with 1.0mm Pb equivalency protection and a 2" belt on the Lap-Guard version. Excellent protection for pregnant patients requiring x-ray.

1.0mm PREGNANCY SEMI-GUARD & LAP-GUARD

SEMI cat. #	LAP cat. #	lead type	medium	large	x-large	2X
SPL	LPL	Light Wt.	\$140.00	\$170.00	\$200.00	\$290.00
SPG	LPG	LMG Lead-Free	\$210.00	\$240.00	\$270.00	\$360.00

COLORS & FABRICS: page 25.
Default Fabric: **Royal Blue nylon.**
2-Year Limited Warranty.

Semi-Guard

Extra Coverage Half Apronette

Convenient hook & loop closure. 0.5mm Pb equivalency protection. Excellent protection for patient and staff.



Individual default apronette color is Royal Blue, unless specified.

size	(width x length)
S	(12" x 12")
M	(18" x 16")
L	(24" x 20")
X	(24" x 24")
2X	(36" x 24")
Half	(24" x 12")
Kids*	(8" x 10")

*See Kids Lap-Guard p.13

Lap-Guard

Buckle Closure Half Apronette

Easy-to-use 1" buckle belt closure. 0.5mm Pb equivalency protection. Excellent protection for patient and staff.



SEMI-GUARD & LAP-GUARD

SEMI cat. #	LAP cat. #	lead type	small	med/half	large	x-large	2X
SGR	LGR	Regular	\$60.00	\$80.00	\$100.00	\$120.00	\$180.00
SGL	LGL	Light Wt.	\$80.00	\$100.00	\$120.00	\$140.00	\$200.00
SGG	LGG	LMG Lead-Free	\$150.00	\$170.00	\$190.00	\$210.00	\$270.00

SEMI-GUARD Set & LAP-GUARD Set

SEMI cat. #	LAP cat. #	lead type	price/set
STR	LTR	Regular	\$400.00
STL	LTL	Light Wt.	\$470.00
STG	LTG	LMG Lead-Free	\$740.00



Semi-Guard & Lap-Guard SETS include a 2-Peg Rack and the following apronettes:

small (yellow)	large (royal blue)
medium (berry)	extra-large (green)

Demi-Guard

Half Apron

Child to Adult loin protection. 0.5mm Pb equivalency front protection. Available with 1" buckle closure or convenient Waist Ring.



DEMI-RACK (only)
Cat.# DDD \$ 215.00

DEMI-GUARD with Waist Rings

cat. #	lead type	small	medium	large	x-large
DRR	Regular	\$120.00	\$130.00	\$140.00	\$150.00
DRL	Light Wt.	\$140.00	\$150.00	\$160.00	\$170.00
DRG	LMG Lead-Free	\$170.00	\$180.00	\$190.00	\$200.00

DEMI-GUARD with Buckles or Hook & Loop

Buckle cat. #	H&L cat. #	lead type	small	medium	large	x-large
DBR	DVR	Regular	\$60.00	\$70.00	\$80.00	\$90.00
DBL	DVL	Light Wt.	\$80.00	\$90.00	\$100.00	\$110.00
DBG	DVG	LMG Lead-Free	\$110.00	\$120.00	\$130.00	\$140.00

DEMI-GUARD (Ring set of 4 w/ rack)

cat. #	lead type	price	size	(width x length)
DSR	Regular	\$720.00	S	(8" x 10")
DSL	Light Wt.	\$790.00	M	(12" x 12")
DSG	LMG Lead-Free	\$910.00	L	(16" x 14")
			XL	(20" x 16")

DEMI-GUARD (Buckle or Hook & Loop Set of 4 w/ Rack)

Buckle cat. #	H&L cat. #	lead type	price	size	(width x length)
DGR	DYR	Regular	\$490.00	S	(8" x 10")
DGL	DYL	Light Wt.	\$560.00	M	(12" x 12")
DGG	DYG	LMG Lead-Free	\$680.00	L	(16" x 14")
				XL	(20" x 16")

Gonadal-Guard

Adult Gonadal Shields

Attached belt with 1" buckle. 0.5mm Pb equivalency protection. Available in three sizes.



GONADAL-GUARD (individual)

cat. #	lead type	small	medium	large
GNR	Regular	\$35.00	\$45.00	\$55.00
GNL	Light Wt.	\$40.00	\$50.00	\$60.00
GNG	LMG Lead-Free	\$45.00	\$55.00	\$65.00

GONADAL-GUARD (set of all 3)

cat. #	lead type	price	size	(width x length)
GSR	Regular	\$110.00	S	(4" x 4.5")
GSL	Light Wt.	\$140.00	M	(6" x 7")
GSG	LMG Lead-Free	\$160.00	L	(7.5" x 9")

Diaper-Guard

Adult Protective Diaper

Specifically designed radiation protective diapers for adults. 0.5mm Pb equivalency protection. Convenient tie strings at each hip.

#GDR

\$105.00



Gonadal-Guard

Pediatric Gonadal Shields

Specifically designed sets will accommodate neonatal up to small children. Gender appropriate gonadal shielding. 0.5mm Pb equivalency protection.



MALE (Blue/set of 4)

cat. #	lead type	price/set
GMR4	Regular	\$140.00
GML4	Light Wt.	\$180.00
GMG4	LMG Lead-Free	\$220.00



FEMALE (Pink/set of 5)

cat. #	lead type	price/set
GFR5	Regular	\$180.00
GFL5	Light Wt.	\$240.00
GFG5	LMG Lead-Free	\$300.00

Diaper-Guard

Pediatric Protective Diaper

Specifically designed radiation protective diapers accommodate infants up to small children. Color-coded regular lead diapers: small (red), medium (green), large (yellow). 0.5mm Pb equivalency protection.



cat. #	back x length x front	size	price
GBR-S	9.5"x 11"x 5"	small	\$75.00
GBR-M	12.5"x 12"x 7.5"	medium	\$85.00
GBR-L	15.5"x 13"x 9"	large	\$95.00
GBS-3	set of all 3	set of 3	\$245.00

1-Front Up

2-Flap Over

3-Final Flap



Child-Guard

Front Protection for Kids

EZ Fit, Front Protection

Front protection apron with buckle closure. 0.5mm Pb equiv. front protection. Shoulder Pads: concealed for comfort and durability. Ideal for our younger patients, creating safer and more secure protection than lead drapes or shields.

CHILD-GUARD

cat. #	lead type	child
KGR	Regular	\$180.00
KGL	Light Wt.	\$240.00
KGG	LMG Lead-Free	\$290.00

size (width x length)
Child (18" x 24")



Kids Lap/Semi-Guard

Half Apronette

Excellent protection for small patients. 0.5mm Pb equivalency protection. LAP uses 1" buckle closure and SEMI uses hook & loop.



KIDS LAP/SEMI-GUARD

LAP cat. #	SEMI cat. #	lead type	child
LKR	SKR	Regular	\$50.00
LKL	SKL	Light Wt.	\$60.00
LKG	SKG	LMG Lead-Free	\$90.00



Kiddie Kovers

Kid-Friendly Protection

What makes Kiddie Kovers Unique:


Decrease anxiety and fear with a "hug" from a friendly buddy. All purpose lead blanket protects children from radiation. 0.5mm Pb equivalency protection (regular lead only). Protective panel is 18" x 18" and safe for all modalities.

cat. #	description	price
KK-LBUG	Ladybug	\$170.00
KK-BBEE	Bumble Bee	\$170.00
KK-TURT	Turtle	\$170.00
KK-BFLY	Butterfly	\$170.00
KK-PUPP	Puppy	\$170.00
KK-FROG	Frog	\$170.00
KK-CATT	Cat	\$170.00
KK-MNKY	Monkey	\$170.00

Ladybug



Bumble Bee



Turtle



Butterfly




Puppy



Frog



Cat



Monkey



Hand-Guard

Five fingered protection for full dexterity. Durable, soft navy vinyl covering (15"). 0.5mm Pb equiv. dip-molded core.



HAND-GUARD 5-Finger

cat. #	hand	price
HGP	Pair	\$270.00
HGR	Right Hand	\$152.00
HGL	Left Hand	\$152.00

Palm-Guard

Provides hand protection with option of exposing fingers for dexterity. Durable, soft navy vinyl covering (15"). 0.5mm Pb equiv. dip-molded core.



PALM-GUARD Slitted Mittens

cat. #	hand	price
PSP	Pair	\$365.00
PSR	Right Hand	\$195.00
PSL	Left Hand	\$195.00

Finger-Guard

Radiation Reducing Gloves

Finger-Guard offers sterile hand protection from secondary radiation.



#PLF offers superior tactile sensitivity with excellent attenuation.



#PFC offers excellent tactile sensitivity with good attenuation.



#PRL offers reasonable tactile sensitivity with the best attenuation.



PLF: Lead-Free/Powder-Free/ Latex-Free

PFC: Lead-Free/Latex-Free

PRL: Lead-Free/Powder-Free/
contains Natural Latex/beaded
cuff/ textured finger tips.

cat. #	description	sizes	qty/pkg	price
PLF	.18mm Glove	6½, 7, 7½, 8, 8½, 9	5 Pair	\$365.00
PFC	.23mm Glove	6½, 7, 7½, 8, 8½, 9	5 Pair	\$365.00
PRL	.35mm Glove	6½, 7, 7½, 8, 8½, 9	5 Pair	\$249.00

Attenuation for Finger-Guard - Broad Beam

Description	Cat.#	60 KvP	80 KvP	100 KvP	120 KvP
.18mm Glove	PLF	56.5%	47.1%	41.6%	37.7%
.23mm Glove	PFC	48.0%	40.0%	32.0%	28.0%
.35mm Glove	PRL	60.0%	53.0%	46.0%	41.0%

("mm" is actual thickness of glove in millimeters, NOT lead equivalency.)
(Test results in accordance with ASTM Test Method F2547-06.)

Omni-Guard

Protective Flexible Shield

Used over a patient or equipment, protects from radiation exposure. Used underneath a patient, prevents scatter radiation. All purpose covering for patients or equipment. 0.5mm Pb equivalency protection. Choose any color fabric.



Custom Sizes Available

COLORS & FABRICS: page 25.
Default Fabric: **Royal Blue nylon.**
2-Year Limited Warranty.



cat. #	lead type	small	medium	large	x-large	size	(width x length)
OMR	Regular	\$115.00	\$155.00	\$180.00	205.00	S	(18" x 18")
OML	Light Wt.	\$160.00	\$200.00	\$255.00	310.00	M	(24" x 24")
OMG	LMG Lead-Free	\$250.00	\$315.00	\$390.00	550.00	L	(24" x 36")
						X	(24" x 48")



ADD GROMMETS: if you would like to hang your Omni-Guard as a curtain, you can add grommets at the time of ordering. Let us know your desired location and the distance between grommets. Maximum width: 48". State on your Omni-Guard order: Cat.#OM-GR at \$25.00 for up to five (5) grommets.

Belt-Guard

Back Support Belt

Regular #BGB
\$45.00

X-Large #BGB-2
\$55.00



Helps reduce back stress.
Can be used with any apron.
6" elastic back support belt with 2" buckle closure.
The same belt as used on the Ortho-Guard apron.

Guard-Tote

Apparel Tote Bag

#ACB-28
\$99.00



Convenient tote for all your protective wear. Durable rip-stop nylon fabric with zippered closure. 28" long to minimize folding your aprons. EZ-Carry handle and shoulder strap included.

Hat-Guard

Specify size:
Small - 21"
Medium - 22"
Large - 23 1/2"



Protection for the top of head. **HAT-GUARD**

Hook & Loop fitted closure.
0.50mm Pb equiv. sides,
0.25mm Pb equiv. top.

cat. #	lead type	price
HCL	Light Wt.	\$140.00
HCG	LMG Lead-Free	\$160.00

PROTECTIVE APPAREL

radiation protection for your eyes



Model#	Color
LER-9SPL	Silver
LER-9KPL	Black
LER-9TPL	Taupe

Plano Lens \$195.00
(Prescription Lens available)



Ultra-Guard

Ultra-light and super stylish.

Curved lens adds to lateral protection.
Regular hinges with wrap around frame.
Universal bridge support. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.

Weight: 60 grams.

Frame: impact resistant polycarbonate.



Model#	Color
LER-ACPL	Carbon/Gray
LER-APPL	Purple/Gray
LER-AYPL	Yellow/Blue

Plano Lens \$195.00
(Prescription Lens available)



Avant-Guard

Two-tone colorful Ultra-Guards.

Curved lens adds to lateral protection.
Regular hinges with wrap around frame.
Universal bridge support. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.

Weight: 60 grams.

Frame: impact resistant polycarbonate.



Model#	Color
LER-CBPL	Blue
LER-CRPL	Red
LER-CGPL	Green
LER-CPPL	Purple
LER-CKPL	Black
LER-COPL	Orange

Plano Lens \$195.00
(Prescription NOT available)

Color-Guard

Bright-tone colorful Ultra-Guards.

Curved lens adds to lateral protection.
Regular hinges with wrap around frame.
Universal bridge support. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.

Weight: 60 grams.

Frame: impact resistant polycarbonate.

PROTECTIVE APPAREL

radiation protection for your eyes



Model# Color
LER-DKPL Black

Plano Lens \$195.00
(Prescription Lens NOT available)

Predator-Guard

Radiation protection and splash protection.
Close fit provide excellent visibility. Extra soft nose pads for superior comfort. Wrapping lens area increases protection. For larger framed faces. Includes premium case.

Lead Equiv.: 0.75mm lens.
Weight: 70 grams.
Frame: impact resistant polycarbonate.



Model# Color
LER-TRPL Red
LER-TKPL Black
LER-TSPL Silver

Plano Lens \$195.00
(Prescription Lens available)

Turbo-Guard

For that sharp and stream-lined look.
Large lens area for excellent coverage. Rubberized temple bars for comfortable grip. Rubber, ultra-soft, fixed nose pad. Available in almost all prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.
Weight: 65 grams.
Frame: impact resistant nylon.



Model# Color
LER-6SPL Silver
LER-6KPL Black
LER-6RPL Red
LER-6BPL Blue

Plano Lens \$215.00
(Prescription Lens available)

Sure-Guard

Stylish design for medium faces.
Large lens area for excellent coverage. Regular hinges with wrap around frame. Rubber, ultra-soft, fixed nose piece. Available in almost all prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.
Weight: 63 grams.
Frame: impact resistant plastic.





Model# Color
 LES-53PLT Tortoise
 LES-53PLC Crimson
 LES-53PLB Blue
 Plano Lens \$215.00
 (Prescription Lens available)



Wrap-Around

The “original” with great side protection. Generous lens area increases protection. Spring hinges and adjustable nose. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.
 Frame: impact resistant plastic.



Model# Color
 LER-VFPL Fancy Flight
 LER-VVPL Cherry Vines
 LER-VSPL Rad Skull
 Plano Lens \$215.00
 (Prescription Lens available)



Viva-Guard

Lively designs for small to medium faces. Smaller frame for excellent fit of smaller faces. Regular hinges with wrap around frame. Ultra soft bridge support. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.
 Weight: 60 grams.
 Frame: impact resistant polycarbonate.



Model# Color
 LER-PBPL Blue/Orange
 LER-PSPL Silver/Pink
 LER-PKPL Black/Gray
 Plano Lens \$215.00
 (Prescription Lens available)



Power-Guard

Colortone accents increase comfort. Large lens area for excellent coverage. Soft colored endpiece for comfort and fit. Padded bridge, temples, and endpiece. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm lens.
 Weight: 64 grams.
 Frame: impact resistant polycarbonate.

PROTECTIVE APPAREL

radiation protection for your eyes



Model# Color
LES-AFYK Black/Yellow
LES-AFPW Purple/White
LES-AFGP Gray/Pink
Plano Lens \$188.00
(Prescription Lens available)



Astro-Guard

Side Shield protection.

Attractive and sturdy women's frames. Flexible temples fit almost any face. Excellent visibility. Available in most prescriptions. Includes premium case.

Lead Equiv.: 0.75mm front and sides.
Weight: 66 grams.
Frame: impact resistant plastic.



Model# Color
LES-CKPL Black
LES-CWPL White
LES-CTPL Tortoise
Plano Lens \$150.00



Cover-Guard

Comfortable protection over regular glasses.

Excellent front and side visibility. Provides limited splash protection. Can be worn over prescription eyewear. Regular hinges with molded nose piece. Includes 0.75mm Pb glass side shields. Includes premium case.

Lead Equiv.: 0.75mm front and sides.
Weight: 85 grams.
Frame: black rubberized nylon.



Model# Color
LES-RKPL Black
Plano Lens \$175.00
(Prescription Lens available)

Secure-Guard

Vented goggle for extreme protection.

Adjustable headband keeps goggles securely on head and prevents slippage. Neoprene foam cushioning for facial comfort. Large lens pocket allows for maximum viewing area. Vented to prevent fogging.

Lead Equiv.: 0.75mm lens.
Weight: 66 grams.
Frame: impact resistant plastic.



Aprons not included.

Mobile Apron Rack

Ten (10) swing-arm mobile garment stand. Durable all-steel construction. 24" x 24" H-Base for stability and easy cleaning. E-Z-Roll swivel locking casters. Convenient handle for easy mobility. Neutral Almond baked enamel finish. Overall height is 55". Some assembly required. Ships @ 110 lb.



Standard Mobile Rack
#ARM-100
\$950.00

Standard WITH Glove Rack
#ARM-G60
\$1020.00

**Easy-Clean
H-Base!**



Arms swing inward to conserve space.

Mobile Apron Z-Rack

Z-Base mobile rack nests with others when not in use. Heavy duty casters on 12G steel base. Chromed-steel posts and hang bar. 60"long x 60"tall bar allows hanging without creasing. Ships @ 35 lbs.

Rack ONLY
#ARM-ZNH
\$280.00

Rack with
Ten Hangers
#ARM-Z10
\$590.00



Apron Hangers

Designed for standard aprons, panoramic, vests, skirts, etc. Top hook swivels for easy hanging. 21" long chromed-steel with anti-slip tips. Glove holders mounted onto back panel. Use on ARM-ZNH rack, peg rack, etc. Single ships @ 2 lbs.

Single Hanger
#ARH-01
\$28.00

Set of 10
#ARH-10
\$265.00



Wall Mounted Apron Racks

Ten or five swing-arm wall-mount apron racks. Mounting holes on 16" centers. All-steel durability with almond enamel finish. 10-Arm Ships @ 25 lb. / 5-Arm Ships @ 17 lb.



Glove Adapter

Retrofits to Mobile Rack.

Holds six (6) 18" gloves.

Ships @ 8 lb.

#ARG-06

\$80.00



10-ARM #ARW-10 \$360.00



5-ARM, SWING LEFT #ARW-L5 \$210.00



5-ARM, SWING RIGHT #ARW-R5 \$210.00

Apron Peg Racks

6" angled pegs keep garments from slipping. All-steel Almond enamel finish.



6-peg
(for 3 aprons)

#ARP-06

\$120.00

Ships @ 10 lb



6-peg and 8-peg racks: mounting holes on 16" and 24" centers.

2-peg
(for 1 apron)

#ARP-02

\$57.00

Ships @ 4 lb



8-peg
(for 4 aprons)

#ARP-08

\$145.00

Ships @ 11 lb



Apron/Glove Racks

Designed for one apron, but will hold more. All steel, heavy duty construction. 16" & 24" holes match most wall studs. Glove holders mounted onto back panel. 28.5" long x 1.75" diameter bar prevents creasing. Ships @ 13 lbs.



Without

Glove Holders

#AGR-105

\$155.00



With Two

Glove Holders

#AGR-110

\$205.00



Aprons not included.

Best-Fit Sizing Charts

Front Protection Apron				Vest (male) & Full Wrap Around			
size	front width	length	chest circumference	size	chest/tummy	vest length	full wrap length
Petite	20"	32"	< 36"	Petite	35"-36"	19.5"	32"
Small	21"	34"	< 40"	Small	37"-38"	20"	34"
Medium	22"	36"	< 44"	Medium	39"-41"	21"	36"
Large	24"	38"	< 48"	Large	42"-44"	22"	38"
X-Large	26"	40"	< 52"	X-Large	45"-49"	23"	40"
2X	28"	40"	< 56"	2X	50"-54"	24"	40"
3X	30"	40"	< 60"	3X	55"-59"	24"	40"

Sizes shown for apparel other than Front Protection are actual body measurements.


All lengths are from the sternal notch or waist downward, depending on the garment.

Vest (female)					Skirt / Kilt			
size	bust	waist	hip	length	size	waist	hip	length
Petite	33"-34"	< 31"	35"-36"	19"	Petite	25"-27"	32"-34"	19"
Small	35"-36"	< 33"	37"-38"	19.5"	Small	28"-29"	35"-36"	19"
Medium	37"-38"	< 35"	39"-40"	20"	Medium	30"-33"	37"-40"	20"
Large	39"-41"	< 39"	41"-43"	21"	Large	34"-37"	41"-44"	21"
X-Large	42"-44"	< 41"	44"-46"	22"	X-Large	38"-41"	45"-48"	22"
2X	45"-49"	< 46"	47"-51"	23"	2X	42"-45"	49"-52"	22"
3X	50"-54"	< 51"	52"-56"	23"	3X	46"-50"	53"-57"	22"

Matching the Apron to You

Techno-Aide Sizing Charts are based on averages. But let’s face it... nobody is average! So, in order to get the Best Fit for You, please consider the following when measuring yourself for either a Standard or a Custom garment:

- PROTECTION IS PRIMARY: Why are you wearing an apron? Protection from radiation! Some companies promise lighter weight garments, but their garments just don’t measure up to proper specs. At Techno-Aide, we’re much more concerned about your health and safety. Comfort is important, and while our garments are comfortable, comfort will always take a back seat to protection. Isn’t that why you’re purchasing Protective Apparel?
- CIRCUMFERENCE CONSIDERATIONS: when measuring for a Front Protection Apron, measure from one side of the widest part of your body, around to the other side - the actual distance you want protected. Then, tell us if that measurement is your chest, tummy, or hips.
- HEIGHT CONSIDERATIONS: in relation to your circumference, do you have a long or short torso? We want to make sure that you are adequately protected!



Getting the Best-Fit

If you don't think that our standard sized garments will fit you properly, we are pleased to offer you the chance to create your own size, to fit your specific body measurements. Just complete the form below and we will create your customized protective wear. Please note, that due to the nature of customized apparel, these garments cannot be returned. Please be accurate. Thanks!

NAME _____ FACILITY _____

ADDRESS 1 _____ ADDRESS 2 _____

CITY _____ STATE _____ ZIP _____

DAYTIME PHONE _____ EMAIL _____

APRON CAT.# _____ DESCRIPTION _____ LEAD TYPE _____

FABRIC (p25) _____ BINDING COLOR (p24) _____


EMBROIDERY (OPTIONAL) LINE 1, UP TO 24-LETTERS (FREE) _____

EMBROIDERY LINE 2, UP TO 24-LETTERS (\$10.00) _____

EMBROIDERY THREAD COLOR (p24) _____ LETTER STYLE (p24) _____

ADD ATTACHED THYROID COLLAR? YES NO CAT.# _____ (COLOR WILL MATCH APRON)

Getting Accurate Measurements (all in inches)



BUST/CHEST: _____
Circumference around the widest part of your chest.

SIDE-to-SIDE: _____
The widest part of your body that you want covered by a Front Protection Apron. Please specify the this area (circle): chest / tummy / hips

WAIST: _____
Standard Apron: measure your waistline.
Skirt/Kilt: measure where you want the skirt/kilt top (may be higher than your regular waistline).

HIP: _____
Circumference around the widest part of your hips.

LENGTH:
FRONT APRON or WRAP-AROUND: _____
Sternal notch to just above the knee.

VEST/TORSO: _____
Sternal Notch to desired length below your waist.

SKIRT/KILT: _____
Desired skirt/kilt top (may be above your regular waistline) to just above your knee.

Please photocopy this page to submit your order.

Embroidery Sewn-In

Selected garments come with first line FREE embroidery! Add flexibility with an embroidered patch. Please specify what you want embroidered and which type style you prefer (see below). Embroidered items are non-returnable, non-refundable.

First Line up to 24-Letters: \$FREE!



Select a Style:

ALL CAPS
(default)

Block

Script

Additional lines up to 24-letters: \$10.00

Embroidery Add-On Patch

This flexible option has your information embroidered onto a patch. Hook & Loop on the back of the patch and on the garment secure the patch in place.

#EMP (Patch) \$18.00



(Patch above pocket)



Accent Colors

Optional BINDING Colors:

Black (default)	Red
Lt. Gray (default)	Yellow
Navy	Orange
Royal Blue	Tan
NEW! →	Pink

(The black or gray stripe on each fabric square is the default binding color for that fabric.)

Optional EMBROIDERY Colors:

<i>White (default)</i>	<i>Red</i>
<i>Black (default)</i>	<i>Yellow</i>
<i>Lt. Gray</i>	<i>Orange</i>
<i>Charcoal</i>	<i>Tan</i>
<i>Royal Blue</i>	<i>Brown</i>

Things to Remember

CARE AND MAINTENANCE:

Always hang your garment on large diameter hangers. Do not wad or fold your garment, as that may weaken the level of protection on the folded seam. Do not store sharp objects in the pockets or store the garment on top of a sharp object.

WARRANTY:

Two-Year Warranty on materials and workmanship.

Shoulder Strap (FREE Option!)



#EMS

\$-No Charge!

Simply ask for it with your order. Hook & loop Strap is attached to one shoulder and reaches across to the other. Helps keep shoulders high and comfortable.



Recycle

Techno-Aide will gladly accept your old, leaded Techno-Aide Brand aprons for recycling. Call for a "Recycle Return Authorization." There will be no charge to you; we only ask that you pay for your own shipping. "Lead-Free" apparel is environmentally friendly and does not need to be recycled.

PROTECTIVE APPAREL

fabric / color selections

Nylon Solids

Liquid Repellent
Light & Supple



11-Royal Blue



12-Navy



20-Green



21-Teal



22-Forest



23-Slate



33-Black



41-Sunshine



43-Orange



51-Red



52-Hot Pink



53-Purple



54-Berry



55-Burgundy



56-Brown

Colorful Prints

Liquid Repellent
Medium Weight



70-Blue Plaid



71-Leopard



72-Lots•O•Dots



73-Sleepy Sheep



77-Life Vines



78-Waterfall



83-Cool Waves



82-Pink Camo



86-Camo Night



85-Digital Camo



94-Big Dipper



95-Andromeda



74-Playful Pups



76-Blue Paws



80-Red Paws

(Prints are shown **SMALLER** than actual size)

Diamond Ripstop

Medium Weight Nylon
Filament Reinforced



61-Sapphire



62-Ruby



63-Emerald



64-Onyx



65-Amethyst

Feather Ripstop

Light Weight Nylon
Filament Reinforced



R701-Royal Blue



R702-Red



R703-Forest



R704-Black



R705-Chocolate



R706-Navy



R707-Burgundy



R708-Kelly



R709-Charcoal



R710-Gray

Natura

Smooth, Soft, Flexible
with a Natural Feel



N911-Stratos



N921-Ocean



N933-Midnight

- "White" and "Black" denote the default embroidery color.
- The Black and Gray lines denote the default binding trim color.

TECHNO-SHIELD MOBILE BARRIERS

Tilted Panel Barrier



This ergonomically unique barrier features a tilted leaded viewing panel that allows the technologist or physician to be closer to a lower placed patient position.

NEW!



#PBMB-TILT
(Call for Price)
Leaded Glass: 30"w x 24"h
Overall Size: 34"w x 70"h
Ships @ 450 lbs.

Viewing Panel Barriers

Bring a full **2.0mm Pb Equivalency** of radiation protection anywhere you need it with our x-ray protective mobile barriers.

Each barrier features transparent/clear upper section of x-ray protective leaded glass; Opaque lower section of solid, leaded x-ray protection; Locking caster system for easy mobility (it is recommended that two people move the barrier for safety); Crated for secure shipment.



#PBMB-3024
(Call for Price)
Leaded Glass: 30"w x 24"h
Opaque Panel: 30"w x 48"h
Overall Size: 33"w x 75"h
Ships @ 450 lbs.




#PBMB-3060
(Call for Price)
Leaded Glass: 30"w x 60"h
Opaque Panel: 30"w x 12"h
Overall Size: 33"w x 75"h
Ships @ 450 lbs.



#PBMB-3630
(Call for Price)
Leaded Glass: 36"w x 30"h
Opaque Panel: 36"w x 42"h
Overall Size: 39"w x 75"h
Ships @ 450 lbs.

2.0mm Pb Eq STOPPING POWER!

COLORS for ALL Window & Panel Barriers



Neutral
#000N
(Default)




Pink
#3015



Rose
#3014



Rust
#3012



Light Gray
#7035



Azure
#5012



Dark Blue
#5026



Sage
#1000



Teal
#6000



Graphite
#7024

NEW!

Match your facility decor or just brighten up your exam room with your choice of *colorful barriers!

If a color is not specified, our popular NEUTRAL will be the default color.

Order by including the COLOR NAME and the FOUR-DIGIT Code# with your Barrier selection.

*(Barrier color may vary slightly from colors printed on this page.)

Viewing Window Barriers

Bring a full **2.0mm Pb Equivalency** of radiation protection anywhere you need it with our x-ray protective mobile barriers.

Each barrier features a **10" w x 12" h** clear window made of x-ray protective leaded glass; Window center is 60" from the bottom of the panel. Opaque barrier section is made of solid, leaded x-ray protection; Locking caster system for easy mobility (it is recommended that two people move the barrier for safety); Crated for secure shipment. **NOTE:** 84" barriers may not fit through some doorways.



#PBMW-2472

(Call for Price)

Leaded Glass: 10" w x 12" h

Overall Size: 24" w x 72" h

Ships @ 350 lbs.



#PBMW-3072

(Call for Price)

Leaded Glass: 10" w x 12" h

Overall Size: 30" w x 72" h

Ships @ 450 lbs.



#PBMW-3672

(Call for Price)

Leaded Glass: 10" w x 12" h


Overall Size: 36" w x 72" h

Ships @ 450 lbs.

Standard Shield

0.5mm Pb equivalency.

Adjustable height: 36" to 60". Stable, five-legged star base. Default color: royal blue.



0.5mm - SHIELDING

cat. #	media	width x height	price
LPS-S2424	lead	24" x 24"	\$520.00
LPS-S2436	lead	24" x 36"	\$555.00
LPS-S2448	lead	24" x 48"	\$590.00
LPSG-S2424	lead-free	24" x 24"	\$780.00
LPSG-S2436	lead-free	24" x 36"	\$835.00
LPSG-S2448	lead-free	24" x 48"	\$885.00

See Fabric Selection for all Shields on page 25.

Max-Shield

0.5mm Pb equivalency.


Adjustable height: 8" to 60". All steel H-base for maximum stability. Hardened chrome steel post. Posi-Lock height adjustable arm with spring-loaded pin. Default color: royal blue.

0.5mm - SHIELDING

cat. #	media	width x height	price
LPS-H3636	lead	36" x 36"	\$850.00
LPS-H4848	lead	48" x 48"	\$950.00
LPSG-H3636	lead-free	36" x 36"	\$1275.00
LPSG-H4848	lead-free	48" x 48"	\$1425.00

1.0mm - SHIELDING

cat. #	media	width x height	price
LPS-H3636-1	lead	36" x 36"	\$1275.00
LPS-H4848-1	lead	48" x 48"	\$1425.00
LPSG-H3636-1	lead-free	36" x 36"	\$1910.00
LPSG-H4848-1	lead-free	48" x 48"	\$2130.00



Solid Panel Deluxe Shield

0.5mm Pb equivalency. Adjustable height: 8" to 60". Steady all steel T-base. Hardened chrome steel post. Posi-Lock height adjustable arm with spring-loaded pin. Default color: royal blue. Also available in 1.0mm lead and lead-free in both thicknesses.

0.5mm - REGULAR LEAD

cat. #	w x h	price
LPS-T2424S	24" x 24"	\$580.00
LPS-T2430S	24" x 30"	\$595.00
LPS-T2436S	24" x 36"	\$615.00
LPS-T2442S	24" x 42"	\$640.00
LPS-T2448S	24" x 48"	\$650.00
LPS-T3024S	30" x 24"	\$595.00
LPS-T3030S	30" x 30"	\$645.00
LPS-T3036S	30" x 36"	\$675.00
LPS-T3042S	30" x 42"	\$695.00
LPS-T3048S	30" x 48"	\$745.00



Pleated Panel Deluxe Shield

0.5mm Pb equivalency. 24" has four panels, 30" has five. Adjustable height: 8" to 60". Steady all steel T-base. Hardened chrome steel post. Posi-Lock height adjustable arm with spring-loaded pin. Default color: royal blue. Also available in 1.0mm lead and lead-free in both thicknesses.

0.5mm - REGULAR LEAD

cat. #	w x h	price
LPS-T2424P	24" x 24"	\$630.00
LPS-T2430P	24" x 30"	\$645.00
LPS-T2436P	24" x 36"	\$665.00
LPS-T2442P	24" x 42"	\$690.00
LPS-T2448P	24" x 48"	\$700.00
LPS-T3024P	30" x 24"	\$645.00
LPS-T3030P	30" x 30"	\$695.00
LPS-T3036P	30" x 36"	\$725.00
LPS-T3042P	30" x 42"	\$795.00
LPS-T3048P	30" x 48"	\$825.00



Standing Shield

Bring 0.5mm Pb equivalency of radiation protection anywhere you need it.

The 48" tall x 21½" wide x 16" deep unit comes fully assembled on a sturdy H-Base frame of powder coated steel. Grip handle on the back of the unit facilitates easy mobility. Default color: royal blue.

Cat.# LPS-HT-2448
\$1950.00



Table Drape Shielding

Protect staff from scatter radiation without obstructing the C-Arm.

0.5mm Pb equivalency. Attaches with hook & loop (all included).

Order two for double-sided use. Available in three widths: 24", 36" and 48". Mix sizes to fill desired width of coverage. Available in three heights: 24", 30" and 36". Default color: royal blue.



Main Drape

cat. #	panel count	width x height	price
TDP-2424	4-panels	24" x 24"	\$295.00
TDP-2430	4-panels	24" x 30"	\$395.00
TDP-2436	4-panels	24" x 36"	\$495.00
TDP-2442	4-panels	24" x 42"	\$595.00
TDP-3024	5-panels	30" x 24"	\$349.00
TDP-3030	5-panels	30" x 30"	\$429.00
TDP-3036	5-panels	30" x 36"	\$449.00
TDP-3042	5-panels	30" x 42"	\$549.00
TDP-3624	6-panels	36" x 24"	\$395.00
TDP-3630	6-panels	36" x 30"	\$449.00
TDP-3636	6-panels	36" x 36"	\$549.00
TDP-3642	6-panels	36" x 42"	\$649.00
TDP-4824	8-panels	48" x 24"	\$495.00
TDP-4830	8-panels	48" x 30"	\$595.00
TDP-4836	8-panels	48" x 36"	\$695.00
TDP-4842	8-panels	48" x 42"	\$795.00
TDP-CUST	Custom	Choose	Call

Anesthesia Drape

cat. #	panel count	width x height	price
TDP-2024	3-panels	20" x 24"	\$250.00
TDP-2030	3-panels	20" x 30"	\$350.00
TDP-2036	3-panels	20" x 36"	\$450.00
TDP-2042	3-panels	20" x 42"	\$550.00
TDP-2224	4-panels	22" x 24"	\$275.00
TDP-2230	4-panels	22" x 30"	\$375.00
TDP-2236	4-panels	22" x 36"	\$475.00
TDP-2242	4-panels	22" x 42"	\$575.00
TDP-CUST	Custom	Choose	Call

A DRAPE for every need!

Lead Blockers

Maximize efficiency and optimize use of receptors.
Divides receptor for multiple exposures. Covered with blue vinyl laminate. Flexible vinyl-lead positions easily.
Cleans easily with most common disinfectants. Provides 1.0mm Pb Equivalency shielding.



cat. #	blocker size	fits receptor	price
FLB-58	5" x 8" (13x18cm)	8" x 10"	\$35.00
FLB-52	5" x 12" (13x30cm)	10" x 12"	\$45.00
FLB-61	6" x 10" (15x24cm)	10" x 12"	\$45.00
FLB-71	7" x 11" (18x28cm)	11" x 14"	\$60.00
FLB-77	7" x 17" (18x43cm)	14" x 17"	\$85.00
FLS-78	set of all FIVE	Save 10%	\$243.00



Lead-Lined Wall Bin

Protect your investment when plates are not in use.
Solid steel construction with vinyl-lead lining. Provides 1.0mm Pb Equivalency shielding. Specifically designed for receptor plates. Opening is 22"w x 5"d x 14"h.



Lead-Lined with shielded body.

#80LL
\$995.00

Wall Bin overall size:
22"w x 5³/₈"d x 22"h
Shipping weight: 25 lbs.

Lead-Lined Table-Top Bin

Store multiple receptors without worry of scatter radiation. 14-gauge welded steel construction with vinyl-lead lining. Provides 1.0mm Pb Equivalency shielding from scatter. Bumper guards protect receptors from impact damage. Opening is 16"w x 19"d x 20"h.



Lead-Lined with shielded body.

#CPH-LL12
\$1585.00

Unit above with shielded door:
#CPH-LD12
\$2095.00


Mobile Image Receptor Holders

ALL MOBILE UNITS INCLUDE:

Vertical range from bottom of image receptor: floor-to-54". Posi-Lock height adjustable arm with spring-loaded pin. Features 60" chromed steel post. 24" H-Base with locking casters. Broad base prevents tipping. Platinum baked enamel finish. 1½" deep channels accept grids and caps. Units ship in two cartons: 1@15#, 1@45#.



SERIES-B
Tilt & Rotate Models
MCH-33 & MCH-37



SERIES-B
Non-Tilt, Non-Rotate Model
MCH-40




Tilt & Rotate Range of Motion:

Straight Vertical
Incline Forward
Incline Backward
Rotate or Oblique Tilt
Face Up
Face Down
Anywhere in Between

Ordering Guide:

- **SERIES-B (Vertical Clamp)**
 - height limit up to 20"
 - no width limit
- **SERIES-C (Bilateral Clamps)**
 - no height limit
 - width limit 13½" to 20½"
- **WEIGHT LIMITS:**
 - MCH-33/MCH-37 up to 12-lbs.
 - MCH-40 up to 22-lbs.



All models feature
the safe and secure
Posi-Lock Handle

Tabletop Lateral Receptor Holders

Accepts grids and caps (channel is 1¼"d). Upper and lower channels slide on 24" post. Holds receptors up to 24" inches tall. Base plate slides under patient or pad. 35LCH upper channel includes a R & L marker. 40LDH features an upper channel that fits into handle and an extra large base for stability. Upper and lower channels are lined with no-skid vinyl. Shipping weight: 35LCH @ 7 lbs / 40LDH @ 9 lbs.

24" Variable Height
Standard Holder
#35LCH
\$230.00



24" Variable Height
DR-Plate Holder
#40LDH
\$265.00



Mobile Tilt & Rotate Image Receptor Holders

Standard arm length provides maximum stability.

Fully articulating holder for a variety of receptors.

Weight limit for both Series-B and Series-C is 12-lbs.

See additional specs on page 32.



Mobile Non-Tilt, Non-Rotate Holder for Heavy Receptors

Features a fixed position holder (non-articulating) for stability of heavier receptors.

Weight limit for both Series-B and Series-C is 22-lbs.

See additional specs on page 32.



SERIES-B

#MCH-40B

\$1150.00



SERIES-C


#MCH-40C

\$1150.00

Mobile Tilt & Rotate Extending Arm Receptor Holders

Extends image receptors 11" to 23" from post. See additional specs on page 32.


Heavy DR receptors may pose a tipping risk when fully extended. Extended weight limit for both Series-B and Series-C is 12-lbs.



Arm Extended
23" from post



SERIES-B
#MCH-37B
\$1195.00



Arm In & Locked
11" from post



SERIES-C
#MCH-37C
\$1195.00



2-Step Weight Bearing Platform

Allows for weight bearing imaging of the knees, ankles, and feet. Three receptor slots on the top platform and two slots on the bottom platform. 52" high safety rail provides for patient stability. Choose from two platforms: Comfort Top: gray, high-density crosslink closed cell foam for comfort and support. Durability Top: black, high-density, solid polyethylene with anti-slip matting. 350-lb. weight capacity. Platform heights: 11½" and 16¼". All steel frame with 30"l x 20"w base.

Comfort Top
#WMP-16
\$1650.00
Actual Wt: **64 lbs**
Ships@ 105 lbs

Durability Top
#WMP-22
\$2325.00
Actual Wt: **105 lbs**
Ships@ 160 lbs

AP Weight Bearing Receptor Protectors

Provides supportive platform when patient is required to stand on the receptor. Accommodates DR, CR, and analog panels. Standard models support patients up to 425-lbs. Frame is resilient aluminum with a clear polycarbonate support surface. Includes an easy-lift handle.



EXTRA CAPACITY CARBON FIBER PROTECTORS
Perfect for Ortho, Podiatry, and Bariatrics.
800-lb weight capacity.

STANDARD PROTECTORS			
cat. #	receptor size	receptor style	price
WBC-0810W	8" x 10"	cassette or wireless panel	\$595.00
WBC-1012W	10" x 12"	cassette or wireless panel	\$875.00
WBC-1417W	14" x 17"	cassette or wireless panel	\$895.00
WBC-1417T	14" x 17"	tethered DR panel	\$1060.00
WBC-1717T	17" x 17"	tethered DR panel	\$1230.00

EXTRA CAPACITY CARBON FIBER PROTECTORS			
cat. #	receptor size	receptor style	price
WBC-1417CF	14" x 17"	carbon fiber wireless	\$1060.00
WBC-1417CFT	14" x 17"	carbon fiber tethered DR	\$1060.00

3-Step Weight Bearing Platform

Designed for tube stands that don't go low enough. Three receptor slots on the top platform and two slots on the middle platform. 59" high safety rail provides for patient stability. Choose from two platforms: Comfort Top: gray, high-density crosslink closed cell foam for comfort and support. Durability Top: black, high-density, solid polyethylene with anti-slip matting. 350-lb. weight capacity. Platform heights: 8½", 19½", and 23". All steel frame with 42"l x 20"w base.

Comfort Top
#WMP-14
\$2365.00
Actual Wt: **92 lbs**
Ships@ 158 lbs

Durability Top
#WMP-28
\$3145.00
Actual Wt: **133 lbs**
Ships@ 200 lbs



Clear coated crosslink foam comfort top.



Solid HD polyethylene durability top.



Forward and side handrails for patient stability.



Rear wheels for easy tilt & go.



Lower Channel is 1½" wide to accommodate different receptors and Upper Channel is 1½" wide to accommodate GE Flash Pad.

Weight Bearing Receptor Holder

Economical solution for weight bearing examinations. Crosslink closed cell foam is coated for durability and easy cleaning. Overall size: 17½" x 11½" x 4".

Cat.#YCWB
\$99.00



Panel Protectors



Protect your investment with our DR Panel Protectors. Each protector contains three layers of your choice of impact-resistant material or carbon fiber for increased weight capacity. These protectors are designed for both tethered and wireless panels, and also include custom utility ports. Includes rounded corners for safe handling and the patented positive lock mechanism to prevent accidental panel pop-outs. Lightweight design makes for easy carrying and one-handed operation.

cat. #	description	price
WBP-LA1417	handle on long side, plus accessory handle	\$1535.00
WBP-LA1417-CF*	handle on long side, plus accessory handle	\$1775.00
WBP-SA1417	handle on short end, plus accessory handle	\$1535.00
WBP-SA1417-CF*	handle on short end, plus accessory handle	\$1775.00
WBP-S1417	handle on long side	\$1415.00
WBP-S1417-CF*	handle on long side	\$1655.00
WBP-T1417	handle on short end	\$1375.00
WBP-T1417-CF*	handle on short end	\$1625.00

*(Includes **Carbon Fiber** Insert)
NOTE: Not for use with panels with built-in handles.



Mobile Imaging Caddy

Keep everything at your fingertips for mobile imaging. Conveniently hangs from mobile imager mast. Large pouch for Receptor Covers. Medium pouch for accessories. Two small pouches for markers, pens, and other small items. Heavy gauge, reinforced vinyl with clear outer pockets. Caddy is 11" w x 16" h with a strap that extends to 24".

#PMA-C4
\$95.00



Other accessories are sold separately.



No Image Compromise!

Protect Your Investment!



Receptor Covers

Disposable; prevents cross-contamination. Keeps receptors clean and dry from exposure to bodily fluids, solutions, plaster, etc. Specifically designed to fit your imaging receptors. Non-sterile. Latex-Free! 2 mil polyethylene is transparent to x-ray. Can also be used to cover positioning sponges. Available in Standard Overlock and Ziplock™ closures.

Standard Flip-Lip Closure

Allows safe handling and positioning of: CR-Receptors, Cassettes, and DR-Receptors. *Fits 23" x 25" DR-Receptors with Grids!

cat. #	fits receptor size	quantity	price
CCS-08	8" x 10"	1 box of 100	\$65.00
CCS-10	8" x 10"	1 case of 500	\$295.00
CCM-10	10" x 12"	1 box of 100	\$65.00
CCM-12	10" x 12"	1 case of 500	\$295.00
CCL-14	14" x 17"	1 box of 100	\$65.00
CCL-17	14" x 17"	1 case of 500	\$295.00
CCC-DR	20" x 20"	1 box of 100	\$65.00
CCD-DR	20" x 20"	1 case of 500	\$295.00
CCXL-24*	23" x 25"	1 box of 100	\$65.00
CCXL-26*	23" x 25"	1 case of 500	\$295.00

Ships Next Business Day!



Ziplock™ Closure

Disposable bag securely seals in image receptors. Provides Maximum Protection from cross contamination. *Fits 23" x 25" DR-Receptors with Grids!

cat. #	fits receptor size	quantity	price
CCZ-14	14" x 17"	1 box of 100	\$65.00 (2 mil)
CCZ-17	14" x 17"	1 case of 500	\$295.00 (2 mil)
CCXZ-24*	23" x 25"	1 box of 100	\$65.00 (2 mil)
CCXZ-26*	23" x 25"	1 case of 500	\$295.00 (2 mil)
CCH-24*	23" x 25"	1 box of 100	\$80.00 (4 mil)
CCH-26*	23" x 25"	1 case of 500	\$375.00 (4 mil)



Shown:
five 100-pc boxes
inside one 500-pc case

Compact Dispenser
Box of 100 covers

Plastic Bag Dispensing

We realize that some facilities cannot accept products in corrugated cartons. This assortment of Receptor Covers is provided in quantities of 100 covers per dispensing bag to meet those requirements and keep your supply clean and organized. These covers are the very same styles and quality as our originals shown above.

cat. #	fits receptor size	quantity	price
STANDARD BAGS			
CCL-14-B	14" x 17"	1 bag of 100	\$75.00 (2 mil)
CCC-DR-B	20" x 20"	1 bag of 100	\$75.00 (2 mil)
CCXL-24-B*	23" x 25"	1 bag of 100	\$75.00 (2 mil)
ZIPLOCK™ BAGS			
CCZ-14-B	14" x 17"	1 bag of 100	\$75.00 (2 mil)
CCXZ-24-B*	23" x 25"	1 bag of 100	\$75.00 (2 mil)
CCH-24-B*	23" x 25"	1 bag of 100	\$95.00 (4 mil)

Standard Step Stools

All steel 350-lb. capacity. Skid resistant black mat. 9" high stool with 11" x 13 $\frac{3}{4}$ " top; 34" tall with handrail. Handrail is $\frac{7}{8}$ " dia. chrome-plated tubular steel. Easy assembly. Shipping weight: stool 8 lbs.; handrail 4 lbs.



#BSH-65
\$99.00

**Most Popular
Best Selling
Stool!**

Bariatric Step Stools

Bariatric 1000-lb. capacity. Step has anti-slip gripper strips. Flared base legs to prevent tipping. Step Only: 19" long x 15" wide x 8" high. Handrail model height is 36". Shipping weight: stool 11 lbs.; handrail 8 lbs.



#BBH-69
\$415.00

Heavy Duty Step Stools

Bariatric 600-lb. capacity. Steps have a skid resistant black mat surface. Height: Single Step: 8 $\frac{3}{4}$ " high step; 34 $\frac{1}{2}$ " with handrail. Double Step: 8" & 16" high steps; 43" with handrail. Handrail is 1" dia. chrome-plated 16g steel. Shipping weight: Single Step: stool 8 lbs.; handrail 6 lbs. Double Step: stool 20 lbs.; handrail 5 lbs.

#BRF-45
\$17.00



Set of four (4) rubber tips for standard and bariatric step stools.



#BBS-61
\$178.00




#BBD-63
\$375.00



#BBD-66
\$475.00

IV Stands

NEW 5-Leg base poles safety approved 1" and 7/8" telescoping, chromed steel tubes. Height range: 46"-83". Base diameter: economy 22"; 5-Leg 23". 5-Leg bases have scuff guard tips. Two or four rams-horn style IV hooks.



Economy Twist Lock

2-Hook
#BIV-70
\$99.00

Ships @
17 lbs.



5-Leg Base Knob Lock

2-Hook
#BIV-76
\$228.00

4-Hook
#BIV-77
\$249.00

Ships @
25 lbs.

Instrument Stands



Tableside

#BMS-80
\$235.00
Ht: 34"-54"
Ships @ 19 lbs.

Provides ready access to equipment or instruments. Leg posts are 1" chrome plated steel. Steel, easy-slide extension posts. Stainless steel tray is 12" x 19" x 3/4" deep.



Mobile

#BMS-95
\$252.00
Ht: 31"-50"
Ships @ 17 lbs.
in two boxes

Standard Exam Lamp

Flexible gooseneck arm with 360° adjustable steel shade. Chrome steel stand has adjustable height range: 50" to 74". Uses one 60-watt incandescent bulb (not included). Weighted stationery base or 23" diameter mobile base. Weight: 7 lbs.

#BEL-90
\$155.00



Inner shade
keeps outer
shade cool

Halogen Exam Lamp

Flexible gooseneck arm with 540° adjustable head. Beige enamel coated steel stand is adjustable to 65½". Uses one high intensity 35-watt halogen lamp (incl.). 1500 Footcandles @ 24"; 3000° Kelvin. Weighted stationery base or 18" diameter mobile base.

#BEH-93
\$412.00



#BEH-91
\$405.00



Heat sink
technology keeps
light housing cool

Pneumatic Lift Stool with Locking Casters

Pneumatic gas cylinder easily adjusts height with convenient under-cushion ring. Adjustable height: 16" to 21". 16" x 3.5" thick seat cushion. 5-leg, chrome base. Reverse casters auto-lock when seated. 23" base provides excellent stability. Weight: 27 lbs.; Ships @ 33 lbs.

#BPS-30
\$275.00



Pneumatic Lift Stool with Backrest

#BPS-32
\$355.00

(BPS-30 with adjustable
backrest and
paddle height lever)



BPS-30 & BPS-32 ONLY:
choose **VINYL** color
on page 117

Spin-Lift Adjustable Stool 5-Leg

Solid steel, machined thread screw adjusts height by rotation of the stool cushion. Adjustable height: 18" to 22". 16" x 4" thick seat cushion. 5-leg, die cast aluminum base. 23" base provides excellent stability. Weight: 17 lbs.; Ships @ 30 lbs.



#BAS-40
\$325.00



Spin-Lift Adjustable Stool 4-Leg

Solid steel, machined screw adjusts height by rotation of the stool cushion. Adjustable height: 19" to 25½". 14½" x 4" thick seat cushion. 4-leg chrome steel base with footring. Weight: 18 lbs.; Ships @ 30 lbs.

#BAS-30
\$210.00



Pneumatic Lift Stool

Pneumatic gas cylinder easily adjusts height with convenient under-cushion ring. Adjustable height: 19" to 24". 16" x 4" thick seat cushion. 5-leg, die cast aluminum base. 23" base provides excellent stability. Weight: 20 lbs.; Ships @ 30 lbs.

#BPS-20
\$335.00



Pneumatic Lift Stool with Backrest


#BPS-22
\$510.00

(BPS-20 with
adjustable backrest
and D-ring height lever)




Department Signs


All signs are constructed of durable 1/8" Trovicel™, a plastic specific to the sign industry. Double-faced tape is attached to the back for mounting. Large Signs are 8" x 10". Small Signs are 4" x 10". The large MRI "Danger" Signs are 12" x 18".




#PS-08B \$20.00



#PS-08R \$20.00



#PS-04EB \$20.00



#PS-04ER \$20.00



#CXR-13EB \$15.00




#CXR-13ER \$15.00




#XRR-19B \$15.00



#XRR-19R \$15.00



#NCP-24B \$15.00



#NCP-24R \$15.00

Caution / MRI Signs

Signs give critical information to patients & technologists. Enhances safety with a front line written warning. Many states adopt certain signs as law as part of a good quality assurance program. Check your local laws or program requirements.



#CRA-15 \$20.00



#CRA-16 \$20.00



#CRA-17 \$15.00



#CRA-18 \$15.00



#CMM-06 30.00



#MBN-13E \$50.00



#MBN-14S \$50.00



#CMF-01 \$20.00



#MBN-15E \$30.00



#MBN-16S \$30.00

Aquasonic® 100

industry standard ultrasound transmission gel

Recommended for all procedures where a viscous gel is required.

Acoustically correct for the broad range of frequencies used. Aqueous formula will not stain clothing or damage transducers. Hypoallergenic, bacteriostatic, non-sensitizing, non-irritating.

cat. #	size per unit	units/box	price/box
USO-025	250ml bottle	12	\$41.00
USO-100	1 liter bottle	6	\$74.00
USO-500	5 liter jug w/bottle	1	\$30.00
USO-048	20g sterile pouch	48	\$95.00



Aquasonic Clear®

same as Aquasonic® 100, but with no fragrance or color

Recommended for all procedures where a viscous gel is required.

Acoustically correct for the broad range of frequencies used. Aqueous formula will not stain clothing or damage transducers. Hypoallergenic, bacteriostatic, non-sensitizing, non-irritating.

cat. #	size per unit	units/box	price/box
USC-025	250ml bottle	12	\$42.00
USC-100	1 liter bottle	6	\$76.00
USC-500	5 liter jug w/bottle	1	\$31.00



Empty Dispenser Bottles

Standard 250ml bottles can be used for ready refills or as a back-up supply when filled from a larger container. Comes with a black tip cap and can be used in ultrasound gel warmers.

#USB-012
\$45.00/dozen



Ultrasound Probe Covers

Acoustically transparent, minimizing distortion during scanning. Available in Latex or Latex-Free and Sterile or Non-Sterile. A variety of style and sizes for many applications and specialties: Radiology, Vascular, Cardiology, Urology, Phlebology, Emergency Medicine, OB/GYN, etc.

cat. #	description	dia. x length	qty/case	price/cs
LATEX PROBE COVERS, NON-STERILE				
USP-LN01CS	condom size, bulk pack	1.25" x 8"	1200	\$585.00
USP-LN03CS	condom size, individually wrapped	1.25" x 8"	1200	\$926.00
USP-LN04CS	extra-long, bulk	1.18" x 11.8"	600	\$1052.00
LATEX-FREE PROBE COVERS, NON-STERILE				
USP-FN03CS	rolled, individually wrapped	1" x 9"	600	\$915.00
USP-FN04CS	extra long, individually wrapped	1" x 12"	600	\$966.00
LATEX-FREE PROBE COVERS, STERILE				
USP-FS01CS	semi-rolled and tapered, bulk pack	1" d x 9"	300	\$1350.00
USP-FS02CS	semi-rolled, wide abdominal	3.5" flat x 12"	75	\$478.00



Now sold by the CASE!

Elemental Ultrasound Stool


Sit in comfort during extended ultrasound procedures. Back height adjustment. Tilting backrest. Pivot backrest. Tilting seat. Warranty: 7-year limited. Seat: 19"w x 19"d. Back: 19"w x 17.5"h. Seat height range: 23"-33". Overall height: 42"-55". 20" chrome footring. Ships @ 68# (actual wt. 46#)

cat. #	price
ZOES-MX03	\$395.00



*Fabrics
on page
117*

*See page 61
for Ultrasound
positioning
sponges.*



Anti-Fatigue Floor Mat for Sonographers

Provides comfort for standing during long exams. Blown vinyl-nitrile sponge base provides superior cushioning, while a solid vinyl-nitrile surface provides sure-footed protection from wear and tear. This 2' x 3' x 1/2" black mat comes with a 1-year warranty.

#TM-51323
\$75.00



Ultrasound Gel Warmers

fits standard 250ml gel bottles

Steel Gel Warmers

"Uni-Sleeve" heating element adjusts up to 122°F. Illuminated On/Off switch with hospital grade cord. Stainless Steel with white powder coat finish. 5"w/8"w/11"w x 6"d x 5"h.

cat. #	rating	capacity	price/box
UGW-UO1	UL/CSA approved	1	\$175.00
UGW-UO2	UL/CSA approved	2	\$235.00
UGW-UO3	UL/CSA approved	3	\$265.00



Composite Gel Warmers

Maintains constant temperature. Impact-resistant composite housing. Hospital grade power cord. Three Bottle Unit ONLY: - Illuminated On/Off switch. - Energy-saving lid. - 3-position thermostat (UGP-UL3): Low: 97°F / Medium: 102°F / High: 109°F. - Digital display thermostat (UGP-ULD): 1° increments from 97°F to 109°F. Single bottle unit is factory set at 104°F with a heat indicator lamp.

cat. #	rating	capacity	price/box
UGP-UL1	UL/CSA approved	1	\$215.00
UGP-UL3	UL/CSA approved	3	\$340.00
UGP-ULD	UL/CSA approved	3	\$380.00



Counter
Top or
Wall Mount

MRI Furniture

Each piece is non-magnetic and safe to use in the MRI environment.
Some weights do not reflect OverSize UPS/FedEx Rate.



= MR Conditional (no hazards up to 3-Tesla);



= MR Safe (no known hazards);



= MR Unsafe (danger/known hazard).

MRI Step Stool w/ Handrail

Stainless steel construction. 500# weight capacity. Step: 16"w x 12"d x 8½"h. Handrail height: 41". "MRI Green" mat; rubber tips. Ships @ 13 lbs (+oversize) UPS/FedEx.



#MBS-65
\$375.00



MRI Two-Step Platform Stool w/ Handrail

Stainless steel. 400# weight capacity. Top Step: 36"w x 12"d x 16"h, Bottom Step: 36"w x 12"d x 8"h. Handrail height: 41". "MRI Green" mat; rubber tips. Ships @ 48 lbs, motor freight.



#MBD-68
\$1320.00



MRI Chairs & Stools

Screw-post height adjustment. 4" padded seat. Seat upholstered in "MRI Green." Ships via UPS/FedEx.



cat. #	height range	features	base	ship wt.	price/box
CHAIRS					
MBC-30CC	17"-23"	chair only	casters	22	\$830.00
MBC-31CC	17"-23"	chair+back	casters	24	\$1030.00
MBC-30CT	15"-21"	chair only	tips	22	\$760.00
MBC-31CT	15"-21"	chair+back	tips	24	\$910.00

STOOLS

cat. #	height range	features	base	ship wt.	price/box
MBC-30SC	23"-29"	stool only	casters	22	\$830.00
MBC-31SC	23"-29"	stool+back	casters	24	\$1030.00
MBC-30ST	21"-27"	stool only	tips	23	\$800.00
MBC-31ST	21"-27"	stool+back	tips	25	\$920.00

MRI IV Pole

Extends from 48" to 92½". 24" diameter base. 2" swivel casters. Ships via UPS/FedEx: Pole @ 3 lbs, Base @ 7 lbs.

cat. #	description	price
MBI-70	2-hook iv pole	\$520.00
MBI-75	4-hook iv pole	\$540.00



MRI Wheel Chairs

24" diameter wheels, 8" front casters. Burgundy upholstery. Smooth ride, rust-proof hand rims. Retractable wheel locks. Quick release seat belt. Sling style seat (bariatric: solid padded seat). Ships via UPS/FedEx (bariatric: motor freight).



cat. #	seat width	overall	weight cap.	ship wt.	price/box
MBW-18	18"	26"	250#	53	\$2685.00
MBW-20	20"	28"	250#	55	\$3100.00
MBW-22	22"	30"	350#	62	\$3750.00
MBW-24	24"	30"	350#	62	\$4150.00
MBW-24B	24"	35"	550# bari.	67 & 13	\$5650.00
MBW-26B	26"	35"	550# bari.	67 & 13	\$5400.00
MBW-10	I.V.Pole & O ₂ E-tank mount			7	\$240.00



Bariatric
Wheel Chairs



"MR" is on the back
of every chair.

Standard MRI
Wheel Chairs



MRI PVC Transport Gurney

Economic PVC construction. Three weight capacities. Features flip-down side rails. Three position backrest. 79"l x 28"w x 33½". Four locking casters, two or four stationary. Either ships @ 150 lbs, motor freight.



350# CAPACITY
#MBG-3PV
\$1955.00

600# CAPACITY
#MBG-6PV
\$1595.00



MRI Transport Gurney

600# weight capacity. Stainless steel construction. Features flip-down side rails. Backrest adjusts from 0° to 75°. Non-marring bumper pads. 70"l x 30"w x 39"h. Four locking casters. Ships @ 300 lbs, motor freight.



#MBG-6SS (HEAVY DUTY)
\$6850.00

Options for MBG-6SS:

IV Pole Attachment	O ₂ Tank Holder
#MBG-IV	#MBG-O2
\$330.00	\$290.00



MRI Positioning Sponges

LEG REST WITH KNEE CUTOUTS


Provides leg and knee support while maintaining hip flexion. Reduces back stress. Cleanable coating.

Additional MRI positioners and table pads can be found on page 69.

#YMCP (Large)
(patients > 6' tall)

18 x 24 x 10


\$450.00



#YMCM (Medium)
(patients 5'-6' tall)

12½ x 16¼ x 7

\$325.00



#YMCS (Small)
(patients < 5' tall)

14½ x 18 x 4

\$265.00



MRI COMFORT PILLOW

Gently elevates head during MRI & CT exams. Includes handhold straps. 16 x 25 x 8¼. Nylon covering.

#YVFB

\$250.00



ADJUSTABLE LEG POSITIONER

Specifically designed for RADIUS table tops. Removable pads allow three different knee heights. Cleanable coating.

#MCP-AP

\$570.00



MRI Warning Signs

Durable plastic silk-screened signs. Double-faced tape is on the back for mounting. Additional signs are available on pages 50 & 51.

MR	cat. #	sign #	description	size	price
	CMF-01	A	yellow, 1 symbol	8"x 10"	\$19.00
	CMM-06	B	yellow, 6 symbols	8"x 10"	\$27.00
	MBN-13E	C	restricted access (Eng)	12"x 18"	\$53.00
	MBN-14S	D	restricted access (Sp)	12"x 18"	\$53.00



MRI Earplugs

Economical protection. Latex-free. Use w/ headset or alone. Bulk packaged as 1000 pcs.



NRR: 29db

#MBE-29

\$185.00



NRR: 33db

#MBE-33

\$170.00



Metal Detector


Save patients and staff from embarrassing and hazardous situations. Features large scan surface w/ sensitivity reduction switch. Ships @ 15 lbs, UPS/FedEx.



#MBM-MD

\$325.00





MRI Transfer Boards

Allows safe and efficient transfer of patients for general imaging, MRI, and CT. Durable, semi-rigid, slick material makes positioning under the patient easy. Transfer boards are radiolucent can be used in the imaging field. Conveniently located 1½\" x 4½\" hand-holds allow optimal movement. Regular Boards are made of ¾\" thick HDPE (high density polyethylene). Regular board color is white. Transfer boards ship at OverSize rates.




REGULAR BOARDS (white)

cat. #	description	price
TBS-2372	standard	\$185.00
TBN-1872	narrow	\$165.00
TBB-3272	bariatric	\$235.00
TBH-2332	half	\$90.00

ANTI-STATIC BOARDS (black)

cat. #	description	price
TBS-2372AS	standard	\$305.00
TBN-1872AS	narrow	\$285.00
TBB-3272AS	bariatric	\$650.00
TBH-2332AS	half	\$150.00



MRI Non-Ferromagnetic Sandbags



3-lb
8\" x 8\"
#MSB-NF03
\$48.00

5-lb
5½\" x 15\"
#MSB-NF05
\$54.00

7-lb
9\" x 12\"
#MSB-NF07
\$60.00

10-lb
9\" x 14\"
#MSB-NF10
\$68.00

12½-lb
12\" x 12\"
#MSB-NF12
\$75.00

Custom upholstered sandbags with supportive handles help with patient positioning. Non-ferrous, non-magnetic, MRI Safe bags come in five weights in order to provide the ideal set-up. Order individual bags or one set that includes one of each size bag for a total of five bags.



Set of All
Five Bags
#MSB-NFST
\$295.00




Types of Positioning Sponge Foam

Techno-Aide manufactures and provides FOUR different types of sponge foam that meet a wide range of needs to ensure proper patient positioning, imaging efficiency, and accuracy. All sponges are radiolucent and fire-retardant.



Standard Foam (STD) : standard open-cell foam for general positioning. Color: light green. Shown on pages 54-65.



Firm/Bariatric Foam (FRM) : a higher density open-cell foam than standard foam, and is used to support larger/heavier patients for positioning. Color: light green. Shown on pages 60 & 64.



Closed Cell Foam (CLSD) : extra firm closed-cell foam resists fluid absorption. It can only be custom ordered and is limited to a gray rectangular block up to 4" thick. Typically used as weight-bearing.



Radiolucent Gel Pads (GEL) : radiolucent gel pads that aid in pressure distribution. Made of non-conductive cross lined polymer that will not leak, ooze, or evaporate. It is also tear and puncture resistant. Shown on pages 70-71.



MORE INFORMATION!

To compare Techno-Aide sponge imaging to other brands, watch the slide show and video from this link.

COMPARE TO CFI

Techno-Aide offers compatible options for most sponges and pads formerly offered by CFI.


Positioning Sponges are Available in these Finishes:

NON-COATED



Non-Coated Standard

Basic, open-cell positioning sponge shown above (STD). Pages 54-65.



Non-Coated Firm/Bariatric

Firm, open-cell positioning sponge shown above (FRM). Pages 60 & 64.



Non-Coated Closed-Cell

Closed-cell positioning sponge (CLSD) shown above. Custom only.

COATED



Coated Standard

Standard foam (STD) with cleanable, fluid-resistant coating. Minimal shadowing with low attenuation. Light green. Pages 54-65.



Coated Firm/Bariatric

Firm foam (FRM) with cleanable, fluid-resistant coating. Minimal shadowing with low attenuation. Light green. Pages 60 & 64.



Coated Closed-Cell

Closed-cell positioning sponge (CLSD) shown above with a clear coating for added fluid resistance and easier cleaning. Custom only.

COVERED



Vinyl Covered

Firm support foam (FRM) sewn into a covering of durable royal blue upholstery vinyl. Pages 66-68.



Heavy Duty Vinyl Covered

Firm support foam (FRM) sewn into a covering of nylon reinforced heavy vinyl. Blue or black. (Not shown in catalog.)



Healthcare Fabric Covered

Custom covering that features black 4-way stretch fabric that is anti-bacterial, anti-fungal, soft, and comfortable. Custom only.

CUSTOM Positioning Sponges

NEW CUSTOM CONFIGURED Positioning Sponges! Copy this form with a dimensioned diagram of the shape, size, and thickness of your sponge for a quick and easy quote:

Contact Name	
Address	
City/State/Zip	
Phone	
Email	

This image shows a full page of blank graph paper. The grid consists of small, uniform squares formed by thin, light gray lines. There are no margins, text, or other markings on the page.

Please show as many views as you can: front, side, top, oblique view, etc. Write length, width, height/thickness, and critical angle dimensions directly on your sketch or in the labeled spaces provided below.

Basic Shape: ☐ Rectangle ☐ Wedge ☐ Cylinder

☐ Other: _____

Foam Type: ☐ Standard Foam (STD) ☐ Closed Cell Foam (CLSD)
 ☐ Firm/Bariatric Foam (FRM)

Desired Finish: ☐ Non-Coated ☐ Coated ☐ Healthcare Fabric
☐ Vinyl Covered* ☐ Heavy Duty Vinyl Covered* (black or blue)

***Vinyl Color:**_____ (see vinyl options on p.78)

Size Factors :

- ☐ Use flared sides (“Stealth” style) in order to minimize line artifacts (understanding this will increase the overall size).
- ☐ Use exact measurements shown because line artifacts are not an issue for the intended use of this sponge.

Usage: please describe the clinical use for this sponge:

Write your dimensions directly on your sketch or label your lines and angles with these letters and write your dimensions in the middle column.

Dimension A		Inches
Dimension B		Inches
Dimension C		Inches
Dimension D		Inches
Dimension E		Inches
Dimension F		Inches
Dimension G		Inches
Dimension H		Inches

Revolutionary Design Minimizes Ghosting and Line Artifacts!



Stealth-Core™ sponges have been designed with unique angles to offer artifact-free viewing and increased stability.



Stealth-Cote™ sponges include a proprietary coating to add cleanability and durability to every sponge and includes antimicrobial protection.
(The **Stealth-Core™** and **Stealth-Cote™** sponges are Patent Pending.)

Stealth Rectangle Positioning Sponges

Overall dimensions: width X length X height.
May be larger than traditional style.



#YFBL
\$12.00



3" RECTANGLE
BASE: 10.25 X 12.5 X 3
TOP: 7.5 x 10



#YCBL
\$55.00

#YFBN
\$32.00



3" RECTANGLE
BASE: 16.5 X 20.5 X 3
TOP: 14 x 18




#YCBN
\$135.00

#YFBO
\$23.00



4" SQUARE SPONGE
BASE: 13.5 X 13.5 X 4
TOP: 10 x 10




#YCBO
\$73.00

#YFBY
\$56.00



4" SQUARE SPONGE
BASE: 21.5 X 27.5 X 4
TOP: 18 x 24



#YCBY
\$208.00

Non-Stealth Traditional Rectangle Sponges

Overall dimensions: width X length X height.

Standard Uncoated

Standard Coated

#YFGG
\$12.00



1.5" RECTANGLE
8 X 18 X 1.5




#YCGG
\$80.00

#YFEO
\$40.00



1.5" RECTANGLE
14 X 21 X 1.5




#YCEO
\$115.00

#YFEM
\$26.00



2.5" RECTANGLE
15 X 24 X 2.5




#YCEM
\$180.00

#YFBR
\$3.50



1" RECTANGLE
4 X 6 X 1



#YCBR
\$11.00

#YFBS
\$4.50



1.5" RECTANGLE
4 X 6 X 1.5



#YCBS
\$20.00

#YFBT
\$5.50



2" RECTANGLE
4 X 6 X 2



#YCBT
\$25.00

#YFBV
\$16.00




2" RECTANGLE
6 X 6 X 2




#YCBV
\$30.00

#YFEH
\$18.00



4" RECTANGLE
5.5 X 20 X 4



#YCEH
\$60.00

Stealth Wedge Positioning Sponges



Overall dimensions: width X length X height.
May be larger than traditional style.



#YFBE
\$13.00/set



34° WEDGE (SET OF 2)
5.5 X 9.75 X 3



#YCBE
\$48.00/set

#YFBH
\$18.00/set



16° WEDGE (SET OF 2)
8.5 X 11 X 3.5



#YCBH
\$86.00/set

#YFBI
\$16.00/set



15° WEDGE (SET OF 2)
12.5 X 10 X 3



#YCBI
\$108.00/set

#YFBJ
\$14.00/set



12° WEDGE (SET OF 2)
10.5 X 9.25 X 2.25



#YCBJ
\$82.00/set

Non-Stealth Traditional Wedge Sponges

Overall dimensions: width X length X height.

Standard Uncoated

Standard Coated

#YFFD
32.00



45° SPINAL WEDGE
7 X 21 X 7



#YCFD
\$108.00

#YFDA
\$174.00



45° CARDIAC WEDGE
14 X 28 X 14



#YCDA
\$260.00

#YFDB
\$390.00



30° CARDIAC WEDGE
28 X 27 X 16



#YCDB
\$458.00

#YFCL
\$78.00



21° EPIDURAL WEDGE
24 X 27 X 7



#YCCL
\$195.00

#YFEG
\$58.00



65° WEDGE
16 X 25 X 8.25



#YCEG
\$175.00

Stealth Wedge Positioning Sponges



#YFEA
\$38.00



45° SPINAL WEDGE
UNCOATED: 13.75 X 24 X 9.5
COATED: 13.75 X 33 X 9.5



#YCEA
\$123.00

#YFEB
\$44.00



45° SPINAL WEDGE
10.5 X 30 X 7.5



#YCEB
\$118.00

#YFEE
\$38.00



45° ECHOCARDIOGRAPHY WEDGE
UNCOATED: 13 X 30 X 9
COATED: 13 X 38.25 X 9



#YCEE
\$130.00

#YFBD
\$32.00



45° SPINAL WEDGE
10 X 27.75 X 7



#YCBD
\$112.00

#YFBB
23.00



45° WEDGE
10 X 20.25 X 7



#YCBB
\$80.00

#YFBC
\$18.00



45° WEDGE
10 X 13.5 X 7



#YCBC
\$50.00

#YFBU
\$29.00



30°-60°-90° MULTI-ANGLE
10 X 29.5 X 6



#YCBU
\$175.00

#YFBF
\$20.00



27° WEDGE
12.75 X 15.5 X 6



#YCBF
\$66.00

#YFBG
\$14.00



21° WEDGE
12.25 X 14.25 X 4.5



#YCBG
\$68.00

Stealth Extremity Positioning Sponges

Overall dimensions: width X length X height.
May be larger than traditional style.



#YFAA
\$12.00



ADULT OBLIQUE
FINGER BLOCK

#YCAA
\$30.00



#YFAB
\$24.00



BILATERAL ADULT
OBLIQUE FINGER BLOCK

#YCAB
\$52.00



#YFAC
\$12.00



BILATERAL PEDIATRIC
OBLIQUE FINGER BLOCK

#YCAC
\$27.00



#YFAD
\$22.00



ADULT LATERAL
FINGER BLOCK

#YCAD
\$68.00



#YFAE
\$20.00



PEDIATRIC LATERAL
FINGER BLOCK

#YCAE
\$34.00



#YFAH
\$42.00



KNEE ARTHROGRAM BLOCK
14.5 X 19 X 8

#YCAH
\$120.00



#YFAJ
\$98.00



CT KNEE WEDGE SUPPORT
22 X 20 X 7

#YCAJ
\$289.00



#YFCB
\$36.00



MEDIUM EXTREMITY BLOCK
11.5 X 19.75 X 5

#YCCB
\$110.00



#YCNZ
\$175.00



15° PODIATRY AXIAL / SESAMOID
8 X 13 X 4
CLOSED CELL FOAM with CLEAR COAT
(For viewing metatarsals, sesamoids,
and plantar condyles)



Stealth Head/Neck Positioning Sponges

Overall dimensions: width X length X height.
May be larger than traditional style.



#YFAK
\$12.00



CIRCULAR HEAD BLOCK
BASE: 9dia. X 1.5
TOP: 7dia.

#YCAK
\$28.00



#YFAL
\$26.00



HEAD & NECK SUPPORT
14.5 X 16 X 6.5

#YCAL
\$78.00



#YFAM
\$42.00



AP/PA ADULT HEAD HOLDER
14 X 15 X 6.75

#YCAM
\$130.00



#YFAN
\$25.00



AP/PA CHILD HEAD HOLDER
9.75 X 11 X 5

#YCAN
\$80.00



#YFAP
\$50.00



ADULT HEAD IMMOBILIZER
9 X 15.25 X 7

#YCAP
\$155.00



#YFAQ
\$35.00



CHILD HEAD IMMOBILIZER
6.75 X 11.5 X 5.25

#YCAQ
\$88.00



#YFAR
\$20.00



CERVICAL HEAD REST
7.5 X 16 X 4.5

#YCAR
\$74.00



Bariatric Stealth Positioning Sponges

Overall dimensions: width X length X height.

May be larger than traditional style.

HIGH-DENSITY SPONGE: These positioning blocks are specifically designed to provide extra firm support for bariatric and athletic patient positioning.



#YBAT
\$52.00



BARIATRIC TORSO BOLSTER
18.5 X 16.5 X 7



#YCAT
\$145.00

#YBBA
\$40.00



BARIATRIC 45° SPINAL WEDGE
10 X 27.5 X 7



#YCBA
\$119.00

#YBBM
\$20.00



BARIATRIC RECTANGLE
BASE: 10.5 X 14.5 X 3
TOP: 8 x 12




#YCBM
\$60.00

#YBBP
\$85.00



DECUBITUS PAD
BASE: 21.5 X 27.5 X 4
TOP: 18 x 24



#YCBP
\$194.00

#YFEU
\$40.00



BARIATRIC 30°-60°-90° WEDGE
10 X 29.5 X 6



#YCEU
\$195.00

Non-Stealth Traditional Bariatric Wedges

Overall dimensions: width X length X height.


Standard Uncoated

Standard Coated

#YFCG
\$42.00



45° BARIATRIC WEDGE
9 X 30 X 9



#YCCG
\$135.00

#YFHB
\$48.00



BARIATRIC SPINAL WEDGE
7 X 24 X 5



#YCHB
\$165.00

Specialty Stealth Positioning Sponges

Overall dimensions: width X length X height.
May be larger than traditional style.



#YBUZ
\$99.00



ENDO-ULTRASOUND WEDGE
20 X 20 X 6



#YCUZ
\$250.00

#YFEQ
\$70.00



PERINEAL WEDGE
20 X 20 X 3



#YCEQ
\$160.00

#YFAY
\$50.00



ANGULAR BOLSTER
14 X 22.5 X 7.25



#YCAY
\$154.00

Coated-Only Stealth Specialty Sponges


Overall dimensions: width X length X height.

#YCET
\$190.00



CT/CATH LAB HEADREST
12.75 X 13 X 4.5
(Neck Cradle 3"h)

CT HEADREST
12.5 X 12 X 9.75
(Neck Cradle 3.75"h)




#YCBW
\$225.00

#YCEV
\$250.00



CONVERTIBLE
CT HEADREST
14 X 13.5 X 6.75
(2-Pc Flat and
Convex Bases)

PHILIPS CT HEADREST
12 X 12 X 6.5



#YCEW
\$225.00

#YCEZ
\$89.00



PHLEBOTOMY WEDGE
5.5 X 25 X 5.5

CONTOURED TORSO PAD
18.5 X 17 X 5



#YCMG
\$195.00

General Sponge Kits

Overall dimensions: width X length X height.
May be larger than traditional style.



#YKGP
\$66.00



GENERAL KIT P

Uncoated OR Coated Kit includes:


uncoated	description	coated
YFBD	45° Spine Wedge	YCBD
YFBI	15° Wedge	YCBI
YFAK	Circular Head Block	YCAK
YFAA	Oblique Finger Block	YCAA



#YSGP
\$256.00




#YKGA
\$66.00



GENERAL KIT A

Uncoated OR Coated Kit includes:

uncoated	description	coated
YFAK	Circular Head Block	YCAK
YFBC	45° Wedge	YCBC
YFBE	(2) 34° Wedges	YCBE
YFBI	(2) 15° Wedges	YCBI
YFBL	Rectangle	YCBL



#YSGA
\$266.00

#YKGB
\$86.00



GENERAL KIT B


Uncoated OR Coated Kit includes:

uncoated	description	coated
YFAA	Adult Obl Finger Block	YCAA
YFAK	Circular Head Block	YCAK
YFBC	45° Wedge	YCBC
YFBE	(2) 34° Wedges	YCBE
YFBI	(2) 15° Wedge	YCBI
YFBJ	(2) 12° Wedges	YCBJ
YFBL	Rectangle	YCBL



#YSGB
\$366.00

#YKGD
\$126.00



GENERAL KIT D

Uncoated OR Coated Kit includes:

uncoated	description	coated
YFAL	Head & Neck Support	YCAL
YFAN	AP/PA Child Head Holder	YCAN
YFAR	Cervical Head Rest	YCAR
YFBF	27° Wedge	YCBF
YFBJ	(2) 12° Wedges	YCBJ
YFBL	Rectangle	YCBL



#YSGD
\$396.00

#YKGE
\$176.00



GENERAL KIT E


Uncoated OR Coated Kit includes:

uncoated	description	coated
YFAA	Adult Obl Finger Block	YCAA
YFAK	Circular Head Block	YCAK
YFAM	AP/PA Adult Head Holder	YCAM
YFBB	45° Wedge	YCBB
YFBG	21° Wedge	YCBG
YFBI	(2) 15° Wedges	YCBI
YFBL	Rectangle	YCBL
YFBN	Rectangle	YCBN
YFBU	Multi-Angle Wedge	YCBU



#YSGE
\$716.00


#YKGX
\$226.00



GENERAL KIT X

Uncoated OR Coated Kit includes:

uncoated	description	coated
YFAK	Circular Head Block	YCAK
YFBC	45° Wedge	YCBC
YFBE	(2) 34° Wedges	YCBE
YFBJ	(2) 12° Wedges	YCBJ
YFBL	Rectangle	YCBL
YFBN	Rectangle	YCBN
Pad-51	Standard Table Pad	Pad-51



#YSGX
\$446.00

General Sponge Kits

Overall dimensions: width X length X height.
May be larger than traditional style.



#YKWA
\$56.00




WEDGE KIT A
Uncoated OR Coated Kit includes:

uncoated	description	coated
YFBC	45° Wedge	YCBC
YFBF	27° Wedge	YCBF
YFBG	21° Wedge	YCBG
YFBI	(2) 15° Wedges	YCBI

#YSWA
\$266.00




#YKWB
\$76.00



WEDGE KIT B
Uncoated OR Coated Kit includes:

uncoated	description	coated
YFBD	45° Spine Wedge	YCBD
YFBE	(2) 34° Wedges	YCBE
YFBF	27° Wedge	YCBF
YFBG	21° Wedge	YCBG

#YSWB
\$266.00



#YKRA
\$36.00




RECTANGLE KIT
Uncoated OR Coated Kit includes:

uncoated	description	coated
YFBQ	(3) Rectangles	YCBQ
YFBR	(3) Rectangles	YCBR
YFBS	(3) Rectangles	YCBS
YFBT	(3) Rectangles	YCBT

#YSRA
\$180.00




#YKOA
\$106.00




ORTHO KIT A
Uncoated OR Coated Kit includes:

uncoated	description	coated
YFAA	Adult Obl Finger Block	YCAA
YFAM	AP/PA Adult Head Holder	YCAM
YFAR	Cervical Head Rest	YCAR
YFBD	45° Spine Wedge	YCBD
YFBL	Rectangle	YCBL

#YSOA
\$366.00




#YKOB
\$116.00



ORTHO KIT B
Uncoated OR Coated Kit includes:

uncoated	description	coated
YFAL	Head & Neck Support	YCAL
YFAN	AP/PA Child Head Holder	YCAN
YFAR	Cervical Head Rest	YCAR
YFBC	45° Wedge	YCBC
YFBG	21° Wedge	YCBG
YFBJ	(2) 12° Wedges	YCBJ
YFBL	Rectangle	YCBL

#YSOB
\$436.00




#YKOC
\$206.00



ORTHO KIT C
Uncoated OR Coated Kit includes:

uncoated	description	coated
YFAA	Adult Obl Finger Block	YCAA
YFAD	Adult Lateral Finger Block	YCAD
YFAF	Leg Immobilizer	YCAF
YFAH	Knee Arthrogram Block	YCAH
YFAK	Circular Head Block	YCAK
YFAL	Head & Neck Support	YCAL
YFAR	Cervical Head Rest	YCAR
YFBU	Multi-Angle Wedge	YCBU

#YSOC
\$706.00



Specialty Sponge Kits

Overall dimensions: width X length X height.
May be larger than traditional style.



#YKCA
\$56.00




CLINIC KIT A		
Uncoated OR Coated Kit includes:		
uncoated	description	coated
YFBB	45° Wedge	YCBB
YFAA	Adult Obl Finger Block	YCAA
YFAK	Circular Head Block	YCAK
YFBI	(2) 15° Wedges	YCBI



#YSCA
\$226.00



#YKCB
\$116.00




CLINIC KIT B		
Uncoated OR Coated Kit includes:		
uncoated	description	coated
YFAA	Adult Obl Finger Block	YCAA
YFAK	Circular Head Block	YCAK
YFBC	45° Wedge	YCBC
YFBE	(2) 34° Wedges	YCBE
YFBI	(2) 15° Wedges	YCBI
YFBJ	(2) 12° Wedges	YCBJ
YFBL	Rectangle	YCBL
YFBD	45° Wedge	YCBD



#YSCB
\$466.00

#YKPA
\$146.00



PEDIATRIC KIT		
Uncoated OR Coated Kit includes:		
uncoated	description	coated
YFAC	Bilateral Child Obl Finger Block	YCAC
YFAE	Pediatric Lateral Finger Block	YCAE
YFAK	Circular Head Block	YCAK
YFAN	AP/PA Child Head Holder	YCAN
YFAQ	Child Head Immobilizer	YCAQ
YFAX	Child Myelogram Block	YCAX
YFBT	Rectangle	N/A



#YSPA
\$346.00
(YSPA does NOT include the rectangle shown in the YKPA.)

#YKBA
\$176.00
(YKBA uses the high-density bariatric foam as standard.)




BARIATRIC KIT		
Uncoated OR Coated Kit includes:		
uncoated	description	coated
YBAT	Torso Bolster	YCAT
YBBA	Spinal Wedge	YCBA
YBBM	Rectangle	YCBM
YBBP	Decubitus Pad	YCBP




#YSBA
\$466.00

#YKTB
\$176.00



TRAUMA KIT		
Uncoated OR Coated Kit includes:		
uncoated	description	coated
YFAH	Knee Arthrogram Block	YCAH
YFAP	Adult Head Immobilizer	YCAP
YFAY	Angular Bolster	YCAY
YFBB	45° Wedge (7x14x7)	YCBB
YFBH	(2) 16° Wedge (7x8x3.5)	YCBH
YFBL	Rectangle (7.5x10x3)	YCBL



#YSTB
\$586.00

#YKHA
\$186.00



HEAD & NECK KIT		
Uncoated OR Coated Kit includes:		
uncoated	description	coated
YFAK	Circular Head Block	YCAK
YFAL	Head & Neck Support	YCAL
YFAM	AP/PA Adult Head Holder	YCAM
YFAN	AP/PA Child Head Holder	YCAN
YFAP	Adult Head Immobilizer	YCAP
YFAQ	Child Head Immobilizer	YCAQ
YFAR	Cervical Head Rest	YCAR



#YSHA
\$566.00

Coated-Only Specialty Sponge Kits

Overall dimensions: width X length X height.

#YSCT
\$656.00

CT IMAGING KIT 1
Coated Kit includes:

description	coated
45° Wedge	YCBB
45° Wedge	YCBC
(2) 34° Wedges	YCBE
27° Wedge	YCBF
21° Wedge	YCBG
(2) 15° Wedges	YCBH
(2) 12° Wedges	YCBJ
CT Headrest	YCBW



#YSCT2
\$715.00

CT IMAGING KIT 2
Coated Kit includes:

description	coated
CT Knee Wedge	YCAJ
Adult Head Immobilizer	YCAP
(2) 34° Wedges	YCBE
(2) 12° Wedges	YCBJ
CT Headrest	YCBW



#YSCT3
\$795.00

CT IMAGING KIT 3
Coated Kit includes:

description	coated
CT Knee Wedge	YCAJ
Adult Head Immobilizer	YCAP
CT/Cath Lab Headrest	YCET
Convertible CT Headrest	YCEV



#YSCT4
\$650.00

CT IMAGING KIT 4
Coated Kit includes:

description	coated
CT Knee Wedge	YCAJ
45° Wedge	YCBC
27° Wedge	YCBF
21° Wedge	YCBG
(2) 16° Wedges	YCBH
(2) 12° Wedges	YCBJ
45° Spinal Wedge	YCFD



#YSRC
\$496.00

ULTRASOUND KIT
Coated Kit includes:


description	coated
(10) Rectangles	YCBT



#YSAS
\$299.00

**KNEE ARTHROPLASTY
STRESS ASSEMBLY KIT**
Coated Kit includes:

description
45° Wedge (9x24x9)
Contoured Cube (11x9x6)



#YSRD
\$570.00

GENERAL / PERINEAL KIT
Coated Kit includes:

description	coated
Angular Bolster	YCAJ
45° Echocardiography Wedge	YCEE
Endo-Ultrasound Wedge	YCUZ



#YSU8
\$315.00

ULTRASOUND: BREAST OR PARACENTESIS
Coated Kit includes:

description	coated
21° Wedge	YCBG
30°-60°-90° Multi-Angle Wedge	YCBU
45° Spinal Wedge	YCFD



Vinyl-Covered Bolster Blocks & Kits

Reusable blocks, bolsters, and wedges are manufactured from firm, high-density urethane foam. Royal Blue outer covering is durable, heavy gauge, reinforced vinyl. Seams are sewn and tear resistant. Easy to clean Bolsters are RADIOLUCENT, however, the seams may show. Provides excellent support for a wide variety of imaging needs.

#YVF45
Cylinder
26" x 6" dia.
\$124.00



#YVFA
Cylinder
24" x 8" dia.
\$130.00



#YVFB
Semi-Round
8" x 18" x 6"
\$99.00




#YVFJ
Block
8" x 12" x 3"
\$74.00



#YVFG
Block
8" x 12" x 6"
\$85.00



#YVFH
Block
12" x 14" x 3"
\$60.00



#YVFP
Block
18" x 24" x 4"
\$150.00



#YVFK
Block
20" x 13" x 3"
\$85.00



#YVFI
Block
24" x 15" x 3"
\$95.00



#YVF42
Block
12" x 12" x 9"
\$260.00



#YVF43
Block
12" x 12" x 12"
\$305.00



#YVFQ
DEXA Block
20" x 20" x 12"
\$325.00



Bolster Kits



#YSU1
KIT Includes:
YVFB, (2) YVFG, YVFJ
\$289.00



#YSU2
KIT Includes:
(2) YVFG, YVFJ
\$200.00



#YSU3
KIT Includes:
YVFB, YVFD, YVFG
\$330.00

Vinyl-Covered Bolsters Wedges & Kits




#YVFU
30°/ 60°/ 90° Wedge
10" x 24" x 6"
\$180.00



#YVFO
45° Wedge
7" x 21" x 7"
\$170.00



#YVFD
45° Wedge
8" x 16" x 8"
\$195.00




#YVFC
20° Wedge
24" x 18" x 6"
\$170.00




#YVDB
30° CardiacWedge
28" x 27" x 16"
\$348.00



#YVDA
45° CardiacWedge
28" x 14" x 14"
\$215.00



#YVFE
20° Wedge
20" x 32" x 12"
\$239.00



#YVFF
50° Wedge
20" x 32" x 12"
\$310.00



#YVF40
26° Wedge/Knee
13" x 17" x 6"
\$265.00

Bolster Kits



#YSU4
KIT Includes:
YVFA, YVFB, YVFC, YVFD,
YVFG, YVFH, YVFJ, YVFK
\$775.00



#YSU5
KIT Includes:
YVFB, YVFC, YVFD,
YVFJ, YVFO
\$640.00



#YSU6
KIT Includes:
YVFJ, YVFK,
YVFO, YVFO
\$520.00

Vinyl-Covered Specialty Bolsters

Reusable bolsters are manufactured from firm, high-density urethane foam. Royal Blue outer covering is durable, heavy gauge, reinforced vinyl. Seams are sewn and tear resistant. Easy to clean Bolsters are **RADIOLUCENT**, however, the seams may show. Provides excellent support for a wide variety of imaging needs.



#YVFM
Knee Bolster
16" x 25" x 8.25"
\$250.00



#YVFZ
Endo Ultrasound Wedge
20" x 20" x 6"
\$410.00



#YVF41
Tapered Stretcher/Wheelchair Bolster
16"-14"w x 28"long x 4.5"-3"thick
\$495.00

#TBC-2214
Lateral Decubitus Positioner Kit
Pad: 20" x 13" x 3"
\$549.00




This KIT includes a 90° support with hook and loop stability straps and a YVFK decubitus bolster to aid in cross lateral examinations. The pad assures the patient is elevated into the imaging field.

#YVWS
Curved Wheelchair Chest Bolster
18" x 20" x 5"
\$250.00



The YVWS and the YVF41 (sold separately) are both designed to be perfect for chest x-ray positioning of wheelchair patients.



#TI-SET \$575.00
DEXA Positioning Kit
KIT includes:
(1) "A" Frame Foot Splint
(1) YVFJ Head Bolster (8" x 12" x 3")
(1) YVFQ Leg Bolster (20" x 20" x 12")




The large bolster aids in positioning the AP lumbar spine. The small head bolster and 15° "A" frame splint lock the feet into position to abduct both legs to hold full extension at the knees and hips.



MRI Positioners

Unique visco-elastic material molds to patient's body shape and conforms to table platform. Meets specific requirements for MRI compatibility and safety with excellent radiolucent qualities. Antimicrobial cover with anti-slip undersurface - easy to clean and maintain. Warranty: 3-months.

Additional MRI Positioners can be found on page 50.



Knee Wedge
#MFP-10



cat.#	w x l x d	price
MFP-01	14" x 3" x 1"	\$290.00
MFP-02	5" x 3" x 1"	\$290.00
MFP-03	10" x 10" x 3"	\$335.00
MFP-04	14" x 4" x 1.5"	\$335.00
MFP-05	3" x 4" x 2"	\$325.00
MFP-06	5.2" x 5.3" x 3"	\$350.00
MFP-07	3" x 7.3" x 2.3"	\$355.00
MFP-10	16" x 28.2" x 8"	\$685.00

Highest Quality Solution for Patient Comfort and Safety for all MRI Exams!

MRI Table Pads

Warranty: 6-months.

GE 1.5T & 3.0T			GE 0.2T, 0.35 & 0.7T			Hitachi			Philips, Siemens, Toshiba		
cat.#	w x l x d	price	cat.#	w x l x d	price	cat.#	w x l x d	price	cat.#	w x l x d	price
MFP-20	15" x 13" x 0.625"	\$325.00	MFP-30	29" x 14" x 0.625"	\$330.00	MFP-40	31" x 13.5" x 0.625"	\$325.00	MFP-50	19.5" x 14" x 0.625"	\$325.00
MFP-22	15" x 13" x 1.250"	\$325.00	MFP-32	29" x 14" x 1.250"	\$420.00	MFP-42	31" x 27.0" x 0.625"	\$700.00	MFP-52	19.5" x 14" x 1.250"	\$325.00
MFP-24	15" x 27" x 1.250"	\$790.00	MFP-34	29" x 28" x 1.250"	\$1125.00	MFP-44	27" x 13.5" x 0.625"	\$325.00	MFP-54	19.5" x 28" x 1.250"	\$795.00
(MFP-24 is two pcs)			(MFP-34 is two pcs)			MFP-46	31" x 13.5" x 1.250"	\$355.00	(MFP-54 is two pcs)		
						MFP-48	31" x 27.0" x 1.250"	\$955.00			
						(MFP-42 & MFP-48 are two pcs)					

Other Positioners & Pads are available. Call for quote.

MRI Positioner Pad Kits & Burn Pads



Burn Pads (4 kits available)



#MFP-BBK



#MFP-GE2



#MFP-GNK

(Includes: MFP-01, MFP-02, MFP-03, MFP-04, MFP-05, MFP-06)

cat.#	MRI Kit Description (# pieces)	price
MFP-GNK	General MRI Kit (12)	\$1935.00
MFP-GE2	CTL Coil Pad Kit: GE 1.5T (2)	\$810.00
MFP-BBK	Breast Biopsy Positioning Kit (10)	\$1245.00
Burn Pad Kits		
MFP-20P	MFP-20 Table Pads (2)	\$650.00
MFP-30P	MFP-30 Table Pads (2)	\$660.00
MFP-40P	MFP-40 Table Pads (2)	\$650.00
MFP-50P	MFP-50 Table Pads (2)	\$650.00

GEL Patient Positioners

Radiolucent Gel Pads aid in pressure distribution without bottoming out. These pads are made of a solid, non-conductive, cross lined polymer that will not leak, ooze, or evaporate and is resistant to punctures and tears. Each antimicrobial, anti-bacterial pad is hypoallergenic, latex-free, silicone-free, and plasticizer-free, and will not harden over time. Following the instructions, the pads can be heated and cooled. Note: 0.25" thick gel is equivalent to .051mm Aluminum.

FLAT PAD POSITIONERS

Protects skin from pressure sores and ulceration. Great for short term and long term situations.



cat. #	dimensions	wt	price	cat. #	dimensions	wt	price
YG-4008	12" x 4" x 1/4"	5	\$50.00	YG-7005	27" x 11" x 1/2"	8	\$205.00
YG-8001*	17" x 11" x 1/4"	5	\$80.00	YG-1009*	20" x 20" x 1/2"	10	\$190.00
YG-7003	11" x 20" x 1/4"	5	\$110.00	YG-1002*	39 1/2" x 20" x 1/2"	18	\$325.00
YG-8002*	11" x 9" x 3/8"	3	\$65.00	YG-1008*	46" x 20" x 1/2"	20	\$355.00
YG-7002	17" x 11" x 3/8"	4	\$80.00	YG-1005*	72" x 20" x 1/2"	31	\$580.00
YG-1001	9" x 20" x 1/2"	5	\$120.00	YG-8010*	21" x 18 1/2" x 5/8"	11	\$210.00
YG-8003*	26" x 12" x 1/2"	8	\$160.00	YG-8014	29" x 10" x 1"	12	\$115.00
YG-8005	17" x 11" x 1"	8	\$200.00	YG-8006	20" x 20" x 1"	18	\$160.00

(*Excellent as Body & Support Surface Pads)

FLAT PAD w/ CUTOUT

cat. #	dimensions	wt	price
YG-1003	20" x 20" x 1/2"	8	\$170.00
YG-1007	46" x 20" x 1/2"	20	\$340.00

Same as above with cutout for greater access to patient and table/bed.



LATERAL POSITIONER

29" x 19 1/2" x 6" x 1"
#YG-3026
\$275.00

Stable and safe cradle for arm and upper body when positioned laterally.



FLAT BOTTOM CHEST ROLLS

cat. #	dimensions	wt	price
YG-3004	12" x 4" x 4"	8	\$120.00
YG-3006	12" x 6" x 6"	15	\$275.00
YG-3008	16" x 4" x 4"	10	\$165.00
YG-3009	16" x 6" x 6"	20	\$370.00
YG-3010	20" x 4" x 4"	13	\$200.00
YG-3011	20" x 6" x 6"	25	\$490.00



Used for prone, lateral, and supine positioning. Protects chest and prevents rolling.

HEEL-ELBOW PAD

3 1/4" x 1 1/2" dia.
#YG-5025
\$50.00

Protects skin from pressure sores and ulceration in extended cases.



TOBOGGAN ARM & LEG POSITIONERS

cat. #	description	dimensions	wt	price
YG-5030	arm pad only	7 3/4" x 6 1/4" x 3/8"	3	\$65.00
YG-5020	arm pad only	17" x 7 1/2" x 3/8"	4	\$80.00
YG-7020	leg pad only	17" x 7 1/2" x 3/8"	3	\$80.00
YG-5022	shield only	16 1/2" x 7 1/2"	5	\$170.00
YG-5021	arm pad & shield	17" x 7 1/2" x 3/8"	6	\$245.00
YG-7021	leg pad & shield	17" x 7 1/2" x 3/8"	5	\$245.00

Used with toboggan positioner to protect patients' arms and legs during OR procedures.



**All Items at
New Lower
Prices!**

TRAPEZOID POSITIONER

7" x 10" x 3 1/2"
#YG-3035
\$145.00

Used for leg separation. Can be used as a knee stabilizer or under the feet to prevent foot drop.



HEEL CUP-PROTECTOR

6" x 3 1/4" x 1 3/4"
#YG-4004
\$55.00

Helps relieve heel pressure in supine position. Provides for even distribution of weight on the foot and ankle.



POSITIONING

GEL Patient Positioners

**All Items at
New Lower
Prices!**

gel pads & positioners



Cradles and supports head and allows clear air passage.

HORSESHOE OPEN HEAD RING

cat. #	description	dimensions	wt	price
YG-2007	adult	8" x 3" x 2"	4	\$105.00
YG-2009	pediatric	5½" x 2¼" x 1¼"	3	\$85.00
YG-2008	neonatal	3½" x 1½" x ¾"	2	\$50.00
YG-2012	bariatric	8" x 3" x 3"	5	\$115.00



Cradles head and protects the ear. Reduces pressure on occipital while supine.

DONUT CLOSED HEAD RING

cat. #	description	dimensions	wt	price
YG-2000	adult	8" x 3" x 2"	4	\$110.00
YG-2004	pediatric	5½" x 2¼" x 1¼"	3	\$80.00
YG-2005	neonatal	3½" x 1½" x ¾"	2	\$50.00
YG-2015	bariatric	8" x 3" x 3"	5	\$115.00



PRONE-SUPINE HEADREST

11" x 9" x 6"
#YG-2010
\$265.00

Patient can rest comfortably even with respiratory tubes in place.



Unique center dish enhances patient comfort and protection.

SUPINE HEAD PAD w/ CENTER DISH

cat. #	description	dimensions	wt	price
YG-2022	adult	10" x 9" x 2"	8	\$130.00
YG-2024	pediatric	7" x 6" x 2"	3	\$85.00
YG-2026	neonatal	5" x 4" x 2"	3	\$65.00



Conforming pad helps reduce pressure to side of patients' face and ear.

LATERAL HEAD PAD w/ CENTER DISH

cat. #	description	dimensions	wt	price
YG-2023	adult	10" x 9" x 1"	6	\$110.00
YG-2025	pediatric	7" x 6" x 1"	3	\$65.00
YG-2027	neonatal	5" x 4" x 1"	3	\$55.00

SIDE-SUPINE HEADREST

9" x 8" x 3"
#YG-2021
\$235.00



Comfortable ear cutout for lateral position.



GEL ARMBOARDS relieve pressure on the hand, forearm, elbow and bicep. Ensures proper positioning and stability.

FLAT ARMBOARD PADS

cat. #	size	dimensions	wt	price
YG-5000	small	15" x 4" x ¼"	3	\$45.00
YG-5005	large	25" x 5" x ¼"	3	\$100.00



CONTOURED ARMBOARD PADS

cat. #	size	dimensions	wt	price
YG-5026	medium	16" x 6" x 2"	7	\$160.00
YG-5027	large	21" x 6" x 2"	9	\$185.00

ULNAR & BRACHIAL NERVE PROTECTOR

15¾" x 6" x ¾"
#YG-5012
\$85.00



ULNAR & BRACHIAL NERVE PROTECTOR

17" x 11" x ¾"
#YG-5011
\$98.00



AXILLARY ROLL COVER

17" x 12" x ½"
#YG-3032
\$155.00



GEL Sandbags

**All Items at
New Lower
Prices!**

cat. #	dimensions	wt	price
YG-6000	6" x 4½" x 1½"	2	\$55.00
YG-6001	6" x 4½" x 2"	2	\$55.00
YG-6002	8" x 4½" x 1½"	2	\$60.00
YG-6003	10" x 4½" x 1½"	3	\$60.00
YG-6004	12" x 4½" x 1½"	5	\$60.00
YG-6008	10" x 4½" x 2"	5	\$95.00
YG-6010	12" x 4½" x 2"	4	\$95.00
YG-6011	14" x 5" x 2½"	5	\$95.00
YG-6012	16" x 12" x 1½"	12	\$205.00
YG-3015	10" x 4½" x 2"	3	\$80.00
YG-3017	14" x 5" x 2"	5	\$95.00
YG-3025	14" x 9" x 2½"	10	\$260.00

Extremely versatile positioners are soft on the surface, but firm inside where it counts. Each sandbag provides exceptional patient support as well as increased patient stability.



Sandbags

All sandbags are filled with sand and are ready to use. Heat sealed inner sleeve prevents the leakage of sand. All seams are sewn and tear resistant. Deep green outer cover is heavy gauge, reinforced vinyl. Every sandbag includes a reinforced webbing handle.



Custom Sizes Available!

Includes weight ID tag



						
2-lb 4" x 9" #SBP-02 \$23.00	3-lb 5" x 9" #SBP-03 \$23.00	5-lb 7" x 9" #SBP-05 \$30.00	5-lb 6" x 12" #SBF-05 \$30.00	5-lb 9" x 9" #SBC-05 \$30.00	7-lb 10" x 10" #SBG-07 \$35.00	7-lb 6" x 14" #SBF-07 \$35.00
						
10-lb 7" x 16" #SBF-10 \$36.00	10-lb 11" x 11" #SBC-10 \$36.00	10-lb 11" x 11" #SBC-HL \$48.00	10-lb 11" x 11" #SBC-AC \$45.00	15-lb 12" x 12" #SBG-15 \$41.00	20-lb 12" x 14" #SBC-20 \$41.00	

Sandsocks

1-lb 2" x 8"
#SBS-01
\$30.00 ea.



3-lb 4" x 12"
#SBS-03
\$31.00 ea.



5-lb 5" x 15"
#SBS-05
\$35.00 ea.



10-lb 6" x 20"
#SBS-10
\$38.00 ea.



Includes weight ID tag



Sandsocks are loosely filled sandbags that contour to a broad range of shapes and forms. Stabilizes hard-to-hold and irregular anatomy. No need for additional person to hold the patient. Double seams for strength and longevity. Flexible, heavy gauge, reinforced vinyl. Only available in "honey" color.



12-lb 6" x 25"
#SBS-12
\$42.00 ea.



10-lb 6" x 28"
#SBS-28
\$42.00 ea.



15-lb 6" x 40"
#SBS-40
\$52.00 ea.

Sandbag Kits (Save over ordering individually!)



5-lb cervical
#SBK-CA \$55.00 set
Set: (2) 5-lb / 9" x 9"



10-lb cervical
#SBK-CB \$70.00 set
Set: (2) 10-lb / 11" x 11"



20-lb cervical
#SBK-CC \$80.00 set
Set: (2) 20-lb / 12" x 14"



10-lb AC Joint w/wrist grips
#SBK-AC \$85.00 set
Set: (2) 10-lb / 11" x 11"



10-lb w/handles
#SBK-HL \$90.00 set
Set: (2) 10-lb / 11" x 11"



femoral/angio
#SBK-FB \$85.00 set
Set: (1) 5-lb / 5" x 12" (1) 7-lb / 6" x 14"
(1) 10-lb / 7" x 16"



general set (a)
#SBK-GA \$230.00 set
Set: (2) 3-lb / 5" x 9" (2) 10-lb / 11" x 11"
(2) 7-lb / 6" x 14" (1) 5-lb / 7" x 9"
(1) 5-lb / 5" x 12"



general set (b)
#SBK-GB \$180.00 set
Set: (2) 3-lb / 5" x 9"
(2) 7-lb / 6" x 14"
(2) 10-lb / 11" x 11"



pediatric
#SBK-PA \$125.00 set
Set: (1) 1-lb / 4" x 5" (1) 5-lb / 7" x 9"
(1) 2-lb / 4" x 9" (1) 7-lb / 6" x 14"
(1) 3-lb / 5" x 9"



sandsock kit
#SBK-SA \$165.00 set
Set: (1) 1-lb / 2" x 8" (1) 12-lb / 6" x 25"
(1) 10-lb / 6" x 20" (1) 5-lb / 5" x 15"
(1) 3-lb / 4" x 12"



athletic/bariatric
#SBK-AA \$210.00 set
Set: (2) 7-lb / 10" x 10"
(2) 10-lb / 7" x 16"
(2) 20-lb / 12" x 14"

Mobile Sandbag Rack

#SBM-100
\$320.00



Sandbags
not included




Organize & store sandbags for easy retrieval. 2-sided panel with 20 hanging hooks. Base shelf on both sides. 3" swivel casters. 53"h x 25"w x 19"d.


(A) Table Pad Foam Types:



Comfort Foam: standard table pad foam. Great for basic table pads and is available in a variety of sizes and thicknesses.




Firm Foam: higher density foam pad than standard foam in order to provide more patient support.



Memory Foam: a bottom layer of 1" comfort foam is bonded to a top layer of 1" 5.5-lb memory foam for superior comfort. Great for long procedures and sensitive patients.

(C) Table Pad End Finishes:




Tapered End with Grommet : standard closure; side seams may be sewn or fused (depending on material) with a tapered end closure for brass grommets.



Premium Box Corner : completely sewn covering that features squared finished corners.

(B) Table Pad Cover Materials:



Standard Vinyl (light blue only): a durable pad covering that is economical and easy to clean.



Nylon Fabric (see colors on p.25): nylon coverings offer the greatest variety of solid and printed fabrics. Durable and liquid repellent.



Upholstery Vinyl (see colors on p.78): extremely durable and very easy to clean. Choose from six colors.



Heavy Duty Reinforced Mesh Vinyl (blue or black only): nylon filament mesh provides a rip-stop covering for high use pads.



Healthcare 4-Way Stretch Fabric (black only): pliable and resilient, antibacterial, antifungal covering.

Putting the Options Together:

pads	page no.	(A) foam type	(B) cover mat'l	(C) end finish
Basic	76	comfort, firm	standard vinyl, nylon, upholstery vinyl	tapered
Techno-Tuff	77	comfort, firm	reinforced mesh vinyl	tapered
Premium	78	comfort, firm, memory	upholstery vinyl, 4-way stretch fabric	boxed
Stretcher	78	comfort	4-way stretch fabric	boxed
Standard Long & Narrow	79	comfort, firm	upholstery vinyl	tapered
Premium Long & Narrow	79	comfort, firm, memory	4-way stretch fabric	boxed
Specialty Pads	80-83	comfort, memory	4-way stretch fabric	boxed
Custom	75, 77, 78	three type choices	five style choices	tapered or boxed

Techno-Aide offers many table pads previously offered by CFI. Find what you are looking for in these pages or call us for a cross-reference.

CUSTOM CONFIGURED Table Pads! Copy this form with a dimensioned diagram of the shape, size, and thickness of your table pad for a quick and easy quote:



Contact Name	
Address	
City/State/Zip	
Phone	
Email	

[illegible]

Foam Type: ☐ Comfort Foam **Desired Finish:** ☐ Tapered End
☐ Firm/Bariatric Foam ☐ Tapered w/ Grommets
☐ Memory Foam ☐ Boxed End

Cover Mat'l : ☐ Standard Vinyl ☐ Upholstery Vinyl ☐ Nylon Fabric
☐ HD Reinforced Vinyl ☐ Healthcare 4-Way Stretch Fabric

Color:_____ (see options for specific materials on p.74)

Anti-Slip Backing: ☐ Yes ☐ No **Hanging Straps:** ☐ Yes ☐ No

Add Hook & Loop Strips to the Underside of the Pad: ☐ Yes ☐ No


- If "YES," write down whether you want the Hook or the Loop sewn to the pad (This should be the opposite of what is on your table).
 - Please note on the diagram above the exact placement of the Hook or Loop and the distances from the sides and ends of the pad, and the length of the strips.
- [For reference: the HOOK is the soft & fuzzy part, the LOOP is rough & prickly.]

Write your dimensions directly on your sketch or label your lines with these letters and write your dimensions in the middle column.

Dimension A		Inches
Dimension B		Inches
Dimension C		Inches
Dimension D		Inches
Dimension E		Inches
Dimension F		Inches
Dimension G		Inches
Dimension H		Inches

Basic Table Pads

Pads can be used for imaging tables or stretchers. Fabricated for durability and longevity. Helps maintain complete control of the patient. Affordable for each room and table. Three types of covering, depending on usage. Two brass eyelets in each end for hanging. Soft and comfortable. Cleans easy with mild cleaner or germicide.



TWO FOAM CHOICES:

- 1" comfort foam
- 2" firm density foam

THREE COVER MATERIALS:

- Standard light blue vinyl (1 choice)
- Polyester-backed upholstery vinyl (6 choices)
- Nylon fabric (29 choices on page 25)



Light Blue (econ vinyl)



Nylon (solid or print)



Green Upholstery



Red Upholstery



Black Upholstery



Blue Upholstery



Yellow Upholstery



Mint Upholstery

Comfort Foam Pads

72" x 23 1/4" x 1"

cat.#	color	covering	price
PAD-51	Light Blue	Standard Vinyl	\$109.00
PAD-53	Green	Upholstery Vinyl	\$189.00
PAD-54	Red	Upholstery Vinyl	\$189.00
PAD-55	Black	Upholstery Vinyl	\$189.00
PAD-50	Blue	Upholstery Vinyl	\$189.00
PAD-52	Yellow	Upholstery Vinyl	\$189.00
PAD-56	Mint	Upholstery Vinyl	\$189.00
PAD-80	See Page 25	Nylon Fabric	\$159.00

(Edge of pads may show up on some images.)

Firm Density Pads

72" x 23 1/4" x 2"

cat.#	color	covering	price
PAH-51	Light Blue	Standard Vinyl	\$159.00
PAH-53	Green	Upholstery Vinyl	\$239.00
PAH-54	Red	Upholstery Vinyl	\$239.00
PAH-55	Black	Upholstery Vinyl	\$239.00
PAH-50	Blue	Upholstery Vinyl	\$239.00
PAH-52	Yellow	Upholstery Vinyl	\$239.00
PAH-56	Mint	Upholstery Vinyl	\$239.00
PAH-80	See Page 25	Nylon Fabric	\$209.00

(Edge of pads may show up on some images.)

Pad Storage Rack & Set

Two hooks on 12" centers.
Mounts on 16" centers. Steel
with white finish. Ships @ 4 lbs.

#PAD-R2 \$ 53.00 RACK ONLY

#PAD-S1 \$145.00 COMBO SET

(Combo Set includes: Standard Pad & Rack)



Techno-Tuff™ Table Pads

Constructed of durable mesh laminate reinforced vinyl. 72" x 23" x 1" pad is radiolucent. Handles allow patients to be easily moved from stretcher to table. Low-friction, water proof, stain proof, and antimicrobial.

TECHNO-TUFF™ PADS FEATURE:

TWO FOAM CHOICES:

1" comfort foam

2" firm density foam

TECHNO-TUFF™ COVER MATERIAL:

Mesh laminate reinforced vinyl

Colors: black or royal blue. With or without handles.



Techno-Tuff™ Comfort Foam

72" x 23¼" x 1"

cat.#	color	description	price
PAD-H8R	Royal Blue	With Handles	\$275.00
PAD-H8K	Black	With Handles	\$275.00
PAD-X8R	Royal Blue	No Handles	\$245.00
PAD-X8K	Black	No Handles	\$245.00

Techno-Tuff™ Firm Density

72" x 23¼" x 2"


cat.#	color	description	price
PAH-H8R	Royal Blue	With Handles	\$325.00
PAH-H8K	Black	With Handles	\$325.00
PAH-X8R	Royal Blue	No Handles	\$295.00
PAH-X8K	Black	No Handles	\$295.00

Custom Pads

Rectangular pads made to your specs. Submit length" x width" x thickness".

Pads larger than 80"long x 30"wide x 2"thick may require further quote.

Custom Pads are NON-RETURNABLE.



Custom Comfort Foam

Up to: 80" x 30" x 2"

cat.#	color	covering	price
PAD-SP51	Light Blue	Standard Vinyl	\$209.00
PAD-SP53	Green	Upholstery Vinyl	\$289.00
PAD-SP54	Red	Upholstery Vinyl	\$289.00
PAD-SP55	Black	Upholstery Vinyl	\$289.00
PAD-SP50	Blue	Upholstery Vinyl	\$289.00
PAD-SP52	Yellow	Upholstery Vinyl	\$289.00
PAD-SP56	Mint	Upholstery Vinyl	\$289.00
PAD-SP80	See Page 25	Nylon Fabric	\$259.00

*Sizes exceeding 80"x30"x2" in any direction call for quote.

Custom Firm Density

Up to: 80" x 30" x 2"

cat.#	color	covering	price
PAH-SP51	Light Blue	Standard Vinyl	\$259.00
PAH-SP53	Green	Upholstery Vinyl	\$339.00
PAH-SP54	Red	Upholstery Vinyl	\$339.00
PAH-SP55	Black	Upholstery Vinyl	\$339.00
PAH-SP50	Blue	Upholstery Vinyl	\$339.00
PAH-SP52	Yellow	Upholstery Vinyl	\$339.00
PAH-SP56	Mint	Upholstery Vinyl	\$339.00
PAH-SP80	See Page 25	Nylon Fabric	\$309.00

*Sizes exceeding 80"x30"x2" in any direction call for quote.

Custom Techno-Tuff™ Comfort

Up to: 80" x 30" x 2"

cat.#	color	description	price
PAD-SPH8R	Royal Blue	With Handles	\$375.00
PAD-SPH8K	Black	With Handles	\$375.00
PAD-SPX8R	Royal Blue	No Handles	\$345.00
PAD-SPX8K	Black	No Handles	\$345.00

*Sizes exceeding 80"x30"x2" in any direction call for quote.

Custom Techno-Tuff™ Firm Density

Up to: 80" x 30" x 2"

cat.#	color	description	price
PAH-SPH8R	Royal Blue	With Handles	\$425.00
PAH-SPH8K	Black	With Handles	\$425.00
PAH-SPX8R	Royal Blue	No Handles	\$395.00
PAH-SPX8K	Black	No Handles	\$395.00

*Sizes exceeding 80"x30"x2" in any direction call for quote.



Standard Long & Narrow Table Pads

NEW Long & Narrow Table Pads can be used for both long or narrow imaging tables. They are covered in very durable upholstery vinyl that is extremely easy to clean with a mild cleaner or germicide. Six (6) color available.

NEW!

cat.#	length" x width" x thickness"	foam type	price
Long Pad (standard width)			
PAD-LS77	77" x 23.25" x 1"	comfort	\$289.00
PAH-LS77	77" x 23.25" x 2"	firm	\$339.00
Narrow Pad (standard length)			
PAD-SN19	72" x 19" x 1"	comfort	\$289.00
PAH-SN19	72" x 19" x 2"	firm	\$339.00
Both Long & Narrow			
PAD-LN79	77" x 19" x 1"	comfort	\$289.00
PAH-LN79	77" x 19" x 2"	firm	\$339.00

COLORS:

Comfort & Firm Foam:

Standard Long & Narrow Table Pads are available in the six (6) upholstery vinyl colors shown below.



Green Upholstery



Red Upholstery



Black Upholstery



Blue Upholstery



Yellow Upholstery



Mint Upholstery

Premium Long & Narrow Table Pads

NEW Premium Long & Narrow Table Pads offer excellent support and durability. Squared off corners are tailored to provide a nice, finished look that is perfect for any imaging table. Pads clean easily with mild cleaner or germicide.

NEW!

cat.#	length" x width" x thickness"	foam type	price
Long Pad (standard width)			
PPAD-LS77	77" x 23.25" x 1"	comfort	\$440.00
PPAH-LS77	77" x 23.25" x 2"	firm	\$500.00
PPAM-LS77	77" x 23.25" x 2"	memory	\$735.00
Narrow Pad (standard length)			
PPAD-SN19	72" x 19" x 1"	comfort	\$325.00
PPAH-SN19	72" x 19" x 2"	firm	\$400.00
PPAM-SN19	72" x 19" x 2"	memory	\$540.00
Both Long & Narrow			
PPAD-LN79	77" x 19" x 1"	comfort	\$405.00
PPAH-LN79	77" x 19" x 2"	firm	\$455.00
PPAM-LN79	77" x 19" x 2"	memory	\$680.00

COLORS:

Comfort & Firm Foam:

available in the Blue or Black upholstery vinyls shown above.

Memory Foam:

available in the Black 4-way stretch fabric shown below.



Premium Box Corner : completely sewn covering that features squared finished corners.



Healthcare 4-Way Stretch Fabric (black only): pliable and resilient, antibacterial, antifungal covering.

Techno-Aide offers many table pads previously offered by CFI. Find what you are looking for in these pages or call us for a cross-reference.

Premium Table Pads

NEW Premium Table Pads offer excellent support and durability. Squared off corners are tailored to provide a nice, finished look that is perfect for any imaging table. Pads clean easily with mild cleaner or germicide.



Standard or
Firm Density
Foam Core

2" Premium Comfort Pad

72" x 23¼" x 2"

cat.#	color	covering	price
PPAD-55	Black	Upholstery Vinyl	\$395.00
PPAD-50	Blue	Upholstery Vinyl	\$395.00
PPAD-75	Black	4-Way Stretch	\$455.00

Custom Premium Comfort Pad

Up to: 80" x 30" x 2"

cat.#	color	covering	price
PPAD-SP55	Black	Upholstery Vinyl	\$535.00
PPAD-SP50	Blue	Upholstery Vinyl	\$535.00
PPAD-SP75	Black	4-Way Stretch	\$615.00



Squared
Box Corner

2" Premium FIRM Density Pad

72" x 23¼" x 2"


cat.#	color	covering	price
PPAH-55	Black	Upholstery Vinyl	\$420.00
PPAH-50	Blue	Upholstery Vinyl	\$420.00
PPAH-75	Black	4-Way Stretch	\$475.00

Custom Premium FIRM Density Pad

Up to: 80" x 30" x 2"

cat.#	color	covering	price
PPAH-SP55	Black	Upholstery Vinyl	\$560.00
PPAH-SP50	Blue	Upholstery Vinyl	\$560.00
PPAH-SP75	Black	4-Way Stretch	\$635.00

Memory Foam Pads are constructed of a bottom layer of 1" standard foam and a top layer of 1" 5.5-lb Memory Foam to provide superior comfort for long procedures, recovery periods, or sensitive patients. Contouring support foam is radiolucent.



2-Ply
Memory
Foam

2" Premium MEMORY Pad

72" x 23¼" x 2"

cat.#	color	covering	price
PPAM-80	Black	4-Way Stretch	\$595.00

Custom Premium MEMORY Pad

Up to: 80" x 30" x 2"

cat.#	color	covering	price
PPAM-SP80	Black	4-Way Stretch	call for price



All Premium Table Pads are available with an optional **Anti-Slip Bottom Surface** that must be added to the same order as any of the above pads as: Cat. #PAD-ASB for an additional \$75.00.

Stretcher Replacement Pads

Economically priced replacement stretcher pads offer excellent support and durability. The black cover of 4-way-stretch antimicrobial fabric is tailored fitted over our firm 3" polyfoam. Squared off corners provide a nice, finished look with a resealable hook and loop end panel. Pads are CAL 117 compliant and clean easily with mild cleaner or germicide.



Ask us about Custom Pads with the standard four squared corners.

3" Stretcher Pad

cat. #	length x width	price
PMR-7224	72" x 24"	\$575.00
PMR-7226	72" x 26"	\$590.00
PMR-7228	72" x 28"	\$635.00
PMR-7230	72" x 30"	\$685.00
PMR-7424	74" x 24"	\$585.00
PMR-7426	74" x 26"	\$600.00
PMR-7428	74" x 28"	\$655.00
PMR-7430	74" x 30"	\$700.00
PMR-7624	76" x 24"	\$585.00
PMR-7626	76" x 26"	\$605.00
PMR-7628	76" x 28"	\$655.00
PMR-7630	76" x 30"	\$705.00

3" Pad w/ Anti-Slip Bottom

cat. #	length x width	price
PMR-7224ASB	72" x 24"	\$700.00
PMR-7226ASB	72" x 26"	\$715.00
PMR-7228ASB	72" x 28"	\$760.00
PMR-7230ASB	72" x 30"	\$810.00
PMR-7424ASB	74" x 24"	\$710.00
PMR-7426ASB	74" x 26"	\$725.00
PMR-7428ASB	74" x 28"	\$780.00
PMR-7430ASB	74" x 30"	\$825.00
PMR-7624ASB	76" x 24"	\$710.00
PMR-7626ASB	76" x 26"	\$730.00
PMR-7628ASB	76" x 28"	\$780.00
PMR-7630ASB	76" x 30"	\$830.00



Also available: Custom Slant Corner Pad "A"
For stretchers with hand holes at only one end.



Also available: Custom Slant Corner Pad "B"
For stretchers with hand holes at both ends.

Specialty Premium Table Pads

Our collection of Specialty Premium Table Pads are designed to fit on specific manufacturers' imaging tables. Each pad offers excellent support and durability. The black cover of 4-way-stretch antimicrobial fabric is tailor-fitted over your choice of either 2" premium comfort foam or our 2" two-ply foam with memory. Squared off corners and ends provide a nice, finished look that is perfect for any imaging table. Pads clean easily with mild cleaner or germicide. Custom size pads are also available.


**Now Available in
both Comfort and
Memory Foams!**

All pads are shown in their proper proportions and to scale with each other.




COMFORT
#PPMS-ANS
\$725.00

MEMORY
#PMMS-ANS
\$1115.00




COMFORT
#PPMS-CYS
\$350.00

MEMORY
#PMMS-CYS
\$503.00



COMFORT
#PPMS-MEM
\$525.00

MEMORY
#PMMS-MEM
\$770.00



COMFORT
#PPMS-S50
\$996.00

MEMORY
#PMMS-S50
\$665.00



COMFORT
#PPMS-XE8
\$625.00

MEMORY
#PMMS-XE8
\$860.00



COMFORT
#PPMS-XEA
\$835.00

MEMORY
#PMMS-XEA
\$1335.00






TABLE PADS

specialty table pads for GE


COMFORT	MEMORY
#PPGE-98R	#PMGE-98R
\$645.00	\$955.00




COMFORT	MEMORY
#PPGE-AFM	#PMGE-AFM
\$960.00	\$1629.00



COMFORT	MEMORY
#PPGE-98N	#PMGE-98N
\$525.00	\$780.00




COMFORT	MEMORY
#PPGE-ANG	#PMGE-ANG
\$825.00	\$1326.00




COMFORT	MEMORY
#PPGE-NRW	#PMGE-NRW
\$630.00	\$842.00




COMFORT	MEMORY
#PPGE-LUP	#PMGE-LUP
\$950.00	\$1517.00



COMFORT	MEMORY
#PPGE-OM3	#PMGE-OM3
\$815.00	\$1317.00



COMFORT	MEMORY
#PPGE-OM4	#PMGE-OM4
\$855.00	\$1629.00



COMFORT	MEMORY
#PPGE-OM5	#PMGE-OM5
\$1040.00	\$1746.00



**Now Available in
both Comfort and
Memory Foams!**

All pads are shown in their
proper proportions and to
scale with each other.

TABLE PADS

specialty table pads for Siemens

COMFORT	MEMORY
#PPSI-ANG	#PMSI-ANG
\$630.00	\$862.00

Siemens Angio Pad

Overall Dimensions
66" x 21"

COMFORT	MEMORY
#PPSI-ANK	#PMSI-ANK
\$690.00	\$1056.00

Siemens Koordinat Angio Table Pad

Overall Dimensions
79" x 21"

COMFORT	MEMORY
#PPSI-CCL	#PMSI-CCL
\$720.00	\$1098.00

Siemens CCL Cardiac Cath

Overall Dimensions
79" x 23.5"

COMFORT	MEMORY
#PPSI-ECT	#PMSI-ECT
\$660.00	\$976.00

Siemens Orbiter ECT Special Procedure

Overall Dimensions
78" x 18.5"

COMFORT	MEMORY
#PPSI-HIQ	#PMSI-HIQ
\$680.00	\$1023.00

Siemens HiQ

Overall Dimensions
96" x 21"

COMFORT	MEMORY
#PPSI-SDR	#PMSI-SDR
\$520.00	\$746.00

Siemens Somatom DR

Overall Dimensions
82" x 16"

COMFORT	MEMORY
#PPSI-HIS	#PMSI-HIS
\$480.00	\$700.00

Siemens HiQ Short

Overall Dimensions
76" x 15"

COMFORT	MEMORY
#PPSI-SPL	#PMSI-SPL
\$480.00	\$700.00

Siemens Somatom Plus

Overall Dimensions
75" x 16"

COMFORT	MEMORY
#PPSI-HIX	#PMSI-HIX
\$640.00	\$888.00

Siemens HiQ XL

Overall Dimensions
114" x 15"

**Now Available in
both Comfort and
Memory Foams!**

All pads are shown in their
proper proportions and to
scale with each other.

TABLE PADS

specialty table pads for Toshiba / Picker / Nuclear X-Ray

COMFORT
#PPTO-150
\$730.00

MEMORY
#PMTO-150
\$1120.00

Toshiba Cat 150

Overall Dimensions
91" x 24"

COMFORT
#PPTO-ANG
\$790.00

MEMORY
#PMTO-ANG
\$1116.00

Toshiba Angio Table Pad

Overall Dimensions
104" x 21"

COMFORT
#PPTO-250
\$985.00

MEMORY
#PMTO-250
\$1629.00

Toshiba Cat 250

Overall Dimensions
104" x 27.5"

COMFORT
#PPTO-ANC
\$730.00

MEMORY
#PMTO-ANC
\$1121.00

Toshiba Angio Contoured

Overall Dimensions
96" x 24"

COMFORT
#PPTO-FX2
\$960.00

MEMORY
#PMTO-FX2
\$1638.00

Toshiba Cat FX-2

Overall Dimensions
112" x 26"

COMFORT
#PPTO-SCP
\$710.00

MEMORY
#PMTO-SCP
\$1079.00

Toshiba Special Procedure

Overall Dimensions
94" x 21"

COMFORT
#PPPI-ANG
\$730.00

MEMORY
#PMPI-ANG
\$1116.00

Picker Angio

Overall Dimensions
99" x 22"

COMFORT
#PPPI-120
\$520.00

MEMORY
#PMPI-120
\$746.00

Picker 1200

Overall Dimensions
82" x 16"

COMFORT
#PPPI-PCT
\$530.00

MEMORY
#PMPI-PCT
\$785.00

Picker CT

Overall Dimensions
89" x 16"

COMFORT
#PPPI-NXR
\$640.00

MEMORY
#PMPI-NXR
\$916.00

Nuclear X-Ray


Overall Dimensions
71" x 24"

Transfer Boards



Allows safe and efficient transfer of patients for general imaging, MRI, and CT. Durable, semi-rigid, slick material makes positioning under the patient easy. Transfer boards are radiolucent can be used in the imaging field. Conveniently located 1½" x 4½" hand-holds allow optimal movement. Regular Boards are made of ½" thick HDPE (high density polyethylene). Regular board color is white. Transfer boards ship at OverSize rates.


REGULAR BOARDS



#TBS-2372	#TBN-1872	#TBB-3272
Qty 1-4: \$185.00 ea.	Qty 1-4: \$165.00 ea.	\$235.00 ea.
Qty 5-9: \$160.00 ea.	Qty 5-9: \$150.00 ea.	
Qty 10+: \$135.00 ea.	Qty 10+: \$130.00 ea.	

#TBH-2332
\$90.00 ea.

ANTI-STATIC BOARDS



Anti-static boards are made of UHMW (ultra high molecular weight polyethylene). Anti-static boards are black. Standard, Narrow, and Bariatric boards ship at OverSize rates.

cat.#	description	price ea.
TBS-2372AS	standard/half	\$305.00
TBN-1872AS	narrow	\$285.00
TBB-3272AS	bariatric	\$650.00
TBH-2332AS	half-size	\$150.00

Engraved Transfer Option

Personalize your transfer boards for your facility. Each board can be engraved with your name, phone number, etc. to make sure your boards come back to YOU!
(call for pricing)



Board Wall Rack



#TBR-2P
\$53.00

Two 6" pegs are angled to keep board from slipping. Mounting holes are on vertical centers. Steel with white enamel. Ships @ 4 lbs.



Lifting Straps

#TBS-ST2
\$30.00/pair

Convenient nylon straps loop through hand holds. Provides better grip. Provides extra length for better reach. Ships @ 1 lb.

Restraint Straps



Super strong 2" nylon webbing. Secure snap buckle closure.

cat.#	girth	description	price ea.
TBR-S72	72"	standard/half	\$35.00
TBR-N60	60"	narrow	\$30.00
TBR-B90	90"	bariatric	\$40.00

Portable Blanket Warmer Bags



Holds 2-3 Blankets

#BWS-03

\$475.00 ea.

Holds 4-5 Blankets

#BWS-05

\$515.00 ea.

Holds 6-7 Blankets

#BWS-07

\$545.00 ea.

All blanket warmer bags are IEC 60601-1 Medical Accessory Compliant.

Designed for radiology, surgery, outpatient, ER and EMS, these specially insulated blanket warmers keep warmed blankets warm. The durable and washable nylon exterior bag features a handy hook and loop closing flap. The removable interior heating elements keep the temperature at a steady 130° - 140°F when plugged into an optional, low voltage power supply (12vDC, <50w adapters sold separately). These lightweight and portable bags are designed and assembled in the USA and come with a 1-year warranty.

Optional AC Power Supply

#BWS-ACS

\$180.00 ea.



Optional 12v Cord Adapter

#BWS-12A

\$35.00 ea.



Blanket Warmer Cabinets

Designed for radiology departments, interventional suites, ER and OR, mammography, ultrasound, CT outpatient clinics and surgery centers, these stainless steel blanket warmers are adjustable from 98° to 180°F (37° to 82°C) with the help of an intelligent 4-digit LED controller. With no moving parts and an audio-visual temperature alarm, these heavy duty units are designed for safety and built to last. All three units come with hospital grade power cords and the two larger units include locking casters. Operational at either 120VAC or 220VAC, all three units are UL & CSA approved.

Countertop Unit

Short Cabinet

Tall Cabinet

Holds 7 bed sized blankets.
16" x 24" x 20"
70 lbs.

Holds 25 bed sized blankets.
22" x 29" x 40"
170 lbs.

Holds 50 bed sized blankets.
22" x 29" x 70"
350 lbs.

#BW-210

(call for price)

#BW-710

(call for price)

#BW-1400

(call for price)




Leg Support

Safely supports and positions patient leg. Easy to use Posi-Lock height adjustment. Lift height range: 14½" to 19". 6½" x 9½" calf support comes with nylon covered pad. All steel frame with sturdy 9" x 12" base. Shipping weight: 14 lbs.



#LGS-09
\$495.00



Bariatric Leg Support

Same as LGS-09, except: Extra large 15" x 18" base. Shipping weight: 19 lbs.

#LGS-B11
\$535.00



Provides true 90° Lateral View of Hip. Unobstructed Femoral Head and Neck.

Lateral Chest Support



Height adjustable bar provides stability for patients of almost every height. Helps reduce re-takes. Cast aluminum base easily mounts to the wall. Repositionable aluminum bar is comfortable to grip.


#RLC-36 (36" tall support)
\$1090.00
#RLC-48 (48" tall support)
\$1185.00

Caliper

Accurately measures body thickness to be x-rayed, up to 21" in depth. Calibrated in both inches and cm. Made of durable aluminum. 1" wide x 21¼" tall x 10⅛" leg. Ships @ 1 lb.



#CAL-21
\$57.00



Head Paddles

Ideal for a designated person to keep patient's head comfortable and steady. Rigid vinyl discs with comfort-grip handles. Coated foam cushions: 7" dia. x 1½" thick. Shipping weight: 4 lbs.


#RHI-CP
\$84.00



Head Immobilizer

Safely supports and positions patient head. Solid steel angles with soft, vinyl covered cushions. Blocks with handles weigh 15-lbs each. Cushions: 8"w x 6"h x 3"d. Shipping weight: 32 lbs.

#RHI-AB
\$250.00




(Edges may appear on images using low kVp.)

Radiolucent Immobilizer Strap

Comfortably secures patient position and eliminates motion. Helps prevent patient falls without limiting respiration or circulation. Ideal for MRI, CT, and Radiology dept. Available in 3 widths and 5 length for a variety of patients and tables. Strap is heavy gauge reinforced vinyl. Fast and easy hook and loop closure. Easy to remove in case of emergency. Fast and easy hook and loop closure. Easy to remove in case of emergency.



BLUE is the default color.	Available Colors	cat.#	width x length	price
	 Blue (default)	RDH-060	7" x 60"	\$75.00
		RDH-072	7" x 72"	\$85.00
		RDH-084	7" x 84"	\$105.00
		RDH-096	7" x 96"	\$125.00
	 Black	RDH-120	7" x 120"	\$160.00
		RDP-060	14" x 60"	\$115.00
		RDP-072	14" x 72"	\$125.00
		RDP-084	14" x 84"	\$145.00
		RDP-096	14" x 96"	\$165.00
	 Green	RDP-120	14" x 120"	\$210.00
		RDW-060	20" x 60"	\$145.00
		RDW-072	20" x 72"	\$155.00
		RDW-084	20" x 84"	\$175.00
		RDW-096	20" x 96"	\$195.00
		RDW-120	20" x 120"	\$240.00
	 Yellow			
	 Red			
	 Mint			

Colors are for Immobilizer Straps only. Hook and Loop is black.



Universal Restraint

Universal clamp fits most imaging tables. 6" x 52" or 10" x 72" nylon strap firmly and comfortably holds the patient to the table with hook and loop closure. Brackets are stainless steel for long life.



#RDS-12 6" x 52" \$160.00
#RDS-22 10" x 72" \$195.00

Elastic Extremity Immobilizer

6" wide strap is designed for superior stability. Strap extensions then cross over patient where they are secured with "hook" and "loop" closure. Three standard lengths; custom available.



cat.#	width x length	price
RDF-084	6" x 84"	\$50.00
RDF-096	6" x 96"	\$55.00
RDF-120	6" x 120"	\$110.00

Traction Straps

MANUAL #RDT-HN \$99.00

Straightening the legs easily applies traction. Hook and loop hand-hold straps provide secure grip.



Hand-Hold Straps

**Low
ers
shoulders
for cervical
imaging.**

AUTOMATIC #RDT-SH \$160.00

Primarily for elderly or unconscious patients. Rigid foot plate provides stable, steady support.



Buckle Straps

Mammography Biopsy Chair



Ideal for use as a procedural chair or treatment table. The chair's versatility allows it to be used for pre-induction, transport, treatment, and recovery. Easy controls from either side of the chair allow height adjustment for optimal patient positioning. Seat height ranges from 23.5" to 31.5". Premium 5" Tente casters make for easy maneuvering. Backrest and leg section adjust from 90° upright to full spine, and 10° Trendelenburg. Features fold-down handrails and 3" comfort mattress. Weight capacity: 325 lbs.

#CRH-MBC
(Call for Price)



Video Fluoroscopic Chair

Ideal solution for radiology, mammography, and transport. Specially designed back has removable section to accommodate the space confinements of large imaging equipment. Easy controls allow height adjustment and full swivel for perfect patient positioning. Seat height range: 23.5" to 31.5". Premium 5" Tente casters make for easy maneuvering. Backrest and leg section adjust from 90° upright to full supine, and 10° Trendelenburg. Features a removable back section and 3" comfort mattress.



with seat rotated for
increased imaging
capacity.

CRH-VIC
(Call for Price)





Defecogram Chair

Designed for use for patients with evacuation dysfunction. Defecography is a dynamic radiological examination that produces an anatomical depiction of the changes in the rectum and anal canal, and provides physicians with functional information about the movements of the pelvic floor at rest, during voluntary contraction, and with defecation. The chair allows radiography of a patient's bowel movement in radiology, rather than the more costly MRI/CT.

300-lb capacity unit includes locking casters, removable 10-qt basin, wrap around skirt for basin modesty, and a 9" raised heel support. Manufactured of Schedule 40 1½" healthcare-grade PVC, this chair is non-returnable and offers a one-year warranty. Made in USA.

Cat.# CRH-DFC
\$2500.00



Sterile Equipment Drapes

Take the guesswork out of selecting the correct sterile drape for your specific model C-Arm. Let us know which system you have and place your order. Our range of our drapes have been form-fitted for C-Arms to cover the image intensifier, x-ray tube, and flex arm, helping you eliminate the frustration and time of buying the wrong drape. Each drape is manufactured from 2 mil clear polyethylene and is latex-free. Each drape is designed to be durable and budget-friendly.



cat. #	description	qty	price
SED-4002	Elastic-end drape for Mini-C-arms: Officemate, Fluoroscope I and Fluoroscope III with 3/4" FOV	40	\$490.00
SED-4003	Elastic-end drape for Mini-C-arms: Fluoroscanner Premier/Encore, Fluoroscanner/Hologic Insight, Orthoscan, Lunar Orca, OEC with all FOV sizes, Xi Tex/Xi Scan	40	\$490.00
SED-1001	Standard C-arm drape, 42" x 74" (compatible with full-size C-arms: OEC, Ziehm, Philips, Comed)	30	\$240.00
SED-1002	Consumable pack includes: C-arm drape, 24" x 20" footswitch cover, 36" x 36" image intensifier & x-ray tube cover, and 5 clips. Fits OEC 9600/9800/9900.	20	\$480.00
SED-1003B	Banded bag for 9" image intensifier, 30" x 30"	20	\$160.00
SED-1003D	Dome cover w/ latex-free elastic, 18" deep	20	\$120.00
SED-8004F	Foot switch cover, 24" x 20"	100	\$300.00
SED-8003F	Foot switch cover, 17" x 15"	100	\$400.00





Mammography Exam Capes

Maintains patient modesty and confidence. Provides protection against cross-contamination. Comfortable, non-irritating materials. Eliminates laundry and replacement costs.



White Cape

Open front/back. Soft 3-ply tissue. Size: 30" x 21". Color: WHITE. 100 capes per case.

#GMC-21
\$42.00/case

Pink Cape

Open front/back. Soft 3-ply tissue. Size: 30" x 21". Color: PINK. 100 capes per case.

#GMC-22
\$50.00/case



Standard Exam Gown

Open front/back w/ties. Soft tissue/poly/tissue. Regular size: 30" x 42" Color: BLUE. 50 gowns per case.

#GXE-55
\$50.00/case



Extrasize Exam Gown

Open sides w/ties. Tissue/poly/tissue. Size: 40" x 50". Color: BLUE. 25 per case.

Size: XL
#GXS-57
\$125.00/case



Stretcher Sheets

Soft and water resistant. 2-ply tissue/poly. Size: 40" x 90". Color: WHITE/BLUE. 50 sheets per case.

#GSS-49
\$68.00/case



Pillow Cases

Soft and water resistant. 2-ply tissue/poly. Size: 21" x 30". Color: BLUE. 100 cases per case.

#GPC-23
\$62.00/case

True Wipes

A new standard for Lint-Free Wipes

Having passed the USP-797 clean environment testing, True-Wipes are ideal for lint-free cleaning of: CR/DR/CT/mammography screens, x-ray cassettes, work stations, PACS monitors, microscopes, lenses and general work surfaces. These highly absorbent 9" x 9" lint-free wipes can be used with screen cleaner, alcohol, bleach, and peroxide to clean and disinfect. Wipes are packaged in a resealable double bag to avoid contamination from typical box dispensers. Reinforces your vital QA standards.

#TW-200 (200 total wipes)

\$55.00/bag



SANI-CLOTH® AF3 Alcohol-Free Wipes

Contains an excellent range of effectiveness. Cleans, disinfects, deodorizes & fragrance-free. Alcohol content: ALCOHOL-FREE (less than 0.016%). Recommended Contact Time: 3-minutes or less. Commonly used for **MAMMOGRAPHY, ULTRASOUND and most other IMAGING EQUIPMENT**. PLUS daily use items like table pads, sandbags, coated sponges, transfer boards, aprons, shields, vinyl straps... you name it. Meets CDC, OSHA and CMS Tag F441 guidelines.

cat.#	wipe size	quantity	price
GWP-WAF	6" x 6¾"	tub of 160 wipes	\$14.50
GWP-WAFCS	6" x 6¾"	case of 12 tubs	\$170.00



SANI-CLOTH® BLEACH Wipes

Ideal for disinfecting high risk areas endemic with Multidrug-Resistant Organisms, Clostridium difficile spores, Norovirus, MRSA, and 50 other pathogenic microorganisms. Recommended Contact Time: 4-minutes. Compatible with a wide range of surfaces and equipment. Meets CDC, OSHA and CMS Tag F441 guidelines with a 1:10 bleach dilution.

cat.#	wipe size	quantity	case price
GWP-B12	6" x 10½"	12 tubs of 75 wipes	\$230.00
GWP-BB14	7½" x 15"	6 tubs of 65 wipes	\$215.00



EASY SCREEN™ Cleaning Wipes

Touch-Screen cleaner that "works at the speed of technology." Rapidly cleans dirt, grime, fingerprints and smudges from healthcare equipment. Helps prevent cross-contamination with a 70% isopropyl alcohol (IPA) solution. Commonly used for **Touch-Screens, Tablets, Lighting, and most other imaging equipment**. Effectively treats acrylic, polycarbonate, glass, stainless steel, and more.



#GWP-ESW
\$155.00 (12 tubs of 70 wipes)



Radiation Badge Storage Boards

Provides a central location for all badge related activity. Exposure record display area and note clip. Accommodates different styles and sizes of badges. Ring board is also available for Nuclear Medicine and PET.

NARROW BADGES

less than 1" wide



#FB-10N



#FB-20N



#FB-40N

SQUARE, HEX, WIDE,
OBLONG BADGES
1" to 2" wide



#FB-10W



#FB-20W



#FB-40W



HEX



SQUARE



SQUARE



RING



NARROW



WIDE

(badges available from your local provider)

cat. #	style	badges	board size	price
FB-10N	narrow	10	12" x 5"	\$79.00
FB-20N	narrow	20	12" x 9.25"	\$87.00
FB-40N	narrow	40	12" x 18"	\$135.00
FB-10W	square/hex	10	12" x 9.25"	\$105.00
FB-20W	square/hex	20	12" x 18"	\$135.00
FB-40W*	square/hex	40	12" x 18"	\$199.00
FB-10R	ring(nuc.med.)	10	12" x 5"	\$76.00

*(The FB-40W comes as two boards, each holding 20 badges.)


RING BADGES



#FB-10R

Merrill's Pocket Guide to Radiography, 13th Edition

Quick reference, soft cover spiral. Ideal for the clinical setting. Step-by-step illustrations for positioning 200 of the most common projections. Unique digital radiography reference radiographs. Handy technique chart for every projection, with manual and AEC techniques. 4 $\frac{3}{8}$ " x 7 $\frac{1}{4}$ ", ships @ 1 lb.



**Includes digital
collimation sizes!**

#BK-MPGR \$66.00




#1E \$22.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"
Without initials, order:
#4E @ \$14.00/set



#7E \$25.00 ea. set
R & L - 1"
2 initials - 3/8" or 3 initials - 1/4"
Without initials, order:
#5E @ \$16.00/set



#2E \$27.00 ea. set
RIGHT & LEFT - 1/2"
2-3 initials - 1/4"
Without initials, order:
#122E @ \$22.00/set



#6EB \$27.00 ea. set
R & L - 1/2", arrow - 1"
2-3 initials - 1/4"
Without initials, order:
#6EA @ \$20.00/set

Elite Style

Elite markers have the most DURABLE aluminum casing. RED for "R" markers and misc. words / BLUE for "L" markers. Provides the best in QC!

REFERENCE: Merrill's Atlas of Radiographic Positioning and Procedures, 9780323073219, text extract, "Each radiograph must include an appropriate marker that clearly identifies the patient's right (R) or Left (L) side. Medico-legal requirements mandate that these markers be present. R and L markers must be used with CR and DR digital imaging. The development of digital imaging and the use of CR and DR have enabled an environment in which the R and L markers can be placed on the image electronically. This is not recommended because of the great potential for error and legal implications. Anatomic markers should be placed on the CR cassette or DR table similar to screen-film cassettes."

**WORD and ARROW
Markers are on page 102!**

**Personal Touch Markers
are on pages 96-97!**

Marker Mates

A new, reusable sticky strip that adheres markers to receptors or equipment. Simply peel off the protective film and apply to the back of your markers. Each pair of Marker Mates will last approximately 30-days under normal use and care. No more need for tape when you can instantly mate your markers!

#MM-XRAY


1 pair
\$2.00

#MM-XRAY10

10 pair
\$20.00

#MM-XRAY25

25 pair
\$50.00



Position Indicator Markers

Beads Indicate Patient Position



#TA-13 \$32.00 ea. set
Elite Style R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"
Without initials, order:
#TA-14 @ \$27.00/set




#TA-1 \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"
Without initials, order:
#TA-4 @ \$13.00/set

Identifier Style

Identifier markers have high-impact polystyrene casing. Enhances quality control with fewer repeat exams. Provides basic, yet personalized image identification.

LOOK!
It's NEW!



#TA-LK \$22.00 ea. set
Look Right / Look Left
Arrow + "LOOK" - 1/2"



#TA-7 \$24.00 ea. set
R & L - 1"
2 initials - 3/8" or 3 initials - 1/4"
Without initials, order:
#TA-5 @ \$14.00/set



#TA-002 \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"
Without initials, order:
#TA-001 @ \$13.00/set



#TA-2 \$26.00 ea. set
Right & Left - 1/2"
2 or 3 initials - 1/4"
Without initials, order:
#TA-122 @ \$22.00/set

WORD and ARROW
Markers are on page 102!

Identifier Position Indicator Markers

Personal Touch Markers
are on pages 96-97!



#TA-17 \$37.00 ea. set
Identifier Style R & L - 7/8" Position
2 initials - 3/8" or 3 initials - 1/4"



#TA-16 \$32.00 ea. set
Identifier Style R & L - 7/8" Position
No initials

Personal Touch Colors & Shapes

Who said work had to be boring?



RED/BLUE SOLID
#TA-1 \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



RED/BLUE/PURPLE TIE-DYE
#PT-TR \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



PURPLE SOLID
#PT-CP \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



BROWN SOLID
#PT-BR \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



BLACK SOLID
#PT-BK \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



PINK SOLID
#PT-PK \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



TROPICAL FISH
#PT-FRB
\$22.00 ea. set

Ask about other Markers in Fun Colors!

Personal Touch Glitter

Have fun in what you do!



RED/BLUE GLITTER
#PT-GR \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



PURPLE GLITTER
#PT-GP \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



LIGHT BLUE GLITTER
#PT-GB \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



BURNT SIENNA GLITTER
#PT-GS \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



PINK GLITTER
#PT-GA \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



ORANGE GLITTER
#PT-GO \$20.00 ea. set
R & L - 5/8"
2 initials - 3/8" or 3 initials - 1/4"



SKULLS
#PT-SKW
\$22.00 ea. set

Ask about other Markers in Fun Colors!

Image Scaling Marker Device

Scaling a digital image is mandatory in order to assure the highest level of diagnostic accuracy. It is critical to know the difference between the size of the anatomy shown on a digital monitor as compared to the actual anatomy. An easy to use, easy to clean, easy to understand calibration tool is a must, and a Rollr Precision Marker is the answer!

The patent pending radiolucent sphere holder and a size indicator sphere come together to allow perfect visibility of the sphere on the image as taken from any direction. A millimeter marker on the sphere holder indicates the size of the sphere that is being used.

The suction cup base allows the unit to be firmly secured to a flat table or a vertical wall bucky board that is close to the patient. The small footprint of the suction cup does not interfere with the patient, the technologist, nor limit its use.

The flexible gooseneck arm is perfect for positioning the calibration sphere as close as possible to the critical anatomy. The arm is covered with a smooth, crevice-free surface that allows for easy sanitation.

The autoclavable unit comes with the suction stand and one calibration sphere. Additional spheres can be bought separately in packs of six (6).



Full unit with 25mm sphere
#A25DR \$645.00/ea.

Full unit with 25.4mm (1-inch) sphere
#A24DR \$645.00/ea.



Full unit with 30mm sphere
#A30DR \$695.00/ea.



Replacement 25mm sphere (set of 6)
#S25DR \$240.00/set


Replacement 25.4mm sphere (set of 6)
#S24DR \$240.00/set

Replacement 30mm sphere (set of 6)
#S30DR \$575.00/set



CM Ruler Markers

The personalized R & L marker sets or the stand alone 3cm ruler are specifically designed for Teleradiology, CR and DR Imaging, where the image-size does not always match the life-size anatomy. Place a marker on the Receptor Plate and capture the radiographed anatomy. The image of the unique 3cm ruler gives an immediate visual and measureable cue as to the size of the screen image compared to the life size anatomy. Each marker has a lead 3cm ruler. Color coded: L-Blue, R-Red. High density lead characters.



#1DR \$50.00 ea. set
Ruler Style R & L - $\frac{3}{4}$ "
2 or 3 initials - $\frac{1}{4}$ "



#3DR \$30.00 ea.
Ruler only



#2DR \$50.00 ea. set
Ruler Style R & L - $\frac{1}{2}$ "
2 or 3 initials - $\frac{1}{4}$ "

Wafer-Thin Digital Imaging Markers

Perfect for equipment (CR, DR, etc.) with very thin bucky clearance. Exceptionally thin markers: $< \frac{3}{32}$ " (2mm). Copper filtered for a broad range of KV. Elite style and color-coded.

**Super Thin for
Slim Bucky Clearance**

**Other Wafer-Thin markers
are on pages 100 & 102!**



#WTS-DR \$30.00 ea. set
Ruler Style R & L with
3cm spacing

Radi-O-Paque Extremity Rulers

Sturdy $\frac{1}{4}$ " wide x $\frac{1}{8}$ " thick. High density lead letters. Brilliant YELLOW color. Spaced in 0.5cm increments.


45cm
#CPT-R45
\$175.00/ea.



115cm
#CPT-R115
\$275.00/ea.



Enlarged View



Wafer Thin

Super Thin for Slim Bucky Clearance

Perfect for newer equipment (CR, DR, etc.) with very thin bucky clearance. Exceptionally thin markers: $< \frac{3}{32}$ " (2mm). Copper filtered for a broad range of KV. Elite style and color-coded.



#WTS-6B \$32.00 ea.set
R & L - $\frac{3}{8}$ " 2 initials - $\frac{3}{16}$ "

Without initials (with arrow), order:
#WTS-6A @ \$21.00/set



#WTS-6C \$32.00 ea.set
R & L - $\frac{3}{8}$ " 2 initials - $\frac{3}{16}$ "

Without initials, order:
#WTS-6D @ \$21.00/set



#WTS-1 \$32.00 ea.set
R & L - $\frac{3}{8}$ " 2 initials - $\frac{3}{16}$ "

Without initials, order:
#WTS-4 @ \$21.00/set



(See Wafer Thin Words on p.98)



#WTW-A3 ($\frac{1}{2}$ " arrow)
\$10.00

Oncology, CT and CR

Super Thick for High Radiation

$\frac{1}{4}$ " super thick, high-density lead letters. Designed for Oncology and Industry, and other applications. Elite style markers, color-coded and Techno-Aide tough.



#RO-107B \$32.00 ea.set
R & L - $\frac{3}{4}$ " 2 initials - $\frac{3}{8}$ "

Without initials, order:
#RO-107A @ \$24.00/set

Vinyl & Epoxy encased, for longer life!



Ex.: ROW-23 $\frac{3}{4}$ " letters ($\frac{1}{4}$ " thick)

The Oncology markers are $\frac{1}{4}$ " thick and will NOT fit in most low clearance bucky trays.

cat.# (3/4")	description	price
ROW-01	single letter	\$12.00
ROW-23	2-3 letters	\$15.00
ROW-45	4-5 letters	\$18.00
ROW-67	6-7 letters	\$22.00
ROW-80	8+ letters	\$25.00

Copper Filtered

Copper filtered for High KV. High density lead letters. Available in ELITE or IDENTIFIER. Color-coded, impact resistant. (Should not be used for mammography.)

**Filtered for
High KV Exams**



#COP-1 \$42.00 set ELITE
(2 sets of 1E markers: 1 set with copper filter, 1 set without.)



#COP-2 \$40.00 set ELITE
R & L - $\frac{1}{2}$ ", 2-3 initials - $\frac{1}{4}$ "
(Filtered & unfiltered on the same marker.)

Standard Thin Clipper

Constructed of high strength Kydex™. Extra thin profile, double-sided, 1 7/8" wide. R and L: 3/4", with or without initials. 2 Initials: 3/8"; 3 initials or words: 1/4".

#CPP-40
(with initials)
\$28.00

#CPP-42
(no initials)
\$25.00



Side View



Stencil Clipper

Plated Steel Clipper

#XSCP-10
\$15.00 ea.
Stenciled non-
personalized R & L.



Gridcap Thin Clipper

Designed to fit DR, CR receptors with gridcap. Constructed of high strength Kydex™. Extra thin profile, double-sided, 1 7/8" wide. R and L: 3/4", with or without initials. 2 Initials: 3/8"; 3 initials: 1/4".

#CPP-60 (with initials)
#CPP-61 (no initials)
\$38.00 ea.



Elite Clippers

Both R & L are on one double clip. Durable Stainless Steel clips. Color-coded: R side-Red / L side-Blue.

Double Personalized



#CP-35
\$30.00 ea.
(2-3 initials)

Portable & Position Double Clippers



#CP-36
(with initials)
\$37.00

#CP-37
(no initials)
\$30.00



#CP-13
(with initials)
\$37.00

#CP-14
(no initials)
\$30.00

ELITE Words (aluminum casing)

38" Elite



12" Elite



**Personalized
Marker Sets
are on
pages 94-97**

IDENTIFIER Words (polystyrene casing)

38" Identifier



12" Identifier



description	example	cat.# (3/8")	price	cat.# (1/2")	price
single letter	3	2E-Y	\$9.00	2E-Z	\$10.00
2-3 letters	EXT	3E-Y	\$11.00	3E-Z	\$13.00
4-5 letters	25 MIN	4E-Y	\$13.00	4E-Z	\$16.00
6-7 letters	STRESS	5E-Y	\$15.00	5E-Z	\$19.00
8+ letters	PORTABLE	6E-Y	\$18.00	6E-Z	\$22.00

description	example	cat.# (3/8")	price	cat.# (1/2")	price
single letter	3	2TA-Y	\$8.00	2TA-Z	\$9.00
2-3 letters	EXT	3TA-Y	\$10.00	3TA-Z	\$11.00
4-5 letters	25 MIN	4TA-Y	\$12.00	4TA-Z	\$13.00
6-7 letters	STRESS	5TA-Y	\$14.00	5TA-Z	\$15.00
8-15/8-10	PORTABLE	6TA-Y	\$17.00	6TA-Z	\$18.00

WAFER THIN Words

(1/4" letters x <3/32"thin)



cat.# (3/4")	description	price
WTW-01	single letter	\$13.00
WTW-23	2-3 letters	\$16.00
WTW-45	4-5 letters	\$19.00
WTW-67	6-7 letters	\$22.00
WTW-80	8+ letters	\$25.00

ONCOLOGY Words

(3/4" letters x 1/4"thick)



cat.# (3/4")	description	price
ROW-01	single letter	\$12.00
ROW-23	2-3 letters	\$15.00
ROW-45	4-5 letters	\$18.00
ROW-67	6-7 letters	\$22.00
ROW-80	8+ letters	\$25.00

Wall Boards

Walnut grain w/brass hooks.



cat.#	description	price
HK-32	32 Hook Board	\$100.00
HK-16	16 Hook Board	\$80.00
HK-8	8 Hook Board	\$70.00
HK-50	50 Hook Board	\$135.00

These are just a few suggestions from all possible word choices:

2-3 letters	4-5 letters	6-7 letters	8+ letters
AP	FLAT	AP ERECT	ANTERIOR
PA	FLEX	CLOSED	EXTENSION
ANT	INSP	INT ROT	INFUSION
EXP	LEFT	MOBILE	PORTABLE
EXT	OPEN	P EVAC	POST EVAC
INT	POST	POST OP	POST VOID
KUB	TOMO	RECUMB	STANDING
LAO	DECUB	STRESS	VERTICAL W/WEIGHT
LPO	ERECT	SUPINE	WEIGHT BEARING
MAG	HOURS	DELAYED	INJECTION
MIN	PRE OP	FLEXION	INVERSION
OBL	PRONE	LATERAL	SEMI ERECT
RAO	PUREE	MINUTES	VALSALVAR
LOA	RIGHT	OBLIQUE	WT BEARING
25	SCOUT	SUNRISE	COMPARISON
600	60 MIN	UPRIGHT	DECUBITUS



#A-3 \$9.00 each
Arrow - 1" Elite Style



#TA-A3 \$7.00 each
Arrow - 1" Identifier Style

General Marker Set



Identifier Set shown (Board: 12" x 15")
3/8" letters / 1/2" numbers
Elite Identifier
#GME-60 #GMI-65
\$575.00 \$525.00



250 Stickers
per box



Cover Tab Stickers

Ideal for coding mismarked or unmarked film. Color coded: L-Blue / R-Red. Self-adhesive backing. 1¼" diameter with square for initials. Available in greater quantity discounts!

Take Advantage of Quantity Discounts!

cat.# (¾")	1-3 boxes	4-9 boxes	10+ boxes
TAB-R	\$22.00	\$16.00	\$14.00
TAB-L	\$22.00	\$16.00	\$14.00

Marking Pens

Highlights points of interest on a film. Wax markers resist smudging, but are not permanent. NOT to be used in place of lead markers. NOT a legally acceptable marking method.

Order each item per number of packs.

Example: 2ea. #MPY-11

(will be shipped with 10 yellow pens.)



Example of the
MPA-19
Assorted Pack



PENS (5 per pack)

cat.# (¾")	color	price/pk
MPY-11	Yellow	\$30.00/pk
MPR-12	Red	\$30.00/pk
MPW-13	White	\$30.00/pk
MPA-19	2-Y, 2-R, 1-W	\$30.00/pk

REFILLS (72 per pack)

cat.# (¾")	color	price/pk
MPY-21	Yellow	\$20.00/pk
MPR-22	Red	\$20.00/pk
MPW-23	White	\$20.00/pk

Refills come in pkgs of 72 pieces.



**CHOOSE YOUR
OWN COLOR!**

A tool pouch for your markers!



Marker Keeper

Convenient four-pouch vinyl holder for markers. Handy clip attaches the Marker Keeper to belt, lapel, pocket, etc. Keeps markers in easy reach for quick retrieval.

cat.# (¾")	color	price
MKR-11	Red	\$9.90
MKB-12	True Blue	\$9.90
MKG-13	Green	\$9.90
MKY-14	Yellow	\$9.90
MKC-15	Clear	\$9.90
MKA-19	*Assorted (4)	\$35.00

(*Assorted set does NOT include Clear.)



Skin Markers for General X-Ray, Mammography, CT and Oncology

GENERAL IMAGING: used to identify areas of pain or trauma, soft tissue masses, to document points of entry and exit sites, or mark prominent nipples for identification of nipple shadows on chest x-rays.




Skin Dots
Artifact/nipple identifier.
Helps eliminate repeat exams. Lead-free. Easy-Peel tab.

cat.#		BB size
SDM-BB15		1.5mm
SDM-BB20		2.0mm
SDM-BB25		2.5mm
1-5 bxs	6-14 bxs	15+ bxs
\$66.00	\$63.00	\$60.00
(100 per box)		




Skin Wings
Artifact/nipple identifier.
Extra support for larger, more prominent nipples. Lead-free. Easy-Peel tab.

cat.#		BB size
SDM-W15		1.5mm
SDM-W20		2.0mm
SDM-W25		2.5mm
1-5 bxs	6-14 bxs	15+ bxs
\$70.00	\$67.00	\$64.00
(100 per box)		



Skin Dots (Rad/Onc)
Lead-free artifact identifier. Larger BBs for Radiology and Oncology. Easy-Peel tab.

cat.#		BB size
SDM-BB30		3.0mm
SDM-BB40		4.0mm
SDM-BB50		5.0mm
1-5 bxs	6-14 bxs	15+ bxs
\$125.00	\$120.00	\$115.00
(50 per box)		




Skin Arrow (Rad/Onc)
Lead-free arrow pointer. Highlights areas of interest for Radiology and Oncology. Easy-Peel tab.

cat.#		BB size
SDM-A1		7.0mm
1-5 bxs	6-14 bxs	15+ bxs
\$105.00	\$101.00	\$97.00
(100 per box)		



Easy-Peel, adhesive-free tabs allow comfortable removal of the marker.



Skin Line/Scar Marker
Lead-free scar identifier. Line contours to complex scar tissue. Flexible adhesive tabs; no scissors required.

cat.#		BB size
SDW-SS3		1.0mm
1-5 bxs	6-14 bxs	15+ bxs
\$85.00	\$81.00	\$77.00
(120"/305cm per box)		

Economy Nipple-Artifact Markers

Encapsulated 2mm lead ball. Used to mark nipple, lesion, scar, artifact, etc. Neutral color, peel-off, self-adhesive backing. Non-sterile. 100 markers mounted on ten convenient strips.
(This product contains lead.)

#LNA-20
\$50.00



2D/3D Breast Tomo Markers

Used for Digital Breast Tomosynthesis.
Low density markers specifically designed to show up clearly when “in slice” and minimize the “slinky” artifact when they’re “out of slice.”

NEW!



Mole Marker

Lead-free mole marker.
Easy-Peel tab.

cat.# SDG-MSD3
\$70.00 / box of 100
Hole Size: 1/2"



Palpable Mass Marker

Lead-free Mass skin marker.
Easy-Peel tab.

cat.# SDM-TSD3
\$70.00 / box of 100



Non-Palpable Mass Marker

Lead-free arrow skin marker.
Easy-Peel tab.

cat.# SDM-ASD3
\$70.00 / box of 100



Skin Line/Scar Marker

Lead-free line skin marker.
Repeated adhesive tabs allow accurate bending to scar outline.

cat.# SDW-SSD3
\$85.00 / roll of 120"
Line: 1mm

Digital Mammography Skin Markers



Digital Dot

Lead-free nipple identifier.
Highlights areas of interest for digital and analog imaging. Easy-Peel tab.

cat.#	BB size		
SDM-BB23	2.3mm		
1-5 bxs	6-14 bxs	15+ bxs	
\$66.00	\$63.00	\$60.00	
(100 per box)			

CT Skin Markers



CT Dot: 2.3mm

Lead-free 2.3mm marker.
Highlights areas of interest for Radiology and CT. Easy-Peel tab.

cat.#	BB size		
SDM-CT1	2.3mm		
1-5 bxs	6-14 bxs	15+ bxs	
\$89.00	\$84.00	\$80.00	
(50 per box)			



CT Dot: 4.0mm

Lead-free 4.0mm marker.
Highlights areas of interest for Radiology and CT. Easy-Peel tab.

cat.#	BB size		
SDM-CT2	4.0mm		
1-5 bxs	6-14 bxs	15+ bxs	
\$89.00	\$84.00	\$80.00	
(50 per box)			

Multi-Modality Skin Markers



CT / MRI

No center hole; appears as bright object on CT and MRI scans. Also highlights areas of interest for radiology and mammography.

cat.#	BB size		
MM3007	solid, no hole		
1-5 bxs	6-14 bxs	15+ bxs	
\$240.00	\$235.00	\$230.00	
(50 per box)			

Lead Characters

If you have any questions about the size characters you need or are looking for additional sizes and styles, please call our customer service department, or mail us a sample so that we may send you exactly what you're looking for.



UNmounted Lead Characters

style	(3/16")	(1/4")	(3/8")	(1/2")	(5/8")	(1" R&L only)	(1" arrow)
FLAT FACE	XTUF-36 .50	XTUF-14 .50	XTUF-38 .50	XTUF-12 .60	XTUF-58 .60	XTUF-01 1.00	X
DEEP BLOCK	X	XTUB-14 .90	XTUB-38 1.00	XTUB-12 1.00	(see Flat Face)	X	X
ARROWS	X	XTUA-14 1.00	X	XTUA-12 .80	X	X	XTUA-01 1.00

Flat Face | Deep Block | Arrow



"Loose" lead letters & numbers.

Individual Mounted Characters (Flat Face)

cat.#	lead/plastic size	price	comments
XTSF-3658	3/16" on 5/8"	.80	Fits: LLH-58, LLH-2, LLH-3, Veri-Thin
XTSF-3634	3/16" on 3/4"	.80	Fits: Legal (regular)
XTSF-1458	1/4" on 5/8"	.80	
XTSF-1434	1/4" on 3/4"	.80	
XTSF-1478	1/4" on 7/8"	.80	Fits: Eclipse, Perfection
XTSF-1234	1/2" on 3/4"	.80	
XTIF-3658	3/16" on 5/8"	.80	Fits: Ideal (special offset mounting)



Characters are mounted on CLEAR flat plastic tabs to fit a variety of letter holders.

Flat Surface Mounted Words (Flat Face)

description	example	cat.#	price	cat.#	price	cat.#	price	cat.#	price
single letter	3	XVFW-0136	...1.00	XVFW-0114	...1.00	XVFW-0112	...1.00	XVFW-0101	.2.00
2-3 letters	EXT	XVFW-2336	2.00	XVFW-2314	2.00	XVFW-2312	2.00	XVFW-2301	4.00
4-5 letters	25 MIN	XVFW-4536	4.00	XVFW-4514	4.00	XVFW-4512	4.00	XVFW-4501	7.00
6-7 letters	STRESS	XVFW-6736	6.00	XVFW-6714	6.00	XVFW-6712	6.00	XVFW-6701	9.00
8+ letters	PORTABLE	XVFW-8036	7.00	XVFW-8014	7.00	XVFW-8012	7.00	XVFW-8001	11.00
Lead/Plastic Sizes		3/16" on 3/4"		1/4" on 3/4"		1/2" on 3/4"		1" on 1 1/2"	



Words made from Mounted Letters glued to the surface of a single CLEAR vinyl strip.

Lead Letter Case

Injection Molded Box. 36 scoop-bottom compartments (for A-Z and 0-9). Side compartment for holder.

#XSLC-36


\$44.00



Ideal Cushion Marker

Double channel with cushions; imprint in the upper section of each channel.

#XSIC-21 \$55.00 ea.
Holder only. 3½" x 15/16"




#XSIC-22 \$150.00 ea set

Complete Set Includes:
1ea. holder
2ea. A-Z
4ea. 0-9
2ea. A, E, I, O, U
1ea. 36-letter case
(Replacement characters:
3/16" on 5/8" plastic tabs, offset.
Reorder #XTIF-3658)

Eclipse Cushion Marker

Single channel with cushion; imprint above and below the centerline.

#XSEC-11 \$49.00 ea.
Holder only. 3½" x 15/16"



#XSEC-12 \$100.00 ea set

Complete Set Includes:
1ea. holder
2ea. R, L
4ea. 0-9
1ea. 36-letter case
(Replacement characters:
1/4" on 7/8" plastic tabs.
Reorder #XTSF-1478)

Legal Marker

Double aluminum channels; imprint between the two channels.

#XSLG-21 \$50.00 ea.
Holder only. 3½" x 15/8"



#XSLG-22 \$130.00 ea set

Complete Set Includes:
1ea. holder
2ea. A-Z
4ea. 0-9
2ea. A, E, I, O, U
1ea. 36-letter case
(Replacement characters:
3/16" on 3/4" plastic tabs;
Reorder #XTSF-3634)

Veri-Thin Legal

Double aluminum channels; imprint down and in between the two channels. Great for Bucky Trays!

#XSVL-21 \$50.00 ea.
Holder only. 3½" x 113/16"



#XSVL-22 \$130.00 ea set

Complete Set Includes:
1ea. holder
2ea. A-Z
4ea. 0-9
2ea. A, E, I, O, U
1ea. 36-letter case
(Replacement characters:
3/16" on 5/8" plastic tabs.
Reorder #XTSF-3658)

Perfection

Single aluminum channel; imprint above and below the channel.

#XSPF-11 \$40.00 ea.
Holder only. 3½" x 15/16"



#XSPF-12 \$100.00 ea set

Complete Set Includes:
1ea. holder
2ea. R, L
4ea. 0-9
1ea. 36-letter case
(Replacement characters:
1/4" on 7/8" plastic tabs.
Reorder #XTSF-1478)

Illuminated "In Use" Door Signs

Lighted sign focuses attention on exam room status. Enhances Radiation Safety by suggesting the risk of accidental entry into an active exam room. NOT a plug-in sign; must be hardwired @ 120VAC. Each sign measures 11½" x 2½" x 7¾"h. Shipping weight of each sign: 4 lbs.



cat.# (38")	description	price
SLU-XR	X-RAY IN USE	\$169.00
SLU-CT	CT IN USE	\$169.00
SLU-MR	MRI IN USE	\$169.00



L.E.D. "X-RAY IN USE" Sign

Power-saving LED light technology meets newer compliance demands. Hardwired sign is 11½" wide x 7¼" high. Shipping weight of each sign: 4 lbs.

cat.# (38")	description	price
SLU-XL	LED: X-RAY IN USE	\$199.00




Illuminators (View Boxes)

STANDARD

High quality, cost effective illuminators. Includes upper and lower Film Grips. Internal reflectors offer consistent illumination from two 15-watt fluorescent lamps per viewing area. Use as stand alone or wall mounted. Average light intensity is 2000 nits. 6 foot power cord. UL listed.

LED POWER SAVERS

Use **70% LESS POWER** than Standard Illuminators! High quality illuminators include upper and lower Film Grips. Internal LED array offers consistent illumination from a matrix of lights per viewing area. Use as stand alone or wall mounted. Average light intensity is 4250 nits. 6 foot power cord with hospital grade plug. UL listed.



single standard
#VBS-10 \$240.00

single LED
#VBS-10LED \$1100.00
14"w x 3½"d x 21¾"h
Shipping weight: 17 lbs.



double standard
#VBD-20 \$475.00

double LED
#VBD-20LED \$2125.00
28"w x 3½"d x 21¾"h
Shipping weight: 34 lbs.

Mammography Illuminator

Light source specifically designed to meet ACR Mammography Viewing requirements. Unit offers two (2) side-by-side film viewing areas for 24 x 30 film. Light intensity: 4,200 nits. Upper and Lower Film Grips. White powder coat finish. Includes four (4) 15-watt fluorescent lamps. UL listed with 5-foot power cord. Overall size: 20½"w x 3½"d x 16¾"h. Shipping weight: 17 lbs.




#VBM-24
\$530.00




I.D. Cards - Standard

Standard Cards measure 3" x 5". Each card is imprinted with standard legend. Each pack contains 1,000 cards. Shipped from stock.



#IDS-02M \$46.00/pk
Blank with 2-Line Legend



#IDS-03M \$46.00/pk
Blank with 3-Line Legend

I.D. Cards - Custom

Custom Cards measure 3" x 5" with 1" x 3" imprint window (GE and Kodak compatible cards are 3¼" x 7¾" with 1" x 2½" imprint window). Each card is imprinted with: Your Custom Imprint (0, 1, 2 or 3 lines); Custom Legend (0, 1, 2 or 3 lines); (if not specified, standard legend will be used). Each pack contains 1,000 cards.

CUSTOM
CARD
BUILDER

Number of Lines for Your Imprint
Number of Lines for Legend
Style (see examples)

Cat.# IDC-
0 0 M
1 1 C
2 2 D
3 3 Q

style	imprint location	1000 cards	5000 cards*
M	Middle Top	\$55.00/M	\$49.00/M
C	Top Corner	\$55.00/M	\$49.00/M
D	Double Middle	\$63.00/M	\$57.00/M
Q	GE or Kodak (compatible)	\$65.00/M	\$59.00/M


(*5000 quantity is priced per 1000 cards.)

(When ordering, a faxed diagram is most helpful.)

Artwork may not be an exact match to that of other suppliers.


Custom Cards are non-returnable.

If not specified, standard legend will be used.




Example:
#IDC-32M


Most Popular Card!



Example:
#IDC-23C



Example:
#IDC-22D



Example:
#IDC-22Q

I.D. Printer

Fully automatic exposure for perfect image prints. Adjustable exposure timer for most film types. Film positioning guides for middle or corner ID cards. Durable, all-metal housing. Plug-in model only; 110/120 VAC. 8"w x 8"d x 3"h. Shipping weight: 6 lbs.



#IDP-33
\$365.00



Accommodates Middle & Corner Cards

Pin Adapter can be set for use with middle or corner cards.

Standard X-RayFilm Jackets

Choose one or two pockets. Non color-coded (see label selection for sticker styles and colors). Size: 14.5” x 18.25”, including 3/4” scored self expansion. Made of 11pt. manila stock. Black ink copy includes blocks for color-code labels and black loan indicator. Offset double thumbcut for easy access to contents. Double pocket version includes 10” x 15.75” report pocket with 3/4” self expansion. Open side (opens along upper edge).



single pocket

cat#	qty./ctn	wt./ctn	1-9 ctns	10-24 ctns	25+ ctns
JJL-SPS	100	22#	\$110.00	\$105.00	\$100.00



double pocket

cat#	qty./ctn	wt./ctn	1-9 ctns	10-24 ctns	25+ ctns
JJL-DPS	100	25#	\$170.00	\$165.00	\$160.00

Multi-use Jackets

Universal design allows many uses with one folder. Size: 9 9/16” x 12 3/4”, including 3/4” scored self expansion. Made of 11pt. manila stock. Universal copy includes check boxes for different exams. Offset double thumbcut for easy access to contents. Open side (opens along upper edge). (May require long lead-time.)

cat#	qty./ctn	wt./ctn	1-9 ctns	10-24 ctns	25+ ctns
JJM-SPU	200	19#	\$280.00	\$275.00	\$270.00



Color-Code Labels

Features permanent adhesive and release coated paper backing.
Superior flexibility for long life. Other color systems available.


cat#	description	qty./ctn	1-9 ctns	10-24 ctns	25+ ctns
see list	any single	1000/rl	\$33.00	\$31.00	\$29.00
JLL-CNS	0-9 set	10 rolls	\$330.00	\$310.00	\$290.00

cat#	description	qty./ctn	1-9 ctns	10-24 ctns	25+ ctns
see list	any single	1000/rl	\$27.00	\$25.00	\$23.00
JLS-CCS	set of 10	10 rolls	\$270.00	\$250.00	\$230.00

large digit reverse block



(shown actual size: 1 7/8" x 1 7/8")



(shown actual size: 1 7/8" x 1 7/8")

large year

cat#	year/color	qty./ctn	1-9 ctns	10-24 ctns	25+ ctns
JLL-Y17	17/brown	1000/rl	\$29.00	\$27.00	\$24.00
JLL-Y18	18/pink	1000/rl	\$29.00	\$27.00	\$24.00
JLL-Y19	19/green	1000/rl	\$29.00	\$27.00	\$24.00

Cat.#

Color

Cat.#

JLL-ORD	red	JLS-CRD
JLL-1GY	gray	JLS-CGY
JLL-2BL	blue	JLS-CBL
JLL-3OR	orange	JLS-COR
JLL-4PR	purple	JLS-CPR
JLL-5BK	black	JLS-CBK
JLL-6YL	yellow	JLS-CYL
JLL-7BN	brown	JLS-CBN
JLL-8PK	pink	JLS-CPK
JLL-9GN	green	JLS-CGN
JLL-CNS	(set of 10)	JLS-CCS
--	aqua	JLS-CAQ
--	gold	JLS-CGD
--	wallet	JLS-CNW
--	white	JLS-CWT

solid color bar label



(shown actual size: 1 3/4" x 1 1/2")



(shown actual size: 1 3/4" x 1 1/2")


small year

cat#	year/color	qty./ctn	1-9 ctns	10-24 ctns	25+ ctns
JLS-Y17	17/brown	1000/rl	\$20.00	\$18.00	\$16.00
JLS-Y18	18/pink	1000/rl	\$20.00	\$18.00	\$16.00
JLS-Y19	19/green	1000/rl	\$20.00	\$18.00	\$16.00

Mammography Jackets

Pink border identifies as separate from master jacket.
Size: 11" x 14", including 3/4" scored self expansion. Pocket size: 8 1/4" x 11". Illustrated back allows for Dx & patient communication. Offset double thumbcut for easy access to contents. Open side (opens along upper edge).

	cat#	qty./ctn	wt./ctn	1-9 ctns	10-24 ctns	25+ ctns
PINK Border	JJM-PPM	200	30#	\$235.00	\$230.00	\$225.00
BLUE Border	JJM-BBM	200	30#	\$235.00	\$230.00	\$225.00



Category Inserts

Locate specific examinations at a glance using solid-color and printed border stock color systems. Made of 32lb. white stock. Color coding by category permits quick location of films. "Open side" category name is printed on both sides of thumbcut; "open end" category name is printed on top and bottom of insert for easy identification. All inserts are available in open side (long edge) and open end (short edge) styles.

Choose from two sizes:

large: 14.25" x 17.5"

medium: 10.5" x 12.5"



Large Size

size (inches)	qty./ctn	wt./ctn	1-9 ctns	10-19 ctns	20+ ctns
14.25 x 17.5	250	20#	\$135.00	\$130.00	\$125.00
open end#	open side#	category	stock	copy	border
JCL-ABE	JCL-ABS	Abdomen	White	Black	Lt. Green
JCL-ANE	JCL-ANS	Angio	White	Black	Purple
JCL-BEE	JCL-BES	BE*	White	Black	Gold
JCL-BNE	JCL-BNS	Bone	White	Black	Yellow
JCL-CHE	JCL-CHS	Chest*	White	Black	Brown
JCL-CTE	JCL-CTS	CT*	White	Black	Dk. Green
JCL-EXE	JCL-EXS	Extremities	White	Black	Yellow
JCL-GIE	JCL-GIS	GI*	White	Black	Gray
JCL-GUE	JCL-GUS	GU*	White	Black	Orange
JCL-HNE	JCL-HNS	Head & Neck	White	Black	Lt. Blue
JCL-IVE	JCL-IVS	IVP	White	Black	Salmon
JCL-MME	JCL-MMS	Mammography*	White	Black	Dk. Blue
JCL-MRE	JCL-MRS	MRI*	White	Black	Rust
JCL-NME	JCL-NMS	Nuclear Med	White	Black	Red
JCL-SEE	JCL-SES	Spine*	White	Black	Black
JCL-USE	JCL-USS	Ultrasound	White	Black	Purple

*Includes anatomical diagram on back.

size (inches)	qty./ctn	wt./ctn	1-9 ctns	10-19 ctns	20+ ctns
14.25 x 17.5	250	22#	\$180.00	\$175.00	\$170.00
open end#	open side#	category	stock	copy	border
JCL-NCE	JCL-NCS	(no category)	White	Black	Dk. Blue

Medium Size

size (inches)	qty./ctn	wt./ctn	1-9 ctns	10-19 ctns	20+ ctns
10.5 x 12.5	250	25#	\$135.00	\$130.00	\$125.00
open end#	open side#	category	stock	copy	border
JCM-BDE	JCM-BDS	Bone Density	White	Black	Yellow
JCM-MME	JCM-MMS	Mammography	White	Black	Dk. Blue
JCM-NME	JCM-NMS	Nuclear Med	White	Black	Red
JCM-USE	JCM-USS	Ultrasound	White	Black	Purple

Negative Preservers

Available in three sizes to accommodate all films. All include stock copy lines for: file number, name, date, doctor and remarks. Made of 32lb. gold kraft stock. Thumbcut for easy access to films. Open end (opens along short side edge). Red and green available in large size, upon request.

cat#	size	qty./ctn	wt./ctn	1-7 ctns	8-14 ctns	15+ ctns
JNS-OEK	8.5 x 10.5	500	25#	\$155.00	\$150.00	\$145.00
JNM-OEK	10.5 x 12.5	500	25#	\$175.00	\$170.00	\$165.00
JNL-OEK	14.5 x 17.5	250	25#	\$135.00	\$130.00	\$125.00



CD/DVD Sleeves

Standard 5" x 5" CD/DVD protectio sleeves. Made of 60lb. white offset paper stock. Window made of clear polyester.

cat#	qty./ctn	wt./ctn	1-7 ctns	8-14 ctns	15+ ctns
JCD-CSW	1000	14#	\$139.00	\$136.00	\$130.00



CD/DVD Mailers

Features self-adhesive closure and tear-tape tab opening. Made of 20pt. coated white board. Economy; no liner.

cat#	qty./ctn	wt./ctn	1-7 ctns	8-14 ctns	15+ ctns
JCD-CDX	200	15#	\$75.00	\$70.00	\$65.00



Cardboard Insert Mailers

Available in two sizes. Made of 32lb. brown kraft stock. Open end, string & button closure. Chipboard insert for rigidity. Printed with bold black stock copy.

cat#	size (inches)	qty./ctn	wt./ctn	1-7 ctns	8-14 ctns	15+ ctns
JMM-SBC	11 x 13	50	8#	\$80.00	\$75.00	\$70.00
JML-SBC	15 x 18	50	12#	\$100.00	\$95.00	\$90.00



Green Diamond Border Mailers

Available in two sizes. Made of 11pt. manila stock. Open end, string & button closure. Green diamond border with red. "X-Ray Film, please do not bend."

cat#	size (inches)	qty./ctn	wt./ctn	1-7 ctns	8-14 ctns	15+ ctns
JMM-SBD	11 x 13	50	7#	\$100.00	\$95.00	\$90.00
JML-SBD	15 x 18	50	10#	\$100.00	\$97.00	\$95.00



E-Z X-Ray Storage Boxes

Allows for the efficient storage of x-ray film files. (Sold in case lots.) Strong and durable with reinforced handles. One piece, easy fold assembly. Can be stacked without crushing. Ideal for archiving files, temporary storage, or moving files. Accommodates jackets up to 14" x 17". Upright: lid opens up. Side down: lid opens forward. (printed legend provided both ways.) Sized for OSHA lift load capacity. Each box O.D.: 6.25" x 16" x 19.5". Ships as 30 lbs. per case.

6 cases or less via FedEx/UPS.

7+ cases via motor freight.



Cat.#ETX-70
(15 boxes per case)

# cases	price
1-3	\$145.00
4-7	\$140.00
8-15	\$135.00
16-49	\$130.00
50+	\$125.00



Horizontal Filing



Vertical Filing



Records & Mammo Organizers

Ideal for Medical Records, Mammo Film, CT, Ultrasound and Nuclear Medicine. Increases file room and transportation efficiency. 2-pocket design is available in six (6) colors with frosty clear pockets. Overall size: 13.5" x 10.75".

cat#	color	qty./ctn	price
400	Blue	20	\$155.00
410	Red	20	\$155.00
420	Yellow	20	\$155.00
430	Green	20	\$155.00
440	Purple	20	\$155.00
495	Pink	20	\$155.00



Standard Organizers

Allows handling of files without damage to film or jackets. Keeps records and film in one convenient location. Heavy gauge vinyl welded together with heat-seal construction. Sealed chipboard insert for rigidity. 2-pocket design for files and jackets. Holds up to 14" x 17" film. Overall size: 21" x 14.5".

cat#	color	qty./ctn	price
800	Blue	20	\$165.00
810	Red	20	\$165.00
820	Yellow	20	\$165.00
830	Green	20	\$165.00
840	Purple	20	\$165.00



Standard Wall Bins

Non-shielded steel storage solutions. Will hold receptor plates, film jackets, and caddies. Opening is 22"w x 5"d x 14"h.

Standard Version with polycarbonate see-thru front. HIPAA Version with solid steel no-see front.

#80AC
\$299.00

#80HP
\$299.00



Techno-Tuff™ X-Ray Filing Cabinets

The industry standard cabinets for ARCHIVE or ACTIVE film storage.

Fully assembled, tough, and versatile. All steel construction offers maximum strength. Solid back and side panels for security. 15³/₈" vertical shelf clearance. 12-gauge corner posts for vertical strength. Dividers are 6" apart. Full height dividers add vertical strength and stability. Full depth dividers prevent snagging and ripping of jackets.



36" w x 18" d x 84" h
#618A
\$740.00



36" w x 18" d x 36" h
#632XST
\$445.00



SPECIAL NOTE ABOUT MOTOR FREIGHT:

Please provide information about your receiving area to determine standard or special delivery requirements. Trucking companies will not unload shipments to accounts that do not have a commercial height receiving dock. In most other circumstances, however, special arrangements can be made available at a reasonable extra freight charge.

Tall Receptor/Filing Cart

Protect your investment! "CRC" CR-Plate Units have special, impact resistant padding on shelves and vertical surfaces to prevent damage to expensive CR-Plates. Accommodates film jackets, cassettes, and IR plates. See above for seven colors of durable enamel finish.

Eight 6" wide compartments @ 16.5" vertical clearance. Nominal size: 24" x 18" x 42.5"h. Shipping weight: 130 lbs.

#CRC-38
\$840.00



All CRC Units ship:



Mini Receptor/Filing Cart

Four 6" wide compartments @ 16.5" vert. clearance. Nominal size: 24" x 18" x 36"h. Shipping weight: 100 lbs.


#CRC-24
\$840.00



Electric lift, height-adjustable desks

Straight Front

The **straight front** is designed to resemble a standard desk and is usually placed against a flat wall.




Straight Front

cat.# Length x Depth
ZLSU-E 48"- 96" x 30"- 36"

Contour Front

The **contour front** is designed to keep all of the desk surface within easy reach from any position. This desk is usually placed against a flat wall.




Contour Front

cat.# Length x Depth
ZLCU-E 60"- 96" x 36"

Corner View

The **corner view** is designed to keep all of the desk surface within easy reach from any position. This desk is usually placed in a corner.




Corner View

cat.# Left Wing x Right Wing
ZCVU-E 60"- 96" x 60"- 96" (30" wing depth)

Multi-Purpose Extension

The **multi-purpose extension** is designed to support one-to-two separate viewing or computer systems. This desk can be placed in a corner or against a flat wall.



Right Extension

cat.# Left Wing x Right Wing
ZMPE-ER 60"- 96" (30" deep) x 48"- 72" (24" wide)

Left Extension

cat.# Left Wing x Right Wing
ZMPE-EL 60"- 96" (30" deep) x 48"- 72" (24" wide)



ELECTRIC BASE:

Includes 3-position programmable controller.
25" adjustment range: 25.5" - 50.5" (at 1.5"/sec). Load capacity: 300 lbs. Requires less than 3amps of power. Locking casters available. Limited Lifetime Warranty.

WORK SURFACE:

Laminated solid surface with high endurance PVC edging. Top is supported by a beamless frame that offers generous leg room. Locking casters available. Custom sizes available; please call for quote.



Top Surface
Laminates
(many others
available)



Matte Black



Designer White



Gray Mist



Maple



Persian Cherry




Mahogany

Top Edge Band



Black



White



Gray

Leg/Frame



Elemental Chair & Stool

Comfortably meets all multi-task and office requirements. Back height adjustment. Tilting backrest. Pivot backrest. Tilting seat. Warranty: 7-year limited. Seat: 19"w x 19"d. Back: 19"w x 17.5"h. Stool includes 20" footring.

Ideal for use with sit-to-stand desks.

Chair only:

Seat height range: 16"-21"

Overall height: 35"-43"

Ships @ 101# (actual wt. 52#).

Chair No Arms

#ZOEK-MX03

\$355.00

Stool only:

Seat height range: 23"-33"

Overall height: 42"-55"

Ships @ 68# (actual wt. 46#)

Stool No Arms

#ZOEK-MX03

\$395.00



YES™ High-Back

The Perfect blend of support, comfort, and style. Features contoured lumbar area and mesh-back comfort. Easy back height adjustment. Three-stage locking tilt. Synchronous seat/back tilt. Seat slider. Warranty: 12-year limited. Seat height range: 17.5"-22.5". Seat: 20.5"w x 18"-20"d. Back: 20"w x 24"h with headrest. Overall height: 46"-55". Ships @ 101# (actual wt. 63#)

High-Back +Arms

#ZOYS-HA03

\$775.00

Ideal for desks, workstations & everyday tasks.



Universal Stool

Self-skinned urethane seat and back. Height adjustable footring. Seat height: 22"-32". Overall ht: 37"-47". 7-year limited warranty. Ships @ 37#.

#ZBCH-USB

(Black ONLY, non-upholstered)
\$345.00



UPHOLSTERY FABRICS

Polyester (celestial)						
	1200-Oberon	1201-Rhea	1203-Callisto	1204-Titan	1207-Hyperion	1208-Europa
Polyester (soft sit)						
	2570-Graphene	2572-Ash	2575-Tsunami	2576-Magma	2577-Molten	2578-Kipuka
Polyester (freedom)						
	3400-Orbit	3401-Stealth	3402-Flight	3403-Cocoa	3404-Evening	3408-Oasis
Vinyl (stage)						
	V50-Cats	V53-Grease	V55-Aida	V52-Chicago	V59-Fame	V58-Shrek

Workstations

Workstations feature:

- Omni-directional movement
- Dynamic height adjustment

- Components can be easily relocated along support post.
- Includes: 1-31½" wall mount vertical track
1-extension arm keyboard tray
1-75mm/100mm VESA swivel mount
(6.5-lb min / 17.5-lb max monitor capacity)
- Models with CPU Holder include an additional 31½" track (for 63" total) and CPU holder, adaptable for standard or thin CPU casements.

Mobile Workstation includes:
1-intergrated workstation
1-steel post (total ht.: 67")
1-stabilizer H-base
1-cable channel
4-locking casters

Single View Articulating Arm

Mobile Workstation

#ZTMW-ES21

\$1375.00

Wall Mount (no base)

#ZEVW-AS21

\$880.00



Single View Articulating Arm w/ CPU

Mobile Workstation

#ZTMW-ET31

\$1625.00

Wall Mount (no base)

#ZEVW-AT31

\$1140.00



Desktop Articulating Monitor Supports

Economy Articulating Arm

Cat.# ZEMS-SV21

\$280.00



Articulating gas spring directional adjustment.
21" extension range.
13" height range.
Weight capacity:
6.5-lb min./17.5-lb max.
Pivot/tilt/rotate monitor.
Clamp mount.

Dual Economy Articulating Arms

Cat.# ZEMS-SV22

\$480.00



Similar to the ZEMS-SV21 to the left.
Supports two monitors on independent gas spring articulating arms. **COMBINED** weight capacity: 35-lbs.

Fusion-Lift

NO Assembly! NO Tools! Right out of the box you're ready to convert your existing desk into a Sit-to-Stand desk!


The Upper Desk is 35 $\frac{3}{8}$ " x 15 $\frac{3}{4}$ " and will raise your workspace to a height of 20" with infinite height adjustments throughout its vertical range. The Keyboard surface is 26 $\frac{3}{8}$ " x 13" and provides an ergonomic +/-10° tilt. Smooth, gas-assist movement offers a weight capacity of 35-lbs. Easily mobile at only 40-lbs.

#ZKDR-FL

\$395.00



**Make EVERY DESK
a Sit-to-Stand Desk!**



Trackless Keyboard Support

Spring assisted height adjustment. 5 $\frac{3}{4}$ " height adjustment on a trackless direct mount plate: -1 $\frac{1}{4}$ " to -7" under the surface. +10° to -15° tilt and 360-degree swivel. 19 $\frac{1}{2}$ " wide platform accommodates keyboard. Mouse platform can be placed on either side of the keyboard. Integrated gel palm rest.

#ZKKV-TLS1

\$299.00



Keyboard Supports are not necessary on desks with dual surfaces.

Cable Management Channels

Pair of 24" "J" channels accommodate multiple cables and wires. Cables are quickly and conveniently stuffed into the channels for out-of-sight capture. Easily installs to the underside of desk with self-adhesive strips or screws (included).

#ZACM-CC

\$44.00/pair



Universal CPU Caddy

Stores CPU underneath work surface. Features easy-access slide and swivel action. Adjustable width: 2 $\frac{1}{2}$ " to 8 $\frac{1}{4}$ ". Adjustable height: 13" to 21". 100 lb weight capacity.

#ZEDA-CPUU \$165.00



Mobile CPU Caddy

Stores CPU above the floor for easy access. Fits any CPU within a 3 $\frac{1}{2}$ "-8 $\frac{3}{4}$ " width.

#ZEDA-CPUC \$99.00

Locking Casters

#ZAC3-KE (includes 4 casters)

\$110.00/set of 4

Set of four 3" locking casters. Changes desks into mobile platforms. Easy installation; can be retrofitted. Not recommended for three-legged desks.

Add 3" to your desk height!



INDEX

index by product number

122E.....	94	BMS-95.....	41	DYG.....	11	JCD-CSW.....	113	MFP-XX.....	69	RLC-XX.....	87	TCR.....	9	XSPF-11.....	107
1DR.....	99	BPS-20.....	43	DYL.....	11	JCL-???	112	MKA-19.....	103	RO-107?	100	TDP-????	30	XSPF-12.....	107
1E.....	94	BPS-22.....	43	DYR.....	11	JCM-???	112	MKB-12.....	103	ROW-series	100	TI-SET.....	68	XSVL-21.....	107
2DR.....	99	BPS-30.....	42	EBG.....	4	JJL-DPS.....	110	MKC-15.....	103	S24DR.....	98	TM-51323.....	47	XSVL-22.....	107
2E.....	94	BPS-32.....	42	EBL.....	4	JJL-SPS.....	110	MKG-13.....	103	S25DR.....	98	TMG.....	9	XTIF-3658.....	106
2E-Y.....	102	BRF-45.....	40	EBR.....	4	JJM-BBM.....	111	MKR-11.....	103	S30DR.....	98	TML.....	9	XTSF-????	106
2E-Z.....	102	BSG.....	10	EGG.....	6	JJM-SPU.....	110	MKY-14.....	103	SBC-??.....	72	TMR.....	9	XTUA-??.....	106
2TA-Y.....	102	BSH-65.....	40	EGL.....	6	JJM-PPM.....	111	MM-XRAY.....	94	SBF-??.....	72	TNG.....	9	XTUB-??.....	106
2TA-Z.....	102	BSL.....	10	EGR.....	6	JLL-???	111	MM-XRAY10.....	94	SBG-??.....	72	TNL.....	9	XTUF-??.....	106
35LCH.....	32	BSR.....	10	EMP.....	24	JLS-???	111	MM-XRAY25.....	94	SBK-??.....	73	TNR.....	9	XVFW-????	106
3DR.....	99	BW-??.....	85	EMS.....	24	JML-SBC.....	113	MM3007.....	105	SBM-100.....	73	TPG.....	9	YB??.....	65
3E-Y.....	102	BWS-??.....	85	ETX-70.....	114	JML-SBD.....	113	MPA-19.....	103	SBP-??.....	72	TPL.....	9	YC??.....	65
3E-Z.....	102	CAL-21.....	86	EVG.....	4	JMM-SBC.....	113	MPR-12.....	103	SBS-??.....	72	TPR.....	9	YF??.....	65
3TA-Y.....	102	CBG.....	10	EVL.....	4	JMM-SBD.....	113	MPW-13.....	103	SDG-MSD3.....	105	TVG.....	9	YG??.....	70-71
3TA-Z.....	102	CBL.....	10	EVR.....	4	JNL-OEK.....	112	MPY-11.....	103	SDM-A1.....	104	TVL.....	9	YK??.....	65
400series.....	114	CBR.....	10	FB-10N.....	93	JNM-OEK.....	112	MRR-22.....	103	SDM-BB?.....	104-105	TVR.....	9	YS??.....	65
40LDH.....	32	CCC-DR.....	39	FB-10R.....	93	JNS-OEK.....	112	MRW-23.....	103	SDM-CT1.....	105	TW-200.....	92	YV??.....	66-68
4E.....	94	CCD-DR.....	39	FB-10W.....	93	KGG.....	13	MRY-21.....	103	SDM-CT2.....	105	TWG.....	9	ZCVU.....	116
4E-Y.....	102	CCH-24.....	39	FB-20N.....	93	KGL.....	13	MSB-NF??.....	51	SDM-ASD3.....	105	TWL.....	9	ZAC3-KE.....	119
4E-Z.....	102	CCH-26.....	39	FB-20W.....	93	KGR.....	13	MWG.....	8	SDM-TSD3.....	105	TWR.....	9	ZACM-CC.....	119
4TA-Y.....	102	CCL-14.....	39	FB-40N.....	93	KK-???	13	MWL.....	8	SDM-W??.....	104	TYG.....	9	ZBCH-USB.....	117
4TA-Z.....	102	CCL-17.....	39	FB-40W.....	93	KMG.....	7	MXG.....	8	SDW-SS3.....	104	TYL.....	9	ZEDA-CPUC.....	119
5E.....	94	CCM-10.....	39	FBG.....	4	KML.....	7	MXL.....	8	SDW-SSD3.....	105	TYR.....	9	ZEDA-CPUU.....	119
5E-Y.....	102	CCM-12.....	39	FBL.....	4	KMR.....	7	NCP-24?.....	44	SED-XXXX.....	90	UGP-UL1.....	47	ZEMS-SV21.....	118
5E-Z.....	102	CCS-08.....	39	FBR.....	4	LER-6???	17	OGG.....	5	SFG.....	7	UGP-UL3.....	47	ZEMS-SV22.....	118
5TA-Y.....	102	CCS-10.....	39	FLB-series.....	31	LER-9???	16	OGL.....	5	SFL.....	7	UGP-ULD.....	47	ZKDR-FL.....	119
5TA-Z.....	102	CCZ-14.....	39	FLS-78.....	31	LER-A???	16	OGR.....	5	SFR.....	7	UGW-101.....	47	ZKKV-TLS1.....	119
6EA.....	94	CCZ-17.....	39	FVG.....	4	LER-C???	16	OMG.....	15	SGG.....	11	UGW-102.....	47	ZLCU.....	116
6EB.....	94	CGG.....	10	FVL.....	4	LER-D???	17	OML.....	15	SGL.....	11	UGW-UL1.....	47	ZLSU.....	116
6E-Y.....	102	CGL.....	10	FVR.....	4	LER-P???	18	OMR.....	15	SGR.....	11	UGW-UL2.....	47	ZMPE.....	116
6E-Z.....	102	CGR.....	10	GBR-L.....	12	LER-T???	17	ORG.....	5	SKG.....	13	USB-012.....	46	ZOEC.....	117
6TA-Y.....	102	CMF-01.....	45	GBR-M.....	12	LER-V???	18	ORL.....	5	SKL.....	13	USC-025.....	46	ZOES.....	117
6TA-Z.....	102	CMM-06.....	45	GBR-S.....	12	LES-53??	18	ORR.....	5	SKR.....	13	USC-100.....	46	ZOYS.....	117
618A.....	115	COP-1.....	100	GBS-3.....	12	LES-AF??	19	PAD-51.....	76	SLU-CT.....	108	USC-500.....	46	ZKDW-ES21.....	118
632XST.....	115	COP-2.....	100	GDR.....	12	LES-C???	19	PAD-R2.....	76	SLU-MR.....	108	USO-025.....	46	ZTMW-ET31.....	118
7E.....	94	CP-13.....	101	GFG5.....	12	LES-RXPL.....	19	PAD-S1.....	76	SLU-XL.....	108	USO-048.....	46	ZTVW-AS21.....	118
800series.....	114	CP-14.....	101	GFL5.....	12	LGG.....	11	PAD-SP??	77	SLU-XR.....	108	USO-100.....	46	ZTVW-AT31.....	118
80AC.....	114	CP-34.....	101	GFR5.....	12	LGL.....	11	PAH-????	77	SRG.....	5	USO-500.....	46		
80HP.....	114	CP-35.....	101	GMC-21.....	91	LGR.....	11	PBMB-3024.....	26	SRL.....	5	USP-????	46		
80LL.....	31	CP-36.....	101	GMC-22.....	91	LGS-09.....	86	PBMB-3060.....	26	SRR.....	5	VBD-20.....	108		
A24DR.....	98	CP-37.....	101	GME-60.....	102	LGS-B11.....	86	PBMB-3630.....	26	STG.....	11	VBG.....	6		
A25DR.....	98	CPH-LD12.....	31	GMI-65.....	102	LKG.....	13	PBMB-TILT.....	26	STL.....	11	VBL.....	6		
A30DR.....	98	CPH-LL12.....	31	GMG4.....	12	LKL.....	13	PBMW-2472.....	27	STR.....	11	VBM-24.....	108		
A-3.....	102	CPP-40.....	101	GML4.....	12	LKR.....	13	PBMW-3072.....	27	TA-001.....	95	VBR.....	6		
ACB-28.....	15	CPP-42.....	101	GMR4.....	12	LNA-20.....	104	PBMW-3672.....	27	TA-002.....	95	VBS-10.....	108		
AGR-105.....	21	CPP-60.....	101	GNG.....	12	LPG.....	10	PFC.....	14	TA-1.....	95	VFG.....	7		
AGR-110.....	21	CPP-61.....	101	GNL.....	12	LPL.....	10	PLF.....	14	TA-122.....	95	VFL.....	7		
ARG-06.....	20	CPT-R115.....	99	GNR.....	12	LPR.....	10	PMA-C4.....	38	TA-13.....	94	VFR.....	7		
ARH-01.....	20	CPT-R45.....	99	GPC-23.....	91	LPS-????	28-29	PMA-CB.....	38	TA-14.....	94	VGG.....	6		
ARH-10.....	20	CRA-15.....	45	GSG.....	12	LTG.....	11	PMR-????	79	TA-16.....	95	VGL.....	6		
ARM-100.....	20	CRA-16.....	45	GSL.....	12	LTN.....	11	PPAD-???	79	TA-17.....	95	VGR.....	6		
ARM-G60.....	20	CRA-17.....	45	GSR.....	12	LTR.....	11	PPAH-??	79	TA-2.....	95	VMG.....	7		
ARM-Z10.....	20	CRA-18.....	45	GSS-49.....	91	MBC-30??	48	PPAM-??	79	TA-4.....	95	VML.....	7		
ARM-ZNH.....	20	CRC-24.....	115	GWP-B12.....	92	MBC-31??	48	PPGE-???	81	TA-40??	95	VMR.....	7		
ARP-02.....	21	CRC-38.....	115	GWP-BB14.....	92	MBD-68.....	48	PMS-???	80	TA-5.....	95	WBC-XXXX.....	36		
ARP-06.....	21	CRH-DFC.....	89	GWP-ESW.....	92	MBE-29.....	50	PPSI-???	82	TA-7.....	95	WMP-XX.....	36-37		
ARP-08.....	21	CRH-MBC.....	88	GSW-WAF.....	92	MBE-33.....	50	PPTO-???	83	TA-80??	90	WTS-1.....	100		
ARW-10.....	21	CRH-VIC.....	88	GWP-WAFCS.....	92	MBG-3PV.....	49	PRL.....	14	TA-A3.....	102	WTS-6A.....	100		
ARW-L5.....	21	CXR-13E.....	44	GXE-55.....	91	MBG-6PV.....	49	PS-???	44	TA-LK.....	95	WTS-6B.....	100		
ARW-R5.....	21	DBG.....	11	GXS-57.....	91	MBG-6SS.....	49	PSL.....	14	TAB-L.....	103	WTS-DR.....	99		
BAS-30.....	43	DBL.....	11	GXX-58.....	91	MBI-70.....	48	PSP.....	14	TAB-R.....	103	WTW-01.....	102		
BAS-40.....	43	DBR.....	11	HCG.....	15	MBI-75.....	48	PSR.....	14	TAG.....	9	WTW-23.....	102		
BBD-63.....	40	DDD.....	11	HCL.....	15	MBM-MD.....	50	PT-FRB.....	96	TAL.....	9	WTW-45.....	102		
BBD-66.....	40	DGG.....	11	HCR.....	15	MBN-13E.....	50	PT-SKW.....	96	TAR.....	9	WTW-67.....	102		
BBH-69.....	40	DGL.....	11	HGL.....	14	MBN-14S.....	50	PT-??.....	96-97	TBB-3272.....	51, 84	WTW-80.....	102		
BBS-61.....	40	DGR.....	11	HGP.....	14	MBN-15E.....	50	RDF-XX.....	87	TBC-2214.....	68	WTW-A3.....	100		
BEH-93.....	42	DRG.....	11	HGR.....	14	MBN-16S.....	50	RDH-??	87	TBH-2332.....	51, 84	XRR-19?	44		
BEL-90.....	42	DRL.....	11	HOG.....	15	MBS-68.....	48	RDP-???	87	TBN-1872.....	51, 84	XSCP-10.....	107		
BGB.....	15	DRR.....	11	HOL.....	15	MBS-65.....	48	RDS-12.....	87	TBR-???	84	XSEC-11.....	107		
BGB-2.....	15	DSG.....	11	HOR.....	15	MBW-??	49	RDT-??	87	TBR-2P.....	84	XSEC-12.....	107		
BIV-70.....	41	DSL.....	11	HK-series.....	102	MCH-33C.....	33	RDW-??	87	TBR-B90.....	84	XSIC-21.....	107		
BIV-76.....	41	DSR.....	11	IDC-series.....	109	MCH-37C.....	35	RDT-HN.....	87	TBS-ST2.....	84	XSLC-36.....	107		
BIV-77.....	41	DVG.....	11	IDP-33.....	109	MCH-40C.....	34	RDT-SH.....	87	TBS-????	84	XSLG-21.....	107		
BK-MPGR.....	93	DVL.....	11	IDS-series.....	109	MCP-AP.....	50	RHI-AB.....	86	TCG.....	9	XSLG-22.....	107		
BMS-80.....	41	DVR.....	11	JCD-CDX.....	113	MCP-LK.....	50	RHI-CP.....	86	TCL.....	9				

“-series” or “??” in a catalog number denotes multiple products with a common prefix.